

Product datasheet for **SC323903**

LRRC42 (NM_052940) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	LRRC42 (NM_052940) Human Untagged Clone
Tag:	Tag Free
Symbol:	LRRC42
Synonyms:	dj167A19.4
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC (PS100020)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

Fully Sequenced ORF: >OriGene sequence for NM_052940.3
 GGCCGCGGGAGGAGGGAGCCGTCACCGAGGAGCTGCCGCTCGCTGCCCCGGGCAGGGGCA
 CAGTTGCAACCAAGGCAATGTCTTACTACCTCAGCTCAGAAAACCACTGGACCCAGGGC
 CCATCTACATGCGAGAAAATGGGAGCTGCACATGGTCAATCTGGCTCTGGATGGTGTCA
 GGAGTAGCCTGCAGAAGCCAAGGCCTTTAGACTGTTCCCAAGGCTTTCTGTGGAGC
 TTTGCATGAACAGGGAAGACGACACTGCACGAAAGAGAAGACTGATCATTTTCATCTTCA
 CATACACCCGAGAGGGGAATCTTCGGTACTCCGCCAAATCCCTCTTCAGCCTTGTCTGG
 GTTTCATCTCCGACAATGTGGATCACATTGATTCCCTTATTGGCTTTCCTGAGCAGATTG
 CTGAAAAGCTGTTCTCTGCTGCTGAAGCCAGACAGAAATCACTGAGCCAGGTGCAGGGC
 TGAGGGCTTTACAGAAATCACTGAGGCCTATGGAAGTTTGGTGTCTTGTCCCTGTGTT
 TGCGAAACAGGTATCTCGTGATTTAGAAAAGCTTGAGGAGATTAAGTCTTCCGGGAGC
 TGACCTGCCTGGATCTTCTGTTGCAAGCTTGAGATGAGCATGAAGTCTAGAACATC
 TCACCAATGAAGCCCTGTCTAGTGAATCACTCAGCTCCACCTGAAGGATAATTGTTTATCTG
 ATGCTGGGGTGCAGGAAGATGACAGCACCAGTTCGAGTGATGAAAAGAGGCCTTGAGAATC
 TAACATTATTAGACTTATCATGTAACCCTGAGATCACAGATGCAGGCATTGGATACCTCT
 TTTCTTTTAGGAACTAACTGCTTAGATATCTCTGGGACAGGGCTCAAGGACATCAAAA
 CCGTCAAGCACAAGCTCCAGACCCACATAGGCCTTGTTCACTCCAAAGTGCCCTTGAAGG
 AATTTGATCATAGTAAGTCAAGACAGAGGGCTGGGCTGACCAGATCGTTCTGCACTGGG
 AGCGTGTGACTGCGGAAGCTGTGAAGCCACGGGAGACCTCGGAGCCTAGAGCAGCAGCTC
 AGCGCTTCTATGGGAAGCGGTCTCGAGCAGAAGCCCCACTGAAGTGTCCCTGGCAGACA
 CCCACATGAAGTCTTCCGAGAACTCCAGTTCTATAAGAGAAAGCCCCAGATTGCCATG
 GGCCAGTGTTGAAACACGAAGCTATCTCAAGCCAGGAGTCAAAGAAGAGCAAGAAGAGAC
 CTTTGGAGGAGTCAGAGACAGAACAGAACTCTCAACCTTCAACCTTCAAGCAGAAATATG
 TATGTCTTGCTGTGGAAGACTGGGACTTGTTAAATTCCTATTGATTAGTATACAAAGTT
 GACCTTTCTCTGCCCCAGCTCTAGTGTGTTGAGTAAAGGAGACTGAGGATGATTACTT
 TTTGTTTGAATTTACCTATGGCATTAGCATAGTAAGCAGATATTTCTACTTTTGTGGTGT
 GGGAGGGGAATGCTATGGAAGTATGTAATAATTTCTAATGTATATTTCAATACATAGTA
 ATTTTTTCAAATAAACATTTTGTGCTGTCCTTAAAAAAAAAAAAAAAAA

Restriction Sites: ECoRI-NOT

ACCN: NM_052940

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_052940.3](#), [NP_443172.1](#)

RefSeq Size: 1714 bp

RefSeq ORF: 1287 bp

Locus ID: 115353

UniProt ID: [Q9Y546](#)

Cytogenetics: 1p32.3

Domains: LRR