

Product datasheet for **SC324068**

MRFAP1L1 (NM_203462) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MRFAP1L1 (NM_203462) Human Untagged Clone
Tag:	Tag Free
Symbol:	MRFAP1L1
Synonyms:	PP784
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC (PS100020)
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>OriGene sequence for NM_203462.1

```

GCGATTCTTCCCGCAGAGTTGTGAAGCGAAAGGCTTACAATTAAGGAAGAAAAA
TAAAGATAATTCGGGAGTACAATTGACAAAGCGTGTGGTCTGCTCAGCCTCCAGCAGTAA
CTGCTGATCTCCAGTTCTTGGAGGTTCCGGTGAGAAGAACGCCCTACTGCGTACTGA
GGAAGCGCGAGGAGGAGATGCGGCCCTGGACATAGACGAGGTGGAAGCGCCTGAGGAAG
TGGAGGTGCTGGAGCCGAGGAGGATTTGAGCAGTTCCTGCTCCCGGTCATCAACGAGA
TGCGCGAGGACATCGCGTCTTTATACGCGAGCACGGCGGGCGTACCTGCGGACCAGGA
GCAAGCTGTGGGAGATGGACAATATGCTTATCCAGATCAAAACGCAGGTGGAGGCCTCGG
AGGAGAGCGCCCTCAATCACGTGCAGCACCCGAGTGGCGAAGCCGACGAGAGAGTGTGG
AGTTGTGCGAGAAGGCTGAGGAGAAAGCCAAGGAGATTGCGAAGATGGCAGAGATGCTGG
TCGAGCTCGTCTGGCGAATAGAGAGAAGCGAGTCTTCTGAAGGAGGAGATCGGTGGTTC
ACAGCTGAGGACTGTAGTCAACTGGTAAGGATGAAGTACACCTTGAGGAACTGAAACA
GCTTGTCAATTTCTCTGTTTTGTATGTTTTTACCCCATGTAACAGTCTCCCT
TATGGTCAGTGATTGATGACCTCGATATGGATTAGATTATCAATGTGTTTGGTTTGG
AAATATAACTTTTGTCAAAGAACTACTCTCAGAAGAGAAATGGGGCTTAATTAAGTTGT
TTTGTGGTCACGTTTATCTTGTACTTCGCTGTGTTTTGAAATGTTGGGCATGGCCTC
GTATTTTGCTGTTACCTTTGTGACCTGATTGTTTTTGGAAACACGTCAAGACGTGGGATC
AGAATCTTCCAACCTTAGAGGTGCAATGGAAGACACTACGCTACTTGGTTGAGCCTGGTG
AAGAATGTATTAATGAGACTGCTTTGCATAAACTGGGAAGAAAGAGAAGACAGTTGGAG
ATGGAAGATGGTTTTGTATATATTTTGGAACTTTAGTTCTCTGTGAGACGAAAGAGGAG
AGCTATGTTTTGTGCACATTGTCTGATATATATTGTGAACCTGTCAGGTGAGTTGATT
TAGACAACATAGCTGACCTTTTATGACAAGGCAGTTTGAATAGGGACTATTGTAATACC
TCACACATTATAGGGGCATCAGAGAATGGCATGGAAGAGACAGTCTACAGAGAGCTTTAA
GAGGCCGAGAAAGGAAAAGACATTATCAGGGCCTGGAAGTCTCTCCAGTTCATCAGG
GTAGTAGACCTGTCAGATTGGAGGTCAAAGTCAAAGTATAATAAGTGTTAATTTGCCTT
GAAAAACAAAGACCTAAAAAGCCTTTTTTAAGAAATAAATTTTTTGTTCACATGACCAA
GCCCTCTGTGTGTGTTAATAAAGAATCATAAACAAAAA
  
```


[View online »](#)

Restriction Sites:	ECORI-NOT
ACCN:	NM_203462
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_203462.1</u> , <u>NP_982287.1</u>
RefSeq Size:	1616 bp
RefSeq ORF:	384 bp
Locus ID:	114932
UniProt ID:	<u>Q96HT8</u>
Cytogenetics:	4p16.1