

## Product datasheet for **SC324085**

### TCTN1 (NM\_024549) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	TCTN1 (NM_024549) Human Untagged Clone
Tag:	Tag Free
Symbol:	TCTN1
Synonyms:	JBTS13; TECT1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC (PS100020)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >OriGene sequence for NM\_024549.4  
 GCAACGCGCTGTCCATGTCGCGGGCCTCGCTGGGACTCCCTGGGAGATGAGGCCGCGAGG  
 TCTCCCGCCGCTCCTGGTGGTGCCTGGGCTGCTGGGCTCCGTGAGCGCCAGACCGA  
 TGCCACCCCGGCGGTGACGACAGAGGGCCTCAACTCCACCGAGGCAGCCCTGGCCACCTT  
 CGGAACCTTCCCGTCGACCAGGCCCGGGGACTCCCAGGGCTCCAGGGCCCTCCTCCGG  
 CCCCAGGCCTACCCAGTCACGGACGTTGCTGTTCTCTGTGTCTGTGACTTATCCCCAGC  
 ACAGTGTGACATCAACTGCTGTGTGATCCCCGACTGCAGCTCCGTGGATTCAGTGTCTT  
 TTCTGCCTGCTCAGTTCAGTTGTGACGGGCGACAGCCAGTTTTGTAGTCAAAAAGCAGT  
 CATCTATTCAATTGAATTTTACAGCAAACCCACCTCAAAGAGTATTTGAACTTGTGACCA  
 GATTAATCCATCTATTTTCTGCATTATATTACAAACTATAAACCTGCATTATCCTTTAT  
 TAATCCAGAAGTACCTGATGAAAACAATTTTGATACATTGATGAAAACATCTGATGGTTT  
 TACATTGAATGCTGAATCATATGTTTCCTTACAACCAAAGTGGATATTCCTACTGCTGC  
 TAAATATGAGTATGGGGTTCCTCTGCAGACTTCAGATTCGTTTCTGAGATTCCTTCGTC  
 CCTGACATCATCTCTGTGACTGATAATAACCCTGCAGGCCAGGCGTACTGGTTCACACC  
 TGTAATCCAGCACTCTGGGAGGCCGAGGCGAGAGGATCACTTGAGGTACCTGATTCAAG  
 AAAAAAGGTCCCTATCACTGTTCACTCCATCGTCACTTCACTCTCTAAATAAACGCTCAC  
 CCGACGGGAGGACACTGATGTGCTGCAGCCGACTCTCGTCAACGCTGGACACTTTAGCCT  
 TTGCGTGAATGTTGTTCTTGAGGTAAAGTACAGCCTCACATACACAGATGCAGGTGAAGT  
 CACCAAAGCTGATCTCTCATTGTTCTGGGGACAGTTAGCAGCGTAGTGGTCCCCTGCA  
 GCAAAAGTTTGAATTCATTTTCTTCCAGGAAAATACCCAGCCAGTCCCTCTCAGTGGAAA  
 CCCTGGTTATGTCGTGGGGCTCCATTAGCTGCTGGATTCCAGCCTCATAAGGGGTCTGG  
 GATTATTCAGACCACAAATAGATATGGACAGCTTACTATTCTTCATAGCACAAGTGAACA  
 AGACTGCTTAGCACTGGAGGGGCTCCGGACCCAGTATTTTGGTTACACTATGCAATC  
 TGGCTGTAAACTAAGACTGACTGGAGCTCTCCCGTGTCACTCGTAGCACAGAAGGTGAA  
 GAGCCTGCTGTGGGGCCAGGGCTTCCCAGATTACGTGGCCCTTTTGGAAATTCAGGCG  
 CCAGGACATGCTGGACTGGGTGCCATCCACTTCATCACCCAGTCATTCAACAGGAAGGA  
 TTCCTGCCAGCTCCAGGGGCTTTGGTTATAGAAGTGAAGTGGACTAAATACGGATCCCT  
 GCTGAATCCACAGGCCAAAATAGTCAATGTAAGTGAAGTGAAGTGAAGTGAAGTGAAGT  
 TGAGGCCAACTCAGGAAATGAAAGGACGATTCTTATTTCCACTGCGGTTACTTTTGTGGA  
 TGTGTCTGCACCTGCAGAGGCAGGCTTCCAGAGCTCCACCAGCCATCAATGCCAGGCTGCC  
 CTTTAACTTCTTCTCCCGTTTGTGTTGACAATGCTCAGATGCATCAGTTCCTTAATATAC  
 ACGTGAATTTGAAAAGTACATTCGGTGAGATTAATTTTATATACAATAAAAAAAAAA  
 AAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAA  
 AAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAA

**Restriction Sites:** ECoRI-NOT

**ACCN:** NM\_024549

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_024549.4](#), [NP\\_078825.2](#)

**RefSeq Size:** 1935 bp

**RefSeq ORF:** 1722 bp

**Locus ID:** 79600

**UniProt ID:** [Q2MV58](#)

**Cytogenetics:** 12q24.11

**Protein Families:** Secreted Protein

**Gene Summary:** This gene encodes a member of a family of secreted and transmembrane proteins. The orthologous gene in mouse functions downstream of smoothed and rab23 to modulate hedgehog signal transduction. This protein is a component of the tectonic-like complex, which forms a barrier between the ciliary axoneme and the basal body. A mutation in this gene was found in a family with Joubert syndrome-13. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Feb 2016]

Transcript Variant: This variant (3) lacks an exon and contains an alternate in-frame exon in the central coding region and an alternate in-frame splice site in the 3' coding region, compared to variant 1, resulting in a shorter protein (isoform 3).