

Product datasheet for **SC324132**

TMEM80 (NM_174940) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TMEM80 (NM_174940) Human Untagged Clone
Tag:	Tag Free
Symbol:	TMEM80
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC (PS100020)
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>OriGene sequence for NM_174940.2

```
TCCGAAAGTGCCCGAGCGGTGCCGACGGACTAATCGGGCCTCGGCCGTGGCTCGGACGT
CCGCTCCCAGAACGCGGGCGGGTCCGGCCCTTCGTGAGGAGACGCGAAAATGGCCGAAG
GAGCGCGAGCGCGGGCCGAGAGGCTGCCGGGATCGCGACGGACCGGCCGGCGGGCGG
GGAGAGGATCCTCCACAGTGCTCTTTCAGTTCCTTCAAATGCTGTTTTATCTCAGCG
GAACGTAACGCCCTGTATTTCTCGCCACGCTCCTGATGATCAGTATAAAAGTCAGG
TGTTTCAGTATCCTCACCCTACCTGGTCTCGATCTTGCTCTGCTGTTTCTGATGGGGA
TTCTAGAAGCAGTTTCGGTTATACCTGGGACCAGGGGCAACCTGACAGAGGCTGAGAGGC
CGCTGGCCGCCAGCCTGGCCCTCACGGCTGGCACCGCCCTCCTCTGCCCCACTTCCTGC
TTTGGCAGGCCCTAGTGTGTGGGCGGACTGGGCCCTCAGCGCCACGCTCCTGGCCCTTC
ACGGCCTGGAGGCCGCTCCTGCAGTGGTTGCCATCGCGGCCTTACCAGGTAGCTACGGA
CACCCGGGATACCCACACTGGGGCCCTCCTCCTGGGCCTGACCAGTCCCCCAGCTGTCA
CCTCCCCATTCTGGACAGGAAGGGCACTTTTCTAGTGACTGGCCATAGATGGTTTTGG
ATGTTTCCATCTGTTCTGGCAGGAGTGGGAGCAGGAGCCAGGGCAGAACAACACTGCTGGA
GGCCCTGGTGTGGGAACAGCTGCGGGGAGGGTAGGGACCAGACAGAAGTGCCTTCAAGA
TGAGTCCCAGGAGCGCACACTCAGCCCTGTGAGTGGGGTCTGGCTTTAGCAGCCAGGCCCT
CCACAGACCCCCATGGGCCCCAGGGCCGAGAGGGAGGACAGAGCCCTTCAAGACAGAGG
CCTCATCTCACTGCATCCCCATCACCCCTAGTTCCTCAATGGTCTAATTTGTGTTCT
GAGATCCCAGTTTACTCTGTGGCCAGGCCCCACCTGTGTTTCCAAGTCGGGCTGGAGACG
CAGGATGGGGTAGGCCTTGTGCTCTGAGCAACCCAGCTCTGCCTCACAGGCAGGCAGGC
CCGGTGCAAGAGTGGACTCTGGGTTCTAAAGCAATAAATGCAAACAAGCCAACAGCTCT
GCTGCCTAGCAATTTCCATCTTAGCCCACTTCTCCCTTCAGGGGCTTCGGAGGAGAGGT
CAGGGCTAAGGCCGGGGATGAGACTGCAGGAGAGAGAGCAGCGGAGGGCCACATTCGGAG
CCTCCGTCCAATCCAGTTTTATCAGTTTTGCCTTTTGCACGGAGTGCTAAACAAATTTCT
AGCTCTGTGTTTTTTCCATTTCCAGATTTACTATCAGTTCTCCTTAAAAAGTATCTAA
GCTGTTACAGTAGCTTTCCCTTCACTTGATTCTATTGTGTGTTTTCTATGTTTGGAAATA
TTACACCCAAATATCTAGATATTTCTCTTCCACCGCATTTTGTAAATAAAGAGATGTGTA
TGCCCTCCCTGAAAAAAAAAAAAAAAAA
```



[View online >](#)

Restriction Sites:	ECORI-NOT
ACCN:	NM_174940
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_174940.2</u> , <u>NP_777600.2</u>
RefSeq Size:	1613 bp
RefSeq ORF:	483 bp
Locus ID:	283232
UniProt ID:	<u>Q96HE8</u>
Cytogenetics:	11p15.5
Protein Families:	Transmembrane