

Product datasheet for **SC324231**

YTHDC1 (NM_001031732) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	YTHDC1 (NM_001031732) Human Untagged Clone
Tag:	Tag Free
Symbol:	YTHDC1
Synonyms:	YT521; YT521-B
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC (PS100020)
E. coli Selection:	Ampicillin (100 ug/mL)

Fully Sequenced ORF: >OriGene sequence for NM_001031732.2
CGCCGTCGCCATCTGTTTCCCTCCCTCCCCATACCTAGCCCGAGTCTGAGCCCTAACG
AGAAGGCTGGCCCTAGGCCCTGGATGCTGGAGTGAAAGGAAGGGAGAAAGGGAAAAAGC
GGGAAGAGTCGAGAAGGCTGAGTGTTAAGAGGCCAAGTGCACGCGCGTATCCGGGCAGA
CGGACTGACGGACGGGCCGTGCTTCTGCCGCGGCTGCGGCGCCCGCGGAGTTCGCGTCT
AAGCGGCGGCGGCGGTGGCAGCGCGGAAACCGAAGGGGAGCCATGGCGGCTGACAGTCG
GGAGGAGAAAGATGGAGAAGTAAATGTTCTGGATGATATTTAACTGAAGTACCAGAACA
AGATGATGAAGTATAATCCAGAGAGTGAACAAGATAAAAAAGGAAAAAGGGATCAAA
AAGAAAAAGTGAATGGAATCTACTGATACCAAACGACAAAAGCCTTCTGTCCATTC
AAGACAAGTGGTTTCTAAGCCACTGAGCTCATCTGTTAGCAATAACAAAAGAATAGTTAG
TACAAAAGGAAAGTCAGCCACAGAGTATAAAAAAGGAAATATCAAAGATCTGAAAGAAA
CAAGCGTCTAGATGCTGATCGGAAAATTCGTCTATCAAGTAGTGCCTCCAGAGAACCTTA
TAAGAATCAACCTGAAAAACCTGTGTCCGAAAAGGGATCCTGAAAGGAGGGCCAAATC
TCCTACGCCAGATGGTTCTGAGAGAATTGGGCTTGAAGTGGATAGACGTGCAAGCAGATC
CAGCCAGTCTTCTAAGGAAGAAGTGAAGTCTGAAGAATATGGCTCTGACCATGAGACTGG
CAGCAGTGGTTCTTCTGATGAGCAAGGGAACAACACTGAGAATGAGGAGGAAGGAGTGG
AGAAGATGTGGAGGAAGATGAAGAAGTGAAGAAGATGCAGAAGAAGATGAAGAGGTGGA
TGAAGATGGAGAGGAGGAGGAGGAAGAGGAGGAGGAGGAGGAGGAGGAGGAGGAGGAGG
AGAAGAAGAATATGAACAGGATGAGAGAGACCAGAAAAGGAGGAGGAAATGATTATGACAC
TCGAAGTGAGGCCAGTACTCTGGTTCTGAATCTGTTTCCCTTACAGATGGGTCTGTGAG
ATCTGGTTCAAGCACAGATGGATCAGATGAGAAAAAGAAGGAAAGGAAGAGAGCTAGAGG
CATATCTCCAATTGTTTTGATAGAAGTGAAGCTCTGCATCAGAGTCATATGCAGGTTT
AGAAAAGAAGCATGAGAAATTATCATCTTCCGTTCTGCTGTCCGAAAAGATCAAACCG
TAAACTCAAATATGTGCTTCAAGATGCAAGATTTTCTCATAAAGAGTAACAACCATGA
GAATGTGCTCTTGCCAAAGCGAAGGGTGTATGGTCCACGCTCCCTGTAATGAGAAGAA
ATTAATCTTGCAATTTAGATCTGCAAGGAGTGTATCTTAATATTTCTGTGAGAGAGAG
TGGAAAATTTCAAGGTTTGAAGACTTTCTTCAAGTACATCAGGAGGATCTCTAT



[View online »](#)

```

ACACTGGGTGCTTCCAGCAGGAATGAGTGCTAAAATGCTGGGAGGTGTCTTTAAAATTGA
CTGGATTTGCAGGCGTGAATTACCCTTCACTAAGTCGGCTCATCTCACCAATCCTTGAA
TGAACATAAAACCAGTAAAGATCGGACGTGATGGACAGGAAATTGAACTTGAATGTGGAAC
CCAGCTTTGTCTTCTGTTTCCCCCGATGAAAGTATTGACTTGTATCAGGTCATTCATAA
AATGCGTCACAAGAGAAGAATGCATTCTCAGCCCCGATCACGAGGACGTCCATCCCGTCG
AGAACCAGTCCGGGATGTGGGAAGGCGTCGACCAGAAGATTATGATATTCATAACAGCAG
AAAGAAACCAAGGATTGACTATCCCCCTGAGTTTCACCAGAGACCAGGGTATTTAAAGGA
TCCACGATACCAGGAAGTGGACAGACGATTTTCAGGAGTTCGCCGAGATGTGTTTTTAAA
TGGGTCTACAATGATTATGTGAGGGAATTTCAACATGGGACCACCACCACCTTGGA
AGGAATGCCCCCTTACCAGGAATGGAACAACCTCCACACCATCCTTACTATCAGCACCA
TGCTCCACCTCTCAAGCTCATCCCCCTTACTCAGGACATCATCCAGTACCACATGAAGC
AAGATACAGATAAACGAGTACATGATTATGATATGAGGGTGGATGATTTCTTCGTCG
CACACAAGCTGTTGTCAGTGGCCGGAGAAGTAGACCCCGTAAAGAGACCCGGGAACGAGA
GCGAGACCGCCCTAGAGATAACAGACGAGACAGAGAGCGAGATAGAGGACGTGATAGAGA
AAGAGAAAAGAGAGCGATTATGTGATCGAGACAGAGACCGAGGGGAGAGAGGTCGATATAG
AAGATAATGGGCTTTTGAAGCACTGATTGTTTAAAGATACAAAAAATCTTGATTTTTT
TTTTGTGTGTGTTTACAAGTAGTAAATTTATTTTCAGCTGTCTGCCTATGAAGTTCATTG
TGTAGAAGGATTTATTATGACCCCTTTGTTCCAAGCATGCAGTATCATAAGAATCGAA
AAACTCAAATCCGCCAAAAATCCACAGCTGACAGTTGAATTGACACTTTTATTGGGGCAG
AATGGAACAGTCCAAGAATGTAGATACTGATTCTTTCTCCATAAATGTTCTTATAGTGTG
TTCATCTAGAGTTATTTTTTTGTTTTTTTTTTCCTTTTGGATCTTGATTGATAACTG
CCATGATATTTGCTTTGATGTGTTTCTACATGTAGTTGCACACGGTTCAGTAAAAATAA
TGCTGCTATCGAGTATGCAAATATTGAAGTATGATGGTTTGACTGTATGGCAGTGTGTA
GCAGCCTCTGTTTTTTTCCCATTGCCTCTTTTTTAAAAAACTATAAAGTCACTTTT
TATTTTTCTCAGTCTTCAATGACGAGCAATATTAAGAAGACATTGCTATCTAATTTT
TAATCTTTTTAAATGAAAAATTCCTATGTTTCAGTAGCGTGGTTGATGCTATTGTTTAGCC
TTCCCCTCAAATGTATACATTGGCTTGGAAATGTTCACTTGCCTGCGTGGCAGCGG
AAGACGATTCCCATATTCTACATTGCTACTGTTTTGTATAAAAATAAATTGGTAAAGATTG
AAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAA

```

- Restriction Sites:** ECoRI-NOT
- ACCN:** NM_001031732
- OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
- OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
- Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001031732.2](#), [NP_001026902.1](#)

RefSeq Size: 6249 bp

RefSeq ORF: 2184 bp

Locus ID: 91746

UniProt ID: [Q96MU7](#)

Cytogenetics: 4q13.2

Gene Summary: Regulator of alternative splicing that specifically recognizes and binds N6-methyladenosine (m6A)-containing RNAs (PubMed:26318451, PubMed:26876937, PubMed:25242552, PubMed:28984244). M6A is a modification present at internal sites of mRNAs and some non-coding RNAs and plays a role in the efficiency of mRNA splicing, processing and stability (PubMed:26318451, PubMed:25242552). Acts as a key regulator of exon-inclusion or exon-skipping during alternative splicing via interaction with mRNA splicing factors SRSF3 and SRSF10 (PubMed:26876937). Specifically binds m6A-containing mRNAs and promotes recruitment of SRSF3 to its mRNA-binding elements adjacent to m6A sites, leading to exon-inclusion during alternative splicing (PubMed:26876937). In contrast, interaction with SRSF3 prevents interaction with SRSF10, a splicing factor that promotes exon skipping: this prevents SRSF10 from binding to its mRNA-binding sites close to m6A-containing regions, leading to inhibit exon skipping during alternative splicing (PubMed:26876937). May also regulate alternative splice site selection (PubMed:20167602). Also involved in nuclear export of m6A-containing mRNAs via interaction with SRSF3: interaction with SRSF3 facilitates m6A-containing mRNA-binding to both SRSF3 and NXF1, promoting mRNA nuclear export (PubMed:28984244). Also recognizes and binds m6A on other RNA molecules (PubMed:27602518). Involved in random X inactivation mediated by Xist RNA: recognizes and binds m6A-containing Xist and promotes transcription repression activity of Xist (PubMed:27602518). Involved in S-adenosyl-L-methionine homeostasis by regulating expression of MAT2A transcripts, probably by binding m6A-containing MAT2A mRNAs (By similarity).[UniProtKB/Swiss-Prot Function]

Transcript Variant: This variant (1) encodes isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.