

Product datasheet for SC324535

RRP15 (NM_016052) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	RRP15 (NM_016052) Human Untagged Clone
Tag:	Tag Free
Symbol:	RRP15
Synonyms:	CGI-115; KIAA0507
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC (PS100020)
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>OriGene sequence for NM_016052.3

```

CAGAAAAATGGCAGCCGCTCCGGACTCACGTGTGAGTGAGGAAGAAAACCTGAAAAA
GACCCCAAGAAGAAGATGAAAATGGTAACTGGAGCCGTAGCGTCGGTGTGGAAGACGA
GGCCACAGACACTTCTGATAGTGAAGGAAGCTGTGGATCGGAAAAGGACCACTTTTATTC
TGATGATGACGCAATAGAAGCTGACAGTGAGGGTGATGCTGAGCCCTGTGACAAAGAAAA
TGAAAATGATGGAGAATCAAGTGTGGGACTAATATGGGCTGGGCAGATGCTATGGCTAA
AGTCCTCAACAAGAAAACCTCTGAAAGTAAACCTACTATTCTGGTCAAAAATAAGAAGCT
GGAAAAGGAAAAAGAAAGTTAAAGCAAGAAAGACTAGAGAAAATAAAACAGCGTGATAA
GAGGCTGGAGTGGGAAATGATGTGCAGAGTAAAGCCAGATGTTGTCCAAGACAAAGAGAC
AGAGAGAAATCTTCAGAGAATTGCAACAAGGGGTGTGGTGAATTATTTAATGCTGTTCA
GAAACATCAAAAGAATGTTGATGAAAAGTTAAGGAAGCTGGAAGTTCTATGAGAAAGCG
TGCTAAGTTGATATCAACTGTTTCCAAGAAAGATTTTCATCAGTGTGTTGAGAGGGATGGA
TGGAAGTACAAATGAGACTGCTTCAAGCAGGAAGAAACCAAAAGCCAAACAGACTGAAGT
GAAATCAGAAGAAGGCCAGGTTGGACGATCCTACGTGATGATTTTCATGATGGGAGCATC
TATGAAAGACTGGGACAAGGAAAGTGATGGGCCAGATGACAGCAGACCAGAATCTGCAAG
TGACTCTGATACATAAAGCATCATAGGAAATACAATTGCAGTCGTTTTATTTTTCTAGA
AAAATATGTCATCCTCTGATAGTTGGGGAATTATAAGGATACCATTGTGAAGAAAGCCAA
AAGACTTTTGCCAGATTTTCATATTTCCCTTTTCATGTACACTTTATATATACTTCATTA
AAATTATATTTTAAACCCTTGATAATTTTAAGCATTGTTCTCAGAACATTTGTAAAAG
GATATATTTCTGCTTGACCAGCGAGATGTGCATTTTGCCAGGATCATATTGGTCATGTCT
ATTGGTGTATTATTTTCAGTATACCAATGTTTTTCAGAAATACAGTACTAATTCATCATT
AACTCTTTGAAGTTAATATTTTTCTGCCTTCTAACTTATAGACTCAACTATGTATCTGTA
GTTTTTGGGAATGGTGGTGTGTTTTGCTTTGTGTTGGGAAGTTATTGAGAAAACCTATA
TAATAAAATTTAAATTTAAAAA
  
```

Restriction Sites: Please inquire


[View online »](#)

ACCN:	NM_016052
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_016052.3</u> , <u>NP_057136.2</u>
RefSeq Size:	7782 bp
RefSeq ORF:	849 bp
Locus ID:	51018
UniProt ID:	<u>Q9Y3B9</u>
Cytogenetics:	1q41
Gene Summary:	This gene encodes a protein that co-purifies with human nucleoli. A similar protein in budding yeast is a component of pre-60S ribosomal particles, and is required for the early maturation steps of the 60S subunit. [provided by RefSeq, Jul 2008]