

Product datasheet for **SC324653**

UBAP1 (NM_016525) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	UBAP1 (NM_016525) Human Untagged Clone
Tag:	Tag Free
Symbol:	UBAP1
Synonyms:	NAG20; SPG80; UAP; UBAP; UBAP-1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC (PS100020)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

Fully Sequenced ORF: >OriGene sequence for NM_016525.3

```

GGGAAGTAGGACTTCAACATGGCGGCTGCGGCACTGGCGGTGGCTACGGTGACGGCCTGG
CCCGGAGCGGGCAGAGTTGGAGGTGGTGGCGTTTCGCTCTCCCTAGGGGCTGTCCGGAGCT
CAGCGGGGACCGAGCCTGGGAGGCCGCGGTGCCAGCACCTTTTCGGCTTCTGAGACGGC
GGCAGCAGCGGCATTACAGTTCTAAATGGCTTCTAAGAAGTTGGGTGCAGATTTTCATGG
GACTTTTCAGTTACCTTGATGATGTCCTATTAAGACAGGAGACAAATCAAAACACCAGC
TAAAGTTGGTCTACCTATTGGCTTCTCCTTGCCTGATTGTTTGCAGGTTGTCAGAGAAAT
ACAGTATGACTTCTCTTTGGAAAAGAAAACCATTTGAGTGGGCTGAAGAGATTAAGAAAAT
CGAAGAAGCCGAGCGGGAAGCAGAGTGCAAAATTCGGAAGCAGAGCTAAAGTGAATTC
TAAGAGTGGCCAGAGGGCGATAGCAAAATGAGCTTCTCCAAGACTCACAGTACAGCCAC
AATGCCACCTCTATTAACCCCATCCTCGCCAGCTTGACGACACACAGCATCCTCACACC
AACTCGGGTCAGCAGTAGTGCCACGAAACAGAAAGTTCTCAGCCACCTCACATAAAGGC
GGATTTCAATCTTGCTGACTTTGAGTGTGAAGAAGACCCATTTGATAATCTGGAGTTAA
AACTATTGATGAGAAGGAAGAGCTGAGAAATATTCTGGTAGGAACCACTGGACCCATTAT
GGCTCAGTTATTGGACAATAACTTGCCAGGGGAGGCTCTGGGTCTGTGTACAGGATGA
GGAGGTCCTGGCATCCTTGGAACGGGCAACCCTAGATTTCAAGCCTCTTCATAAACCCAA
TGGCTTTATAACCTTACCACAGTTGGGCAACTGTGAAAAGATGTCACTGTCTTCCAAAGT
GTCCCTCCCCCTATACCTGCAGTAAGCAATATCAAAATCCCTGTCTTCCCAAACCTTGA
CTCTGATGACAGCAATCAGAAGACAGCAAGCTGGCGAGCACTTTCATAGCACATCCTG
CCTCCGCAATGGCACGTTCCAGAATCCCTAAAGCCTTCCACCCAAAGCAGTGCCAGTGA
GCTCAATGGGCATCACACTCTTGGGCTTTCAGCTTTGAACTTGGACAGTGGCACAGAGAT
GCCAGCCCTGACATCCTCCAGATGCCTTCCCTCTCTGTTTTGTCTGTGTGCACAGAGGA
ATCATCACCTCCAAATACTGGTCCCACGGTCACCCCTCCTAATTTCTCAGTGTCAACAAT
GCCCCAATGCCCAGCTGTCCCAGGCCCTATTCTGAACTGCAGATGCTGTCCCCCAGCGA
GCGGCAGTGTGTGGAGACGGTGGTCAACATGGGCTACTCGTACGAGTGTGTCTCAGAGC
CATGAAGAAGAAAGGAGAGAATATTGAGCAGATTCTCGACTATCTCTTGCACATGGACA
GCTTTGTGAGAAGGGCTTCGACCCTCTTTAGTGAAGAGGCTCTGAAAATGCACCACTG
TTCAGAAGAAAAGATGATGGAGTTTCTTCAGTTAATGAGCAAAATTAAGGAGATGGGCTT
TGAGCTGAAAGACATTAAGGAAGTTTTGCTATTACACAACATGACCAGGACAATGCTTT
GGAAGACCTCATGGCTCGGGCAGGAGCCAGCTGAGACCAGGCCCTGCCTAGGCCCTGCCG
CAGAACCACCATCCCTGGGAGGCCCTGCAGAGCCACCTGTGGGAAAGAGAAGGGGCAG
CTTCCGGATTTTCTTTTGGGGTTAGAAGGTACAGGTGTGGAGACTGCTCGCCAGTCTCTG
TGAGCCTAGGCCCTGAGCTGGGGAGGTGGGAAGATTGGGCATGTGAGTGCCCCCAGAA
CTGTCTGGCTCCTTCCGTATTAACAGCATTTGCATTTTGAGAAGTGCTTCCCACTTC
AGCCCTCCGGAGAGACTACCCTAGTCTTTCTGGGGTGTATGTCCTCAGCTGAAGCCTG
GCCTAGTTGCTGAGAGGGGCTGGGGAGATGGGGCGGGAGGGCCAGACTCAGTGCTGCTGT
GGAGCTAGGTGCTTCCCCTTCCCCTGAGACTGGTGGACTGAACTCCAGTCAAGTTGAGT
TCAAGTGAAAGATTCTTCAGGGTTTTATTTTTCCCCTCCTAACAAAGTCTCATAGTGT
TAACACTGGTTCTGCAATATCTCTGAGGTGCAAGAATGCACTTTCCCTATGGGGCCCA
GAGTTTGCTTTTCTGCCAGGCAGTCACCACGCTTCCCTACCCAGCCTGTTTCTTTTGG
CTTGTTTTGGACCACAGTCTCTGCTACCCAGGGTTTTAGAGCCCTGCTCTAGGAAACA
GTTTAAGAAATCATTGGCCCTTCCCAGCACATTGAATGGGTAAGCAGACAGGCCATGAT
TTAGTTGGCCAGCACTAACTCCACCTCTGTTCTCCTTGAACAGCTTCCCCTCCAGCCAC
TGCTTTAGGATGACACAATGAATAACACCTAGTCATAGAAATCAGTCTCTCTGGTTTGT
TTGATTATGTTGTACATCATTAAAGATCTAAATACAAAGGAAAAAAAAAAAAAAAAAAAA
AAAAAAA

```

Restriction Sites: Please inquire

ACCN: NM_016525

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_016525.3</u> , <u>NP_057609.2</u>
RefSeq Size:	2757 bp
RefSeq ORF:	1509 bp
Locus ID:	51271
UniProt ID:	<u>Q9NZ09</u>
Cytogenetics:	9p13.3
Domains:	UBA
Gene Summary:	This gene is a member of the UBA domain family, whose members include proteins having connections to ubiquitin and the ubiquitination pathway. The ubiquitin associated domain is thought to be a non-covalent ubiquitin binding domain consisting of a compact three helix bundle. This particular protein originates from a gene locus in a refined region on chromosome 9 undergoing loss of heterozygosity in nasopharyngeal carcinoma (NPC). Taking into account its cytogenetic location, this UBA domain family member is being studies as a putative target for mutation in nasopharyngeal carcinomas. Multiple alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jan 2010] Transcript Variant: This variant (1) is the dominant transcript and encodes the dominant isoform (1). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.