

Product datasheet for **SC324762**

PRELID3A (NM_001142405) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PRELID3A (NM_001142405) Human Untagged Clone
Tag:	Tag Free
Symbol:	PRELID3A
Synonyms:	C18orf43; HFL-EDDG1; SLMO1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC324762 representing NM_001142405. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
 GATCCGGTACCGAGGAGATCTGCCGCC**CGATCGCC**
 ATGAAGATCTGGAGCTCGGAGCAGTGTGGCCACCCGTGGGACACGGTCATCCAGGCGGCCATGCGC
 AAGTACCCGAACCCGATGAACCCGAGCGTGCTGGGCGTGGATGTGCTACAGCGCCGCTGGACGCCGC
 GGCCGCCTGCACAGCTTGCCTGCTCAGCACCAGTGGGGCTGCCAGCCTCGTGAGAGCGATTTTG
 GGAACCACTAGGACATTGACATACATCCGAGAACATTCTGTGGTGGATCCAGTGGAAAAGAAAATGGAA
 CTTTGTCTACCAATATCACACTCACAATTTGGTGTGAGTTAATGAGAGTTGGTGTACACACCTCAT
 CCAGAGAACCCAGAAATGACCGTCTCACACAAGAAGCCATCATCACTGTGAAGGGGATTAGCCTTGGT
 AGTTATTTGGAAAGTTTAAATGGCCAATACGATATCATCCAATGCAAAGAAGGGTGGGCTGCTATCGAG
 TGGATAATTGAACACTCTGAAAGCGCTGTGAGCTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
 TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC

Restriction Sites:	SgfI-MluI
ACCN:	NM_001142405
Insert Size:	519 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).


[View online »](#)

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001142405.1</u>
RefSeq Size:	1734 bp
RefSeq ORF:	519 bp
Locus ID:	10650
UniProt ID:	<u>Q96N28</u>
Cytogenetics:	18p11.21
MW:	19.2 kDa
Gene Summary:	In vitro, the TRIAP1:PRELID3A complex mediates the transfer of phosphatidic acid (PA) between liposomes and probably functions as a PA transporter across the mitochondrion intermembrane space. Phosphatidic acid import is required for cardiolipin (CL) synthesis in the mitochondrial inner membrane.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) represents the more predominant transcript. Transcript variants 1 and 2 encode the same isoform (1). Sequence Note: This RefSeq record was created from transcript and genomic sequence data because no quality transcript was available for the full length of the gene. The extent of this transcript is supported by transcript alignments.