

Product datasheet for SC324776

CITED1 (NM_001144887) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	CITED1 (NM_001144887) Human Untagged Clone
Tag:	Tag Free
Symbol:	CITED1
Synonyms:	MSG1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC324776 representing NM_001144887. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTT TAGTGAACCGTCAGAATTTTGT AATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCC GCGATCGCC
ATGCCAACACGTCGAGGCCTGCACTTGATGTCAAGGGTGGCACCTCACCTGCGAAGGAGGATGCCAAC
CAAGAGATGAGCTCCGTGGCCTACTCCAACCTTGC GGTGAAAGATCGCAAAGCAGTGCCATTCTGCAC
TACCCTGGGGTAGCCTCAAATGGAACCAAGGCCAGTGGGGCTCCCACTAGTTCCTCGGGATCTCCAATA
GGCTCTCTACAACCACCCCTCCCACTAAACCCCATCCTTCAACCTGCACCCCGCCCTCACTTGCTG
GCTAGTATGCACCTGCAGAACTTAATAGCCAGTATCAGGGGATGGCTGCTGCCACTCCAGGCCAACCC
GGGAGGCAGGACCCCTGCAAACTGGGACTTTGGGGCCAGGGGGAGGGGCAGAACTACTCTCTCTCT
TCTGCTGGTGCCAGAGCCCTGCTATCATCGATTCCGACCCAGTGGATGAGGAAGTGCTGATGTCGCTG
GTGGTGGAACTGGGGTTGGACCGAGCCAATGAGCTTCCGGAGCTGTGGCTGGGGCAGAATGAGTTTGAC
TCACTGCCGACTTTCATCTAGCTGCTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: Sgfl-MluI



[View online »](#)

Plasmid Map:


ACCN: NM_001144887

Insert Size: 582 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq:	<u>NM_001144887.1</u>
RefSeq Size:	830 bp
RefSeq ORF:	582 bp
Locus ID:	4435
UniProt ID:	<u>Q99966</u>
Cytogenetics:	Xq13.1
Protein Families:	Transcription Factors
MW:	19.9 kDa
Gene Summary:	<p>This gene encodes a member of the CREB-binding protein/p300-interacting transactivator with Asp/Glu-rich C-terminal domain (CITED) family of proteins. The encoded protein, also known as melanocyte-specific gene 1, may function as a transcriptional coactivator and may play a role in pigmentation of melanocytes. Alternatively spliced transcript variants have been described. [provided by RefSeq, Jan 2009]</p> <p>Transcript Variant: This variant (4) differs in the 5' UTR compared to variant 1. Both variants 1 and 4 encode the same protein.</p>