

Product datasheet for SC324805

PRR20E (NM_001130407) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PRR20E (NM_001130407) Human Untagged Clone
Tag:	Tag Free
Symbol:	PRR20E
Synonyms:	PRR20; PRR20A; PRR20B; PRR20C; PRR20D
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC324805 representing NM_001130407. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTGTAGTAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
 GATCCGGTACCGAGGAGATCTGCCGCC**CGCATCGCC**
 ATGGAGGAACCAAGGCCTTCGAAGCGACTTCGCTCCATGGCCCCCTAATCAAGCCTCAGGTGGGCCTCCT
 CCAGAGCCAGGCTGCTGTGTTCGGGACCCTGAAGGCTCCGTGGAAGCAGATGGGCCCGCACAGCCAGCC
 CAACCCGCAAAACCCATCGCTTACGTGAAACCTTCAGACGGCAGCCCCAGCTCGCCAGAGTCACCC
 CCTCCTGCAGAGAGAGGCCGGCGCCGGGGAGGAAGCCGGCGGCCAGGGCGAGGCCGTGGCAGAAGGGCT
 GGGCCCCGCGGGACGCTGGCCAGAGACAGGGGGCAGAAGGCTTGATGGCACCGACGTGCACATCCAA
 CTGGACCACCATGGAGAGCCAGGCCACCAGGGGAACCGAAATCACGGAGACCGCAGCCTTCTCCCTT
 TCTGAAACAGGTCTCCGCCTGGAAGTGTGCAGGAAGGCCCTGGCCCCGACGTGGCGCAACCTGAGCTG
 GGGTTTCAGGAGCCGCCGCTGCTCCTGGGCCTCAGGCTGTTGACTGGCAACCCGCTTGACCCCTAT
 CCCTGCATCGGGTTTAGGGCTCTGGGTGACTCAGCGGTTTTACAAGTCATTCAAACCCCCAGGGCACC
 TACGTGCAAGGGTCCCAGTGTTCTCACCACATTGCGTATTGA
 ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
 TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC

Restriction Sites:	SgfI-MluI
ACCN:	NM_001130407
Insert Size:	666 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).


[View online »](#)

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001130407.1
RefSeq Size:	1847 bp
RefSeq ORF:	666 bp
Locus ID:	729250
UniProt ID:	P86496
Cytogenetics:	13q21.1
MW:	23.3 kDa
Gene Summary:	This gene is one of five identical loci in a cluster on chromosome 13q21.1. The predicted protein is proline-rich and contains several dopamine D4 receptor signatures and PRINTS domains. [provided by RefSeq, Oct 2008]