

## Product datasheet for **SC324810**

### **RABL2B (NM\_001130922) Human Untagged Clone**

#### **Product data:**

<b>Product Type:</b>	Expression Plasmids
<b>Product Name:</b>	RABL2B (NM_001130922) Human Untagged Clone
<b>Tag:</b>	Tag Free
<b>Symbol:</b>	RABL2B
<b>Vector:</b>	<u>pCMV6 series</u>
<b>Fully Sequenced ORF:</b>	>NCBI ORF sequence for NM_001130922, the custom clone sequence may differ by one or more nucleotides ATGGCAGAAGACAAAACCAACCGAGTGAGTTGGACCAAGGGAAGTATGATGCTGATGAC AACGTGAAGATCATCTGCCTGGGAGACAGCGCAGTGGGCAAATCCAAACTCATGGAGAGA TTTCTCATGGATGGCTTTCAGCCACAGCAGCTGTCCACGTACGCCCTGACCCTGTACAAG CACACAGCCACGGTAGATGGAAGGACCATCCTTGTGGACTTTTGGGACACGGCAGGCCAG GAGCGGTTCCAGAGCATGCATGCCTCCTACTACCACAAGGCCACGCCTGCATCATGGTG TTTGATGTACAGAGGAAAGTCACCTATAGGAACCTGAGCACCTGGTATACAGAGCTTCGG GAGTTTCAGGCCAGAGATCCCATGCATCGTGGTGGCCAATAAAAATTGATGACATAAACGTG ACCCAAAAAAGCTTCAATTTTGCCAAGAAGTTCTCCCTGCCCTGTATTTTCGTCTCGGCT GCTGATGGTACCAATGTTGTGAAGCTCTTCAATGATGCAATTCGATTAGCTGTGTCTTAC AAACAGAACTCCCAGGACTTCATGGATGAGATTTTTCAGGAGCTCGAGAATTCAGCTTG GAGCAGGAAGAGGAGGACGTGCCAGACCAGGAACAGAGCAGCAGCATCGAGACCCATCA GAGGAGGCGGCCTCTCCCCACAGC
<b>Restriction Sites:</b>	Please inquire
<b>ACCN:</b>	NM_001130922
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** NM\_001130922.1, NP\_001124394.1

**RefSeq Size:** 2187 bp

**RefSeq ORF:** 687 bp

**Locus ID:** 11158

**UniProt ID:** Q9UNT1

**Cytogenetics:** 22q13.33

**Protein Families:** Druggable Genome

**Gene Summary:** The RABL2B protein is a member of the RAB gene family which belongs to the RAS GTPase superfamily. RABL2B is located within a subtelomeric region of 22q13.3. Multiple alternatively spliced transcript variants encoding several different isoforms have been found for this gene. [provided by RefSeq, Aug 2008]

Transcript Variant: This variant (6) has multiple differences in the 5' UTR and coding sequence compared to variant 7. The resulting isoform (2) has the same N- and C-termini but is shorter compared to isoform 3. Variants 2, 6, 10, 11 and 12 all encode the same isoform (2).