

Product datasheet for **SC324898**

CTRP7 (C1QTNF7) (NM_001135171) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	CTRP7 (C1QTNF7) (NM_001135171) Human Untagged Clone
Tag:	Tag Free
Symbol:	C1QTNF7
Synonyms:	CTRP7; ZACRP7
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC324898 representing NM_001135171. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAACGACTACTATAGGGCGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGTTTGTCTTGCTCTATGTTACAAGTTTTGCCATTTGTGCCAGTGGACAACCCCGGGTAAATCAGTTG
AAAGGAGAGAAGTACTCCCCAGGTATATCTGCAGCATTCTGGCTTGCCTGGACCTCCAGGGCCCCCT
GGAGCAAATGGTTCCTCCGTCGGCCCATGGTCGCATCGGCCTTCCAGGAAGAGATGGTAGAGACGGCAGG
AAAGGAGAGAAAGGTGAAAAGGGAAGTGCAGGTTTGGAGAGTAAAGTGGACCGCTAGGTCTTGCCGGT
GAGAAAGGGGACCAAGGAGAGACTGGGAAGAAAGGACCCATAGGACCAGAGGGAGAGAAAGGAGAAGTA
GGTCCAATTGGTCTCTCTGGACCAAGGGAGACAGAGGAGAACAAGGGGACCCGGGGCTGCCTGGAGTT
TGCAGATGTGGAAGCATCGTGCTCAAATCCGCCTTTTCTGTTGGCATCACACCAGCTACCCAGAAGAA
AGACTACCTATTATTTAAACAAGTCTCTTCAACGAGGGAGAGCACTACAACCCTGCCACAGGGAAAG
TTCATCTGTGCTTTCCAGGGATCTATTACTTTTCTTATGATATCACATTGGCTAATAAGCATCTGGCA
ATCGGACTGGTACACAATGGGCAATACCGGATAAAGACCTTCGACGCCAACACAGGAAACCATGATGTG
GCTTCGGGGTCCACAGTCATCTATCTGCAGCCAGAAGATGAAGTCTGGCTGGAGATTTTCTTACAGAC
CAGAATGGCCTCTTCTCAGACCCAGGTTGGGCAGACAGCTTATTCTCCGGTTTCTTATACGTTGAC
ACAGATTACCTAGATTCCATATCAGAAGATGATGAATTGTGA
ACGGGTACGCGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```

Restriction Sites: SgfI-MluI



[View online »](#)

Plasmid Map:


ACCN: NM_001135171

Insert Size: 870 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: NM_001135171.1

RefSeq Size: 4477 bp

RefSeq ORF: 870 bp

Locus ID: 114905
UniProt ID: [Q9BXJ2](#)
Cytogenetics: 4p15.32
Protein Families: Secreted Protein
MW: 30.7 kDa