

Product datasheet for **SC324936**

DHRS9 (NM_001142270) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: DHRS9 (NM_001142270) Human Untagged Clone
Tag: Tag Free
Symbol: DHRS9
Synonyms: 3-alpha-HSD; 3ALPHA-HSD; RDH-E2; RDH-TBE; RDH15; RDHL; RDHTBE; RETSDR8; SDR9C4
Vector: pCMV6 series
Fully Sequenced ORF: >NCBI ORF sequence for NM_001142270, the custom clone sequence may differ by one or more nucleotides

```
ATGCTCTTTTGGGTGCTAGGCCTCCTAATCCTCTGTGGTTTTCTGTGGACTCGTAAAGGA
AACTAAAGATTGAAGACATCACTGATAAGTACATTTTTATCACTGGATGTGACTCGGGC
TTTGAAACTTGGCAGCCAGAACCTTTTGATAAAAAGGGATTTTCATGTAATCGCTGCCTGT
CTGACTGAATCAGGATCAACAGCTTTAAAGGCAGAAACCTCAGAGAGACTTCGACTGTG
CTTCTGGATGTGACCCAGCCAGAGAATGTCAAGAGGACTGCCAGTGGGTGAAGAACC
GTTGGGGAGAAAAGTCTCTGGGGTCTGATCAATAATGCTGGTGTCCCGCGTGTGGCT
CCCACTGACTGGCTGACACTAGAGGACTACAGAGAACCTATTGAAGTGAACCTGTTTGG
CTCATCAGTGTGACACTAAATATGCTTCCTTTGGTCAAGAAAGCTCAAGGAGAGTTATT
AATGTCTCCAGTGTGGAGGTCGCCTTGAATCGTTGGAGGGGGCTATACTCCATCCAAA
TATGCAGTGGAAGGTTTCAATGACAGCTTAAGACGGGACATGAAAGCTTTTGGTGTGCAC
GTCTCATGCATTGAACCAGGATTGTTCAAACAACACTTGGCAGATCCAGTAAAGGTAATT
GAAAAAACTCGCCATTTGGGAGCAGCTGTCTCCAGACATCAAACAACAATATGGAGAA
GGTTACATTGAAAAAGTCTAGACAACTGAAAGGCAATAAATCCTATGTGAACATGGAC
CTCTCTCCGGTGGTAGAGTGCATGGACCAGCTCTAACAAGTCTCTTCCCTAAGACTCAT
TATGCCGCTGAAAAGATGCCAAAATTTCTGGATACCTCTGTCTCACATGCCAGCAGT
TTGCAAGACTTTTTATTGTTGAAACAGAAAGCAGAGCTGGCTAATCCCAAGGCAGTG
```

Restriction Sites: Please inquire
ACCN: NM_001142270
OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.



[View online >](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001142270.1</u> , <u>NP_001135742.1</u>
RefSeq Size:	1760 bp
RefSeq ORF:	960 bp
Locus ID:	10170
UniProt ID:	<u>Q9BPW9</u>
Cytogenetics:	2q31.1
Protein Families:	Druggable Genome
Protein Pathways:	Metabolic pathways, Retinol metabolism
Gene Summary:	<p>This gene encodes a member of the short-chain dehydrogenases/reductases (SDR) family. The encoded protein has been identified as a moonlighting protein based on its ability to perform mechanistically distinct functions. This protein demonstrates oxidoreductase activity toward hydroxysteroids and is able to convert 3-alpha-tetrahydroprogesterone to dihydroxyprogesterone and 3-alpha-androstanediol to dihydroxyprogesterone in the cytoplasm, and may additionally function as a transcriptional repressor in the nucleus. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2014]</p> <p>Transcript Variant: This variant (3) differs in the 5' UTR, compared to variant 2. Variants 2-4 encode the same isoform (1).</p>