

Product datasheet for **SC325023**

TMEM229A (NM_001136002) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TMEM229A (NM_001136002) Human Untagged Clone
Tag:	Tag Free
Symbol:	TMEM229A
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC325023 representing NM_001136002. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCC**CGCATCGCC**
ATGGCGGGGAGCGACGTGGACAGCGAGGGCCCCGACGGAGGGGCGGCGGGCGGGCGTCCGGGGGCC
CCTGGCGGGCCAGGAAGCGAGGCGGCAGCCGGCTGCCCGAGCCGCTGTCCACTGCTGAAGCGCCGGCT
GAGAGCGCCACGCTGCCCGCCTGGATGCGCCTCTACTTCTACGGGATGCACGGGATCACCTGGACGTG
CTGGTGTCTCGGCCCGGCGCTTCGCCCGCAGCCCGGACCTGCGGATGCTAGGCTTCTCTCGCCCTAC
CGCTGCCTGCTGCACTCGCTCACCCATTTGCCCTGGAGAAGGTGTACCTGCAGCAGCGGCGCTGTCCC
AACGCCTTCGTCTTCAATTTCTCTCTACCCCTCGGCCACGTGGGGCTGCAGACCTAGCGGGCCAG
GCGCTACTACTCAGCCTGGCGGCGGGGCGGGGTGCGGTGGCGCCAGGGGCGCTGGACCTGGCGCTG
CAGTACGTGCTGGCGCTTACCACTGCCAAGTGTCTGAAGCGCTTCTGCGCTTGGGTACGGGCGA
CAGAGGCGGCGGCAGCAGCAACAGCAGCAGCAGCAGCAGCAGCAGCAGCGGAGGGGCGCGCTCCCGTC
CCTCCCGGCGCCCGGTCCCTACTGCGGCCGAGCCCGGCGGCGACGCCCGTGGCCCCAGGGGCGCC
GGGGGAGCCCCAGCCAGGGGCTGCCGACCTACCCGCTTTCTTTCTCGGAATGCACGGCTTTCTG
GATGAGATCTTCTCACCTTCTTCTTCAACGTACTGGGGCAGGGGACGGGACAACCAGCGGCCACAG
TCGCTCTGGTCTTCTTTATGTACGCGAGCTGCAGTTTCGTGGTGGAAAAGCTCTACTTCCACCTCCAC
TACAGCCGCGGTTGGGGCACTTGAAGCGGTGCCATCTACGTATCTCATCTACGTGTGGGAGCTG
TCTTGGGGTCTGGGACTCCGCACGTGCGGGGCTTGTTCCTGGGACTATTCTCACTACCCGCTCAATTTT
ATGGGCCCTCATACCCCTGATGTATTTACCTGGCTGGATATTCTTAGTGTGTACCAGGACCTAATTTCC
AACGTGTTGTGGAGGGTGCAGTACGTACCAGCTAA**CTAA**
ACGCGTACGCGGCCGCTCGAGCAGAA**CTCATCTC**AGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC

Restriction Sites: SgfI-MluI

ACCN: NM_001136002



[View online »](#)

Insert Size:	1143 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001136002.1</u>
RefSeq Size:	2569 bp
RefSeq ORF:	1143 bp
Locus ID:	730130
UniProt ID:	<u>B2RXF0</u>
Cytogenetics:	7q31.32
MW:	42.4 kDa