

## Product datasheet for SC325067

### SERPINB2 (NM\_001143818) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	SERPINB2 (NM_001143818) Human Untagged Clone
Tag:	Tag Free
Symbol:	SERPINB2
Synonyms:	HsT1201; PAI; PAI-2; PAI2; PLANH2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC325067 representing NM_001143818. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTAGTGAACCGTCAGAATTTGTAAACGACTACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGAGGATCTTTGTGTGGCAAACACACTCTTTGCCCTCAATTTATTCAAGCATCTGGCAAAGCAAGC
CCCACCCAGAACCTCTTCCTCTCCCAATGGAGCATCTCGTCCACCATGGCCATGGTCTACATGGGCTCC
AGGGGCAGCACCGAAGACCAGATGGCCAAGGTGCTTCAGTTTAAATGAAGTGGGAGCCAATGCAGTTACC
CCCATGACTCCAGAGAATTTACCAGCTGTGGTTCATGCAGCAGATCCAGAAGGGTAGTTATCCTGAT
GCGATTTTGCAGGCACAAGCTGCAGATAAAATCCATTCATCCTCCGCTCTCTCAGCTCTGCAATCAAT
GCATCCACAGGAATTATTTACTGGAAAGTGCAATAAGCTGTTGGTGAGAAGTCTGCGAGCTTCCGG
GAAGAATATATTCGACTCTGTGAGAAATATTACTCCTCAGAACCCAGGCAGTAGACTTCCTAGAATGT
GCAGAAGAAGCTAGAAAAAGATTAATTCCTGGGTCAAGACTCAAACCAAAGGCAAAATCCAAACTTG
TTACCTGAAGGTTCTGTAGATGGGGATACCAGGATGGTCTGGTGAATGTGTCTACTTCAAAGGAAAG
TGGAAACTCCATTTGAGAAGAACTAAATGGGCTTTATCCTTTCCGTGTAAGTCTCGGCTCAGCGCACA
CCTGTACAGATGATGACTTGGGTGAAAAGCTAAACATTTGGATACATAGAAGACCTAAAGGCTCAGATT
CTAGAATCCCATATGCTGGAGATGTTAGCATGTTCTTGTGCTCCAGATGAAATTGCCGATGTGTCC
ACTGGCTTGGAGCTGCTGAAAGTGAATAACCTATGACAACTCAACAAGTGGACCAGCAAAGACAAA
ATGGCTGAAGATGAAGTTGAGGTATACATACCCAGTTCAAATTAAGAAGAGCATTATGAAGTCAAGTCC
ATTCTGAGAAGCATGGGCATGGAGGACGCCTTCAACAAGGGACGGCCAAATTTCTCAGGGATGTCGGAG
AGGAATGACCTGTTTCTTCTGAAGTGTCCACCAAGCCATGGTGGATGTGAATGAGGAGGGCACTGAA
GCAGCCGCTGGCAGAGGAGTGTATGACAGGGAGAAGTGGACATGGAGGCCACAGTTTGTGGCAGAT
CATCCTTTTCTTTTCTTATTATGCATAAGATAACCAACTGCATTTTATTTTTCGGCAGATTTTCTCA
CCCTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```



[View online »](#)

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM\_001143818

Insert Size: 1248 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM\\_001143818.1](#)

RefSeq Size: 2180 bp

RefSeq ORF: 1248 bp

Locus ID: 5055

**UniProt ID:** [P05120](#)

**Cytogenetics:** 18q21.33-q22.1

**Protein Families:** Druggable Genome, Secreted Protein

**MW:** 46.6 kDa

**Gene Summary:** Inhibits urokinase-type plasminogen activator. The monocyte derived PAI-2 is distinct from the endothelial cell-derived PAI-1.[UniProtKB/Swiss-Prot Function]  
Transcript Variant: This variant (1) represents the longer transcript. Both variants 1 and 2 encode the same protein.