

Product datasheet for **SC325068**

UNC93A (NM_001143947) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	UNC93A (NM_001143947) Human Untagged Clone
Tag:	Tag Free
Symbol:	UNC93A
Synonyms:	dj366N23.1; dj366N23.2; Unc-93A
Vector:	<u>pCMV6 series</u>
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001143947, the custom clone sequence may differ by one or more nucleotides

```

ATGGACAGAAAGTCTAAGGAACGTCCTTGTGGTTTCCTTTGGGTTCTGCTTCTCTTTACA
GCCTATGGAGGTCTGCAGAGCCTGCAGAGCAGCCTGTACAGCGAGGAGGCCTGGGTGTC
ACAGCGCTCAGCACCTCTATGGAGGCATGCTCCTGTCTCCATGTTCTCTCCACCGCTC
CTCATCGAGAGGTGGGCTGCAAGGGGACCATCATCTCTCCATGTGTGGCTACGTGGCC
TTCTCCGTGGGCAACTTCTCGCCAGCTGGTACACTTTGATCCCCACCTCCATACTGCTG
GGACTCGGGGCCCGCCGCTGTGGTCTGCACAGTGCACATACCTCACGATCACGGGAAAC
ACACATGCAGAGAAGCGGGAAAGCGTGGCAAAGACATGGTGAACCAAGTATTTTGGCATC
TTCTTCTCATATTCCAGTCATCCGGTGTGTGGGGCAACTTGATCTCATCGCTGGTATTT
GGCCAGACTCCCAGCCAAGGGAGTGGTGTCTGGCTGTCTGATGATAGTGCCTTCCTC
CAACCCATACGAGATGTTTCAGCGGGAAAGTGAAGGAGAGAAGAAATCAGTACCTTTCTGG
TCCACTTTACTGTGCACTTTCAAGCTATATAGAGATAAACGTCTGTGCCTCTTAATTCTG
CTGCCGCTGTACAGTGGATTGCAGCAAGGATTCCTCTCCAGCGAATACACAAGGTCTAT
GTCACCTGCACCTGGGCATCCAGTTCGTGGCTACGTGATGATCTGCTTCTCGGCCACT
GACGCGCTGTGCTCCGTGTTGTATGGAAGGTCTCGCAGTACACGGGCAGGGCTGTGCTG
TACGTGCTGGGCGCGGTGACCCACGTGCTGCTGATGATTGCCCTACTGCTGTGGAGACCT
CGTGCTGACCATCTGGCAGTGTTCTTCGTATTCTCTGGCCTGTGGGGCGTGGCAGATGCC
GTCTGGCAGACACAAAACAATGCTCTCTACGGCGTTCTGTTTGAGAAGAGCAAGGAAGCT
GCCTTCGCCAATTACCGCTGTGGGAGGCCCTGGGCTTCGTATTGCCTTCGGGTACAGC
ATGTTTTTGTGCGTGCACGTCAAGCTCTACATTCTGCTGGGGGTCTGAGCCTGACCATG
GTGGCGTATGGGCTTGTGGAGTGCCTGGAGTCCAAGAACCCGATCAGACCCACGCTCCA
GGACAGGTCAACCAGGCAGAGGATGAAGAAATACAAACAAAATG
  
```

Restriction Sites:	Please inquire
ACCN:	NM_001143947
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).


[View online »](#)

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001143947.1 , NP_001137419.1
RefSeq Size:	2006 bp
RefSeq ORF:	1248 bp
Locus ID:	54346
UniProt ID:	Q86WB7
Cytogenetics:	6q27
Protein Families:	Transmembrane