

Product datasheet for **SC325115**

FIBCD1 (NM_001145106) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	FIBCD1 (NM_001145106) Human Untagged Clone
Tag:	Tag Free
Symbol:	FIBCD1
Vector:	<u>pCMV6 series</u>
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001145106, the custom clone sequence may differ by one or more nucleotides

```

ATGGTCAACGACCGGTGGAAGACCATGGGCGGCGCTGCCCAACTTGAGGACCGGCCGCGC
GACAAGCCGCAGCGGCCGAGCTGCGGCTACGTGCTGTGCACCGTCTGCTGGCCCTGGCT
GTGCTGCTGGCTGTAGCTGTACCGGTGCCGTGCTCTTCTGAACACGCCACGCGCCG
GGCACGGCGCCCCACCTGTCTGTCAGCACTGGGGCTGCCAGCGCCAACAGCGCCCTGGTC
ACTGTGGAAGGGCGGACAGCTCGCACCTCAGCATCCTCATTGACCCGCGCTGCCCGAC
CTCACCGACAGCTTCGCACGCTGGAGAGCGCCAGGCCCTCGGTGCTGCAGGCGTGACA
GAGCACAGGCCACGCCAGGCTGGTGGCGGACCAGGAGCAGGAGCTGCTGGACACGCTG
GCCGACCAGCTGCCCGGCTGCTGGCCGAGCCTCAGAGCTGCAGACGGAGTGATGGGG
CTGCGGAAGGGGCATGGCACGCTGGGCCAGGGCCTCAGCGCCCTCAGAGTGAGCAGGGC
CGCCTCATCCAGTTCTCTGAGAGCCAGGGCCACATGGCTCACCTGGTGAACCTCCGTC
AGCGACATCCTGGATGCCCTGCAGAGGGACCGGGGGCTGGGCCGCCCCGCAACAAGGCC
GACCTTCAGAGAGCGCCTGCCCGGGGAACCCGGCCCCGGGGCTGTCCACTGGCTCCCGG
CCCCGAGACTGTCTGGACGTCCTCCTAAGCGGACAGCAGGACGATGGCGTCTACTCTGTC
TTTCCCACCCACTACCCGGCCGGCTTCCAGGTGTACTGTGACATGCGCACGGACGGCGGC
GGCTGGACGGTGTTCAGCGCCGGGAGGACGGCTCCGTGAACCTTCCGGGGCTGGGAT
GGTACCAGAGACGGCTTTGGCAGGCTCACCGGGGAGCACTGGCTAGGGCTCAAGAGGATC
CACGCCCTGACCACAGGCTGCCTACGAGCTGCACGTGGACCTGGAGGACTTTGAGAAT
GGCACGGCCTATGCCCGCTACGGGAGCTTCGGCGTGGGCTTGTCTCCGTGGACCCCTGAG
GAAGACGGGTACCCGCTCACCGTGGCTGACTATTCCGGCACTGCAGGCGACTCCCTCCTG
AAGCACAGCGCATGAGGTTACCACCAAGGACCGTGACAGCGACCATTCAGAGAACAAC
TGTGCCGCTTCTACCGCGTGCCTGGTGGTACCGCAACTGCCACACGTCCAACCTCAAT
GGGCAGTACCTGCGCGGTGCGCACGCTCCTATGCCGACGGCGTGAGTGGTCTCCTCGG
ACCGGCTGGCAGTACTCACTCAAGTTCTCTGAGATGAAGATCCGGCCGGTCCGGGAGGAC
CGC

```

Restriction Sites:	Please inquire
ACCN:	NM_001145106



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001145106.1 , NP_001138578.1
RefSeq Size:	3119 bp
RefSeq ORF:	1386 bp
Locus ID:	84929
UniProt ID:	Q8N539
Cytogenetics:	9q34.12
Protein Families:	Druggable Genome, Transmembrane
Gene Summary:	<p>FIBCD1 is a conserved type II transmembrane endocytic receptor that binds chitin and is located primarily in the intestinal brush border (Schlosser et al., 2009 [PubMed 19710473]). [supplied by OMIM, Apr 2010]</p> <p>Transcript Variant: This variant (2) differs in the 5' UTR compared to variant 1. Both variants 1 and 2 encode the same protein.</p>