

Product datasheet for **SC325116**

Retinoic Acid Receptor alpha (RARA) (NM_001145301) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Retinoic Acid Receptor alpha (RARA) (NM_001145301) Human Untagged Clone
Tag:	Tag Free
Symbol:	RARA
Synonyms:	NR1B1; RAR
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC325116 representing NM_001145301.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCCAGCAACAGCAGCTCTGCCGACACCTGGGGCGGGCACCTCAATGGGTACCCGGTGCCTCCC
TAGCCCTTCTTCTCCCCCTATGCTGGGTGACTCTCCCGCCAGGCGCTCTGACCACTCTCCAGCAC
CAGCTTCCAGTTAGTGGATATAGCACACCATCCCAGCCACCATTGAGACCCAGAGCAGAGTTCTGAA
GAGATAGTGCCAGCCCTCCCTCGCCACCCCTCTACCCCGCATCTACAAGCCTTGCTTTGTCTGTGAG
GACAAGTCTCAGGCTACCACTATGGGGTCAAGCCTGTGAGGGCTGCAAGGGCTTCTCCGCCGACGC
ATCCAGAAGAACATGGTGTACACGTGTACCCGGGACAAGAAGTGCATCATCAACAAGGTGACCCGGAAC
CGTGCCAGTACTGCCGACTGCAGAAGTCTTTGAAGTGGGCATGTCCAAGGAGTCTGTGAGAAACGAC
CGAAACAAGAAGAAGAAGGAGGTGCCAAGCCCGAGTGTCTGAGAGCTACACGCTGACGCCGGAGGTG
GGGGAGCTCATTGAGAAGTGCCAAAGCGCACAGGAAACCTTCCCTGCCCTCTGCCAGCTGGGCAAA
TACACTACGAACAACAGCTCAGAACAACGTGTCTCTGGACATTGACCTCTGGGACAAGTTCAGTGAA
CTCTCCACCAAGTGCATCATTAAAGACTGTGGAGTTCGCCAAGCAGCTGCCCGGCTTACCACCCTCAC
ATCGCCGACCAGATCACCTCCTCAAGGCTGCCTGCCTGGACATCCTGATCCTGCGGATCTGCACGCGG
TACACGCCGAGCAGGACACCATGACCTTCTCGGACGGGCTGACCTGAACCGGACCCAGATGCACAAC
GCTGGCTTCGGCCCCCTCACCGACTGGTCTTTGCCTTCGCCAACCAGCTGTGCCCTGGAGATGGAT
GATGCGGAGACGGGGTCTCAGCGCCATCTGCCTCATCTGCGGAGACGCCAGGACCTGGAGCAGCCG
GACCGGTGGACATGCTGCAAGGACCGTGTGGAGGCGTAAGGTCTACGTGCGGAAGCGGAGGCC
AGCCGCCCCACATGTTCCCAAGATGCTAATGAAGATTACTGACCTGCAAGCATCAGCGCAAGGGG
GCTGAGCGGGTGATCACGCTGAAGATGGAGATCCCGGGCTCCATGCCGCTCTCATCCAGAAATGTTG
GAGAACTCAGAGGGCCTGGACACTCTGAGCGGACAGCCGGGGGTGGGGGGCGGGACGGGGTGGCCCTG
GCCCCCCGCGCAGGCAGCTGTAGCCCCAGCCTAGCCCCAGCTCCAACAGAAGCAGCCCGGCCACCCAC
TCCCCGTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites:

Sgfl-MluI

Plasmid Map:



ACCN: NM_001145301

Insert Size: 1389 bp

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001145301.2</u>
RefSeq Size:	3413 bp
RefSeq ORF:	1389 bp
Locus ID:	5914
UniProt ID:	<u>P10276</u>
Cytogenetics:	17q21.2
Protein Families:	Druggable Genome, Nuclear Hormone Receptor, Transcription Factors
Protein Pathways:	Acute myeloid leukemia, Pathways in cancer
MW:	50.8 kDa
Gene Summary:	<p>This gene represents a nuclear retinoic acid receptor. The encoded protein, retinoic acid receptor alpha, regulates transcription in a ligand-dependent manner. This gene has been implicated in regulation of development, differentiation, apoptosis, granulopoiesis, and transcription of clock genes. Translocations between this locus and several other loci have been associated with acute promyelocytic leukemia. Alternatively spliced transcript variants have been found for this locus.[provided by RefSeq, Sep 2010]</p> <p>Transcript Variant: This variant (3) reflects the use of an alternate promoter and contains a different 5' UTR segment, compared to variant 1. Variants 1 and 3 encode the same isoform (1).</p>