

Product datasheet for **SC325135**

TCHP (NM_001143852) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TCHP (NM_001143852) Human Untagged Clone
Tag:	Tag Free
Symbol:	TCHP
Synonyms:	TpMs
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC325135 representing NM_001143852.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCGCTCCCGACGCTGCCGTCTACTGGTGCAGCCAGCAGCGCCTGAATCAGCAGCTAGCACGACAG
CGAGAGCAGGAGGCCCGCTTCGGCAGCAGTGGGAGCAGAACAGCCGTTACTTCAGGATGTCTGACATC
TGCAGCTCCAAACAGGCAGAATGGAGCTCTAAACCTCCTACCAGCGGAGCATGCATGCCTATCAGCGG
GAGAAGATGAAGGAGGAGAAGAGGAGGAGTCTGGAGGCCGACGGGAAAAGCTCAGGCAGCTCATGCAG
GAGGAGCAGGACCTGCTGGCCAGAGAACTGGAGGAGCTGAGGCTGAGCATGAACTTCAGGAAAAGA
ATCCGGGAGCAGCACGGGAAGCTGAAATCAGCCAAAGAAGAGCAGAGGAACTGATTGCTGAACAACTT
TTGTACGAACTTGAAAAAGAACAACCCGAACTTCGAGAGATGGAGCTGGACCTTCACCAGAAGCAT
GTCGTAACCTCTGGGAAATGCAGAAAGAAGAAAAAACAGCAAGAAGCCACCGCAGAGCAAGAGAAC
AAACGGTATGAAATGAATATGAAAGGGCCGAAGGGAGGCGCTAGAAAGGATGAAAGCTGAAGAGGAG
AGGAGGCAGCTGGAGGACAAGCTCCAGGCCGAGGCACTGCTGCAACAGATGGAGGAGCTGAAGCTGAAG
GAGGTGGAGGCGACAACTAAAGAAGGAGCAGGAGAATCTGTTGAAGCAGCGGTGGGAGCTAGAGAGG
CTGGAGGAAGAGCGAAAGCAGATGGAAGCCTCCGGCAGAAAGCAGAGCTGGGGCGTTTCTTGAGACAT
CAGTATAACGCTCAACTCAGCAGACGCACACAGCAGATCCAAGAGGAGCTGGAGGCAGACAGGCGGATC
CTGCAGGCCCTCCTCGAGAAGGAGGACGAGAGCCAGCGCTCCACCTGGCCAGGCGGGAGCAGGTTCATG
GCCGATGTGGCCTGGATGAAGCAGGCCATTGAGGAGCAGCTGCAGCTGGAGCGGGCGGGAGGCAGAG
CTGCAGATGCTGCTGAGGAGGAGGCCAAGGAGATGTGGGAAAAGAGAGAGGCAGAGTGGGCCGAGAG
CGCAGCGCACGGGACAGACTGATGAGCGAGTTCTGACAGGGAGACAACAGCAAAATACAAGAGAAGATT
GAGCAGAACCACGGGCACAAGAGGAATCCCTGAAACACAGGGAGCAACTTATTCGAAATCTTGAGGAG
GTGAGAGAGTTGGCTCGTCGCGAGAAAGAGGAGAGTGAAGAGCTGAAATCGGCCAGGAAGCAGGAGCTG
GAAGCCCAGTTGCAGAGCGCCGCTGCAGGCATGGGAAGCAGACCAGCAGGAGGAGGAGGAAGAGGAG
GAGGCCCGGGCGGTGAGCAGCTCTCAGATGCCCTGCTGCAGCAGGAGGCGGAGACTATGGCTGAGCAG
GGTACCGGCCTAAGCCTTACGGACATCCAAAAATTGCTTGAACTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites:

SgfI-MluI

Plasmid Map:



ACCN:

NM_001143852

Insert Size:	1497 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001143852.1
RefSeq Size:	3111 bp
RefSeq ORF:	1497 bp
Locus ID:	84260
UniProt ID:	Q9BT92
Cytogenetics:	12q24.11
MW:	61.1 kDa
Gene Summary:	<p>Tumor suppressor which has the ability to inhibit cell growth and be pro-apoptotic during cell stress. Inhibits cell growth in bladder and prostate cancer cells by a down-regulation of HSPB1 by inhibiting its phosphorylation. May act as a 'capping' or 'branching' protein for keratin filaments in the cell periphery. May regulate K8/K18 filament and desmosome organization mainly at the apical or peripheral regions of simple epithelial cells (PubMed:15731013, PubMed:18931701). Is a negative regulator of ciliogenesis (PubMed:25270598).[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (2) differs in the 5' UTR compared to variant 1. Variants 1 and 2 both encode the same protein.</p>