

## Product datasheet for **SC325146**

### **VRK2 (NM\_001130481) Human Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	VRK2 (NM_001130481) Human Untagged Clone
Tag:	Tag Free
Symbol:	VRK2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >SC325146 representing NM\_001130481.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGCCACCAAAAAGAAATGAAAAATACAAACTTCTATCCATTCCAGAAGCAAGTTCTGGATGAT
ATGGAAGCAATCAGTGGTACTGGCAAGAAGATTGGCTCTGGAGATTGGATTGATATATTTAGCT
TTCCCCACAAATAAACCCAGAGAAAAGATGCAAGACATGTAGTAAAAGTGGAAATATCAAGAAAATGGCCCG
TTATTTTCAGAACTTAAATTTTATCAGAGAGTTGCAAAAAAGACTGTATCAAAAAGTGGATAGAACGC
AAACAACCTGATTATTTAGGAATTCCTCTGTTTTATGGATCTGGTCTGACTGAATCAAGGGAAGAAGT
TACAGATTTATGGTAATGGAAAGACTAGGAATAGATTTACAGAAGATCTCAGGCCAGAATGGTACCTTT
AAAAAGTCAACTGTCCTGCAATTAGGTATCCGAATGTTGGATGTACTGGAATATATACATGAAAATGAA
TATGTTTATGGTATATAAAGCAGCAATCTACTTTGGGTTACAAAAATCCAGACCAGGTTTATCTT
GCAGATTATGGACTTTCCTACAGATATTGTCCCAATGGGAACCAAAACAGTATCAGGAAAATCCTAGA
AAAGGCCATAATGGGACAATAGAGTTTACCAGCTTGGATGCCACAAGGGAGTAGCCTTGTCCAGACGA
AGTGACGTTGAGATCCTCGGCTACTGCATGCTGCGGTGGTTGTGTGGGAAACTTCCCTGGGAACAGAAC
CTGAAGGACCCTGTGGCTGTGCAGACTGCTAAAACAAATCTGTTGGACGAGCTCCCCCAGTCAGTGCTT
AAATGGGCTCCTTCTGGAAGCAGTTGCTGTGAAATAGCCCAATTTTGGTATGTGCTCATAGTTTAGCA
TATGATGAAAAGCCAAACTATCAAGCCCTCAAGAAAATTTGAACCTCATGGAATACCTTTAGGACCA
CTGGACTTTTCCACAAAAGGACAGAGTATAAATGTCCATACTCCAACAGTCAAAAAGTTGATTCACAA
AAGGCTGCAACAAAGCAAGTCAACAAGGCACACAATAGGTTAATCGAAAAAAGTCCACAGTGAGAGA
AGCGCTGAGTCTGTGCAACATGGAAAGTGCAGAAAGAGGAGAACTGATTGGATTGATGAACAATGAA
GCAGCTCAGGAAAGCACAAGGAGAAGCAGAAAATATCAAGAGTCTCAAGAACCTTGAATGAAGTAAAC
AGTTTCCACAAAAAATCAGCTATACACAATTCCAAACCTCATTTTATGACCTCATCAAGATTTTACC
AGTCCAGATATATTCAAGAAGTCAAGATCTCCATCTTGGTATAAATACACTTCCACAGTCAAGCAGGGG
ATCACAGACTTAGAAAAGTTCAACTGGACTTTGGCCTACAATTTCCAGTTTACTCTTAGTGAAGAGACA
AACGCAGATGTTTATTATTCGCATCATACACTGCTCTTTTATGATGTTAGTATTTCTGCTTTATTT
TTTCTCTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:**

SgfI-MluI

**Plasmid Map:**



**ACCN:** NM\_001130481

<b>Insert Size:</b>	1527 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001130481.2</a>
<b>RefSeq Size:</b>	1781 bp
<b>RefSeq ORF:</b>	1527 bp
<b>Locus ID:</b>	7444
<b>UniProt ID:</b>	<a href="#">Q86Y07</a>
<b>Cytogenetics:</b>	2p16.1
<b>Protein Families:</b>	Druggable Genome, Protein Kinase, Transmembrane
<b>MW:</b>	58.1 kDa
<b>Gene Summary:</b>	<p>This gene encodes a member of the vaccinia-related kinase (VRK) family of serine/threonine protein kinases. The encoded protein acts as an effector of signaling pathways that regulate apoptosis and tumor cell growth. Variants in this gene have been associated with schizophrenia. Alternative splicing results in multiple transcript variants that differ in their subcellular localization and biological activity. [provided by RefSeq, Jan 2014]</p> <p>Transcript Variant: This variant (3) differs in the 5' UTR compared to variant 1. Variants 1, 2, 3, and 8 encode the same isoform (1, also known as VRK2A).</p>