

Product datasheet for **SC325150**

DELE1 (NM_001142603) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	DELE1 (NM_001142603) Human Untagged Clone
Tag:	Tag Free
Symbol:	DELE1
Synonyms:	DELE; DELE1(L); KIAA0141
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC325150 representing NM_001142603.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTGGCGCCTCCCGGGACTCTGGGCCGAGCTTTCCCGTACACTGGGACCTAGCCTCTGGAGGGTG
ACTCCTAAGTCCACCAGCCAGATGGCCCTCAGACTACCTCCTCACTTTGCTGGTTCTGTGCCTAAC
CTCGACAGGTCAGGTCCCATGGCCAGGCACGAGCGGGGTCCAAGTCCCATGGATGGAAGGATGCC
TTCCAATGGATGTCTTCCCGTGTCTCCCGAACCCCTATGGGATGCCATATCTTGGGGACTCTGGCC
GTGCTGGCCCTGCAGCTGGCAAGGCAGATCCACTTCCAGGCATCCCTGCCAGCAGGACCTCAGCGGGTA
GAACACTGCTCCTGGCACAGTCCCTGGACCGTTTCTTCTCATCTCCCTTGTGGCACCCATGCTCCTCA
CTGCGACAACACATCCTCCCAGCCCCGATGGCCAGCTCCAGGCACACTGGCCTCAGGGAACCCAGG
CTTGGCCAGGAAGAAGCCTCAGCTCAGCCCCGAACTTCTCACAACTCTTTGAGAGGAGCTCGTCTC
CAGGACCCCTCTGAGGAAGTCCCGGTGATTTTGGCTTCTGCATGCCAGTAGTAGCATCGAGTCCGAG
GCAAACAGCCAGCCCTCAGCCACTGGTAAAAGGAACAAGATAAATCAAAAACCTTTCCCTTGAG
GAGGCTGTGACTTCCATTCAGCAGCTTCCAGCTCAGTGTTCATCGCTTTCAACTTCTGGGAACA
GAGAACATGAAGAGTGGCGACCACAGGCAGCCTTTTCTTACTTCCAGAAAGCTGCAGCCCGCGGTAC
AGCAAAGCGCAGTACAATGCGGGCTTGTGTATGAGCATGGCAGAGGCACCCCCAGGGACATTAGCAAG
GCGGTCTTTATTATCAGTTGGCTGCCAGCCAGGGCCACAGCCTGGCTCAGTACCCTATGCCAGGTGC
CTACTACGAGACCCAGCCTTTCGTGGAACCTGAGCGGCAGAGGGCAGTGTCTTGTGAAGCAGGCT
GCAGACTCAGGCTTGAGAGAGGCCAAGCTTTCCTCGGGTGCTTTTACCAAGGAGCCCTACCTGGAT
GAGCAGAGAGCTGTGAAATATCTTTGGCTTGCAGCAACAATGGGACTCACAGAGCAGGTACCACCTT
GGATTTTGTATGAGAAAGGCCTTGGTGTGCAGAGGAATCTGGGAGAGGCCTTGAGATGTTACCAGCAG
TCAGCCGCTCTGGAAATGAGGCCGCCAGGAGAGGCTGCGAGCCCTTTTTCCATGGGGGCTGCAGCC
CCGGGGCCAGCGACCTGACAGTTACAGACTGAAGTCTTCTCCAGCCCTCCCTCTGCAGCTTGAAAC
ACCCTGTAGCAGGAACCTCAGCCTACCACATGCCTCGAGCACAGGCAACCTTGGCCTCCTCTGCAGA
AGTGGGCATCTCGGAGCCAGCCTGGAAGCCTCCAGCAGGGCTATCCCCACACCCTACCCACTGGAA
AGGAGTGTGTAAGACTAGGTTTGGCTAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites:

SgfI-MluI

Plasmid Map:



ACCN:

NM_001142603

Insert Size:	1548 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001142603.2
RefSeq Size:	4195 bp
RefSeq ORF:	1548 bp
Locus ID:	9812
UniProt ID:	Q14154
Cytogenetics:	5q31.3
MW:	55.9 kDa
Gene Summary:	<p>Essential for the induction of death receptor-mediated apoptosis through the regulation of caspase activation.[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (2) is alternatively spliced at the 3' end and is missing a segment of the 3' UTR compared to variant 1. Transcript variants 1 and 2 encode the same protein. Sequence Note: This RefSeq record was created from transcript and genomic sequence data because no single transcript was available for the full length of the gene. The extent of this transcript is supported by transcript alignments.</p>