

Product datasheet for **SC325195**

TBCCD1 (NM_001134415) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TBCCD1 (NM_001134415) Human Untagged Clone
Tag:	Tag Free
Symbol:	TBCCD1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC325195 representing NM_001134415.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGATCAGTCCAGAGTTCTCCTCTGGGTGAAAGCAGAACCCTTTATAGTGGGTGCCTTGCAGGTCCCC
CCTCCATCCAAGTTTAGTCTTCACTATCTCAGGAAGATATCCACCTATGTGCAAATCCGGGCCACAGAA
GGAGCTTACCGCGCCTCTACTGGTCTACATGGAGGCACATCGCTTGTGGGAAGCTGCAGTTGGCCAAG
GACCTGGCGTGGCTTTACTTCGAAATATTTGATAGTCTTTCAATGAAGACACCTGAGGAGCGCCTGGAA
TGGTCTGAGGTTCTGTCCAACCTGCATGTCTGAGGAGGAAGTTGAAAAGCAGAGAAATCAGCTTTCAGTG
GACACGCTACAGTTTCTGCTCTTCTTATACATTCAACAGTTGAACAAGGTCTCCCTAAGGACATCTTTG
ATTGGCGAAGAGTGGCCAGTCCGAAAACAAATCTCAGTCTCCTGACCTGACTGAAAAATCTAATTGT
CATAATAAGAAGTGAATGATTACAGTCACCAAGCTTTTGTCTATGATCATCTGTCTGATCTCCTCGAG
CTGCTTTTAGATCCAAAACAACACTGATCATTTTCATTCAACCCATAGTAGTCTAGTGTCTCGAGAA
GCTGTTGTGGCGCTCAGCTTCTTATTGAAGGTACAATAAGTAGAGCCAGGAAGTCTATCCACTTCAT
GAACTTGCAGTGTGGCAACCACTGCATGCAGATAGTGGCTTCTCAAAGATCTCTAAGACTTCTCTTTC
TACAAACTGGAACCTGGTTGAGGTCTGTTTACTGGGAAATCCATTTGGTACATCAGCTTGCCCTCAAG
TCTGGAAGAAATGGCTTGGGCTCATCAAGTTGAAGGGACCACAAAAGAGCTAAGATTGCTTGTAAAT
ACTCATGTGGCCCTAGGATGCACCGACTGGTAGTGATGAGCCAGGTTTACAAGCAGACACTGGCTAAG
AGCTCAGACACTCTGGCGGGGCACATGTAAGATTATCGTTGCAACGAATCTTTTATATATCTGCTC
TCTCCCTACGATCTGTGACAATTGAGAAGTGCAGGAATAGCATCTTTGTCTTGGGCCCTGTAGGGACT
ACACTTCACTCCACAGTTGTGACAATGTTAAAGTCATTGCTGTTGCCATCGTTTGTCCATCTCTTCT
ACAACAGGTTGCATCTTTCACGTTCTTACGCCTACACGCCCACTTATTCTCTGGGAACCCAGACAGTA
ACTTTTGGCCCTTTTATACACATTACCAATGCTAGAGGACCATATGGCCAGGACTGGCCTTGCTACA
GTGCCAACTATTGGGATAATCCAATGGTTGTGTGCAGAGAGAACAGCGACACAAGAGTCTTCCAGCTT
TTACCACCTTGTGAATCTATGTATTTATTATCCCTTTGAAATGGAAGGGGACACAACAGAGATACCC
GGGGTCTTCCATCTGTATATCAGAAAGCACTGGGTCAAAGAGAACAGAAGATACAGATCTGGCAGAAA
ACTGTGAAGGAGGCTCATTTGACAAAGGATCAAAGGAAGCAGTTCAGGTACTGGTAGAGAAACAAGTTT
TATGAATGGTTGATTAATACAGGACATCGCCAACAGCTGGACAGCCTGTACCCCTGCAGCAGGCTCC
AAACAAGCAGCTGGATAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites:

SgfI-MluI

Plasmid Map:



ACCN:	NM_001134415
Insert Size:	1674 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001134415.1</u>
RefSeq Size:	3080 bp
RefSeq ORF:	1674 bp
Locus ID:	55171
UniProt ID:	<u>Q9NVR7</u>
Cytogenetics:	3q27.3
MW:	63.5 kDa
Gene Summary:	<p>Plays a role in the regulation of centrosome and Golgi apparatus positioning, with consequences on cell shape and cell migration.[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (1) represents the longest transcript and encodes the longer isoform (a). Both variants 1 and 2 encode the same isoform. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>