

Product datasheet for **SC325199**

NT5C2 (NM_001134373) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	NT5C2 (NM_001134373) Human Untagged Clone
Tag:	Tag Free
Symbol:	NT5C2
Synonyms:	cN-II; GMP; NT5B; PNT5; SPG45; SPG65
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC325199 representing NM_001134373.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGTCAACCTCCTGGAGTGATCGGTTACAGAATGCAGCAGATATGCCTGCTAACATGGATAAGCATGCC
CTGAAAAAGTATCGTCGAGAAGCCTATCATCGGGTGTGTGTAACCGAAGTTTAGCAATGGAAAAGATA
AAGTGTTTTGGTTTTGATATGGATTATACCCTTGCTGTGTACAAGTCCCCAGAGTATGAGTCCCTGGT
TTTGAGCTTACTGTGGAGAGATTAGTTTCTATTGGCTATCCCCAGGAGTTGCTCAGCTTTGCTTATGAT
TCTACATTCCCTACCAGGGGACTTGTCTTTGACACACTGTATGGAAATCTTTGAAAGTCGATGCCTAT
GGAAACCTCTTGGTCTGTGCACATGGATTTAACTTTATAAGGGGACCAGAACTAGAGAACAGTATCCA
AATAAATTTATCCAGCGAGATGATACTGAAAGATTTTACATTCTGAACACACTATTCAACCTACCAGAG
ACCTACCTGTTGGCTGCCTAGTAGATTTTTTACTAATTGTCCAGATATACCAGTTGTGAAACAGGA
TTTAAAGATGGGACCTTTCATGTCCTACCGGAGTATGTTCCAGGATGTAAGAGATGCTGTTGACTGG
GTTCAATACAAGGCTCCCTAAGGAAAAGACAGTTGAAAATCTTGAGAAGTATGTAGTCAAAGATGGA
AAACTGCCTTTGCTTCTGAGCCGATGAAGGAAGTAGGAAAAGTATTTCTTGCTACCAACAGTGACTAT
AAATATACAGATAAAATATGACTTACCTGTTGACTTCCCACATGGCCCAAGCCTGGGAGCTCCCAT
CGACCATGGCAGTCTACTTTGACTTGATCTTGGTGGATGCACGGAAACCACTCTTTTTGGAGAAGGC
ACAGTACTGCGTCAGGTGGATACTAAAAGTGGCAAGCTGAAAATGGTACCTACACAGGGCCCTACAG
CATGGTATCGTCTACTCAGGAGTTCTTCTGATACGATCTGTGACCTGTTGGGAGCCAAGGGAAAAGAC
ATTTTGTATATTGGAGATCACATTTTGGGGACATTTTAAATCAAAGAAACGGCAAGGGTGGCGAACT
TTTTGGTGATTCTGAACTCGCACAGGAGCTACATGTCTGGACTGACAAGAGTTCATTTTCGAAGAA
CTTCAGAGCTTGATATTTTCTGGCTGAACTCTACAAGCATCTTGACAGCAGTAGCAATGAGCGTCCA
GACATCAGTTCCATCCAGAGACGTATTAAGAAAAGTAACTCATGACATGGACATGTGCTATGGGATGATG
GGAAGCCTGTTTCGACAGTGGCTCCCGGCAGACCCTTTTGGCAGTCAAGTATGCGTTATGCTGACCTC
TATGCAGCATCTTTCATCAACCTGCTGTATTACCCTTTCAGTACCTCTTCAGGGCTGCCCATGTCTTG
ATGCCTCATGAATCAACGGTGGAGCACACACAGTAGATATCAATGAGATGGAGTCTCCTCTTGCCACC
CGGAACCGCACATCAGTGGATTTCAAAGACACTGACTACAAGCGGCACCAGCTGACACGGTCAATTAGT
GAGATTAACCTCCCAACCTCTCCCACTGGCCCCCAGGAAATTACACACTGCCATGACGAAGATGAT
GATGAAGAGGAGGAGGAGGAGGAAGAAATAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001134373

Insert Size: 1686 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001134373.2](#)

RefSeq Size: 3453 bp

RefSeq ORF: 1686 bp

Locus ID: 22978

UniProt ID: [P49902](#)

Cytogenetics: 10q24.32-q24.33

Protein Pathways:	Metabolic pathways, Nicotinate and nicotinamide metabolism, Purine metabolism, Pyrimidine metabolism
MW:	65 kDa
Gene Summary:	<p>This gene encodes a hydrolase that serves as an important role in cellular purine metabolism by acting primarily on inosine 5'-monophosphate and other purine nucleotides. [provided by RefSeq, Oct 2011]</p> <p>Transcript Variant: This variant (2) differs in the 5' UTR, compared to variant 1. Variants 1-3 encode the same protein (isoform 1).</p>