

Product datasheet for **SC325203**

CDC7 (NM_001134419) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CDC7 (NM_001134419) Human Untagged Clone
Tag:	Tag Free
Symbol:	CDC7
Synonyms:	CDC7L1; HsCDC7; Hsk1; huCDC7
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC325203 representing NM_001134419.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGAGGCGTCTTTGGGGATTCCAGATGGATGAGCCAATGGCTTTTCTCCCAGCGTGACCGGTTTCAG
GCTGAAGGCTCTTTAAAAAAAACGAGCAGAATTTAACTTGAGGTGTTAAAAAGATATTGAGAAG
CTTTATGAAGCTGTACCACAGCTTAGTAATGTGTTAAGATTGAGGACAAAATTGGAGAAGGCATTTTC
AGCTCTGTTTATTTGGCCACAGCACAGTTACAAGTAGGACCTGAAGAGAAAATTGCTCTAAAACACTTG
ATTCCAACAAGTCATCCTATAAGAATTGCAGCTGAACTTCAGTGCCTAACAGTGGCTGGGGGCAAGAT
AATGTCATGGGAGTTAAATACTGCTTTAGGAAGAATGATCATGTAGTTATTGCTATGCCATATCTGGAG
CATGAGTCGTTTTTGGACATTCTGAATTCTTTTCCTTTCAAGAAGTACGGGAATATATGCTTAATCTG
TTCAAAGCTTTGAAACGCATTTCAGTTTGGTATTGTTACCGTGATGTTAAGCCAGCAATTTTTTA
TATAATAGGCGCTGAAAAAGTATGCCTTGGTAGACTTTGGTTTGGCCCAAGGAACCCATGATACGAAA
ATAGAGCTTCTTAAATTTGTCCAGTCTGAAGCTCAGCAGGAAAGGTGTTCAAAAACAAATCCACATA
ATCACAGGAAACAAGATTCCACTGAGTGGCCAGTACCTAAGGAGCTGGATCAGCAGTCCACCACAAAA
GCTTCTGTTAAAAAGCCCTACACAAATGCACAAATTCAGATTAACAAGGAAAAGACGGAAAGGAGGGA
TCTGTAGGCTTTCTGTCCAGCGCTCTGTTTTTGGAGAAAAGAAATTTCAATATACACAGCTCCATTTCA
CATGAGAGCCCTGCAGTAAAACATGAAGCAGTCAAAGACTGTGGATGTACTGTCTAGAAAAGTTAGCA
ACAAAAAGAAGGCTATTTCTACAAAAGTTATGAATAGTGTGTGATGAGGAAAAGTCCAGTCTTGTGC
CCAGCTAGCCTGACCTGTGACTGCTATGCAACAGATAAAGTTTGTAGTATTTGCCTTTCAAGGCGTCAG
CAGGTTGCCCTAGGGCAGGTACACCAGGATTCAGAGCACCAGAGGTCTTGACAAAAGTCCCAATCAA
ACTACAGCAATTGACATGTGGTCTGCAGTGTGCATATTTCTTTCTTTGCTTAGTGGACGATATCCATTT
TATAAAGCAAGTGTGATTTAACTGCTTTGGCCCAAAATATGACAATTAGGGGATCCAGAGAAACTATC
CAAGCTGCTAAAACCTTTTGGGAAATCAATATTATGTAGCAAAGAAGTTCCAGCACAAAGCTTGAGAAAA
CTCTGTGAGAGACTCAGGGGTATGGATTCTAGCACTCCCAAGTTAACAAGTATATAAAGGCGATGCT
TCTCATCAACCAGCTATTTAGAGAAGACTGACCATAAAGCTTCTTGCCTCGTTCAAACCTCCAGGA
CAATACTCAGGGAATTCATTTAAAAGGGGGATAGTAATAGCTGTGAGCATTGTTTTGATGAGTATAAT
ACCAATTTAGAAGGCTGGAATGAGGTACCTGATGAAGCTTATGACCTGCTTGATAAACTTCTAGATCTA
AATCCAGCTTCAAGAATAACAGCAGAAGAAGCTTTGTTGCATCCATTTTTAAAGATATGAGCTTGTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001134419

Insert Size: 1725 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001134419.1](#)

RefSeq Size: 3188 bp

RefSeq ORF: 1725 bp

Locus ID: 8317

UniProt ID: [O00311](#)

Cytogenetics: 1p22.1

Protein Families:	Druggable Genome, Protein Kinase, Stem cell - Pluripotency
Protein Pathways:	Cell cycle
MW:	63.9 kDa
Gene Summary:	<p>This gene encodes a cell division cycle protein with kinase activity that is critical for the G1/S transition. The yeast homolog is also essential for initiation of DNA replication as cell division occurs. Overexpression of this gene product may be associated with neoplastic transformation for some tumors. Multiple alternatively spliced transcript variants that encode the same protein have been detected. [provided by RefSeq, Aug 2008]</p> <p>Transcript Variant: This variant (2) uses a different splice site in the 5' UTR, compared to variant 1. Variants 1, 2, and 3 all encode the same protein.</p>