

Product datasheet for **SC325214**

GAD65 (GAD2) (NM_001134366) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	GAD65 (GAD2) (NM_001134366) Human Untagged Clone
Tag:	Tag Free
Symbol:	GAD65
Synonyms:	GAD65
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >NCBI ORF sequence for NM_001134366, the custom clone sequence may differ by one or more nucleotides

```

ATGGCATCTCCGGCTCTGGCTTTTGGTCTTTCGGGTCGGAAGATGGCTCTGGGGATTCCGAGAATCCCG
GCACAGCGCGAGCCTGGTGCCTGCAAGTGGCTCAGAAGTTCACGGGCGGCATCGGAAACAACTGTGCGCCCT
GCTCTACGGAGACGCGGAGAAGCCGGCGGAGAGCGGGAGCCAAACCCCGGGCCGCCCGCCGGAAG
GCCGCCTGCGCCTGCGACCAGAAGCCCTGCAGCTGCTCCAAAGTGGATGTCAACTACGCGTTTCTCCATG
CAACAGACCTGCTGCCGGCGTGTGATGGAGAAAGGCCACTTTGGCGTTTCTGCAAGATGTTATGAACAT
TTTACTTCAGTATGTGGTAAAAAGTTTCGATAGATCAACCAAAGTGATTGATTTCCATTATCCTAATGAG
CTTCTCCAAGAATAAATTGGGAATTGGCAGACCAACCACAAAATTTGGAGGAAATTTTGTGATTGCC
AAACAACCTCTAAAATATGCAATTAACAGGGCATCTAGATACTTCAACTTTCTACTGGTTTGGAA
TATGGTTGGATTAGCAGCAGACTGGCTGACATCAACAGCAAATACTAACATGTTACCTATGAAATTGCT
CCAGTATTTGTGCTTTTGAATATGTCACACTAAAGAAAATGAGAGAAATCATTGGCTGGCCAGGGGGCT
CTGGCGATGGGATATTTTCTCCCGTGGCGCCATATCTAACATGTATGCCATGATGATCGCACGCTTTAA
GATGTTCCAGAAAGTCAAGGAGAAAGGAATGGCTGCTCTTCCAGGCTCATTGCCTTACGCTCTGAACAT
AGTCATTTTTCTCTCAAGAAGGGAGCTGCAGCCTTAGGGATTGGAACAGACAGCGTGATTCTGATTAAT
GTGATGAGAGAGGGAAAAATGATTCCATCTGATCTTGAAGAAGGATTCTTGAAGCCAAACAGAAAGGGTT
TGTTCTTTCTCGTGAGTGCCACAGCTGGAACCACCGTGTACGGAGCATTTGACCCCTCTTAGCTGTC
GCTGACATTTGCAAAAAGTATAAGATCTGGATGCATGTGGATGCAGCTTGGGGTGGGGGATTACTGATGT
CCCGAAAACACAAGTGGAACTGAGTGGCGTGGAGAGGGCCAACCTGTGACGTGGAATCCACACAAGAT
GATGGGAGTCCCTTTGCAGTGTCTGCTCTCTGGTTAGAGAAGAGGGATTGATGCAGAATGCAACCAA
ATGCATGCCTCCTACCTCTTTCAGCAAGATAAACATTATGACCTGCTCCTATGACACTGGAGACAAGGCCT
TACAGTGGGAGCCACGTTGATGTTTTTAAACTATGGCTGATGTGGAGGGCAAAGGGGACTACCGGGTT
TGAAGCGCATGTTGATAAATGTTTGGAGTTGGCAGAGTATTTATACAACATCATAAAAAACCAGAAAGGA
TATGAGATGGTGTGTTGATGGGAAGCCTCAGCACACAAATGTCTGCTTCTGGTACATTCTCCAAGCTTGC
GTACTCTGGAAGACAATGAAGAGAGAATGAGTCGCCTCTCGAAGGTGGCTCCAGTGATTAAGCCAGAAT
GATGGAGTATGGAACCACAATGGTCAGCTACCAACCCTTGGGAGACAAGGTCAATTTCTCCGCATGGTC
ATCTCAAACCCAGCGCAACTACCAAGACATTGACTTCTGATTGAAGAAATAGAACGCCTTGGACAAG
ATTTATAA
    
```

Restriction Sites: Please inquire

ACCN: NM_001134366

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001134366.1</u> , <u>NP_001127838.1</u>
RefSeq Size:	2419 bp
RefSeq ORF:	1758 bp
Locus ID:	2572
UniProt ID:	<u>Q05329</u>
Cytogenetics:	10p12.1
Protein Families:	Druggable Genome
Protein Pathways:	Alanine, aspartate and glutamate metabolism, beta-Alanine metabolism, Butanoate metabolism, Metabolic pathways, Taurine and hypotaurine metabolism, Type I diabetes mellitus
Gene Summary:	<p>This gene encodes one of several forms of glutamic acid decarboxylase, identified as a major autoantigen in insulin-dependent diabetes. The enzyme encoded is responsible for catalyzing the production of gamma-aminobutyric acid from L-glutamic acid. A pathogenic role for this enzyme has been identified in the human pancreas since it has been identified as an autoantibody and an autoreactive T cell target in insulin-dependent diabetes. This gene may also play a role in the stiff man syndrome. Alternative splicing results in multiple transcript variants that encode the same protein. [provided by RefSeq, Oct 2008]</p> <p>Transcript Variant: This variant (2) differs in the 3' UTR, compared to variant 1.</p>