

Product datasheet for **SC325217**

LRRC15 (NM_001135057) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	LRRC15 (NM_001135057) Human Untagged Clone
Tag:	Tag Free
Symbol:	LRRC15
Synonyms:	LIB
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC325217 representing NM_001135057.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCCTTTGGACAAGGCTATGCCACTGAAGCATTATCTCCTTTTGTGGTGGGCTGCCAAGCCTGGGGT
GCAGGGTTGGCCCTACCATGGCTGCCCTAGCGAGTGTACCTGCTCCAGGGCCTCCAGGTGGAGTGCACC
GGGGCACGCATTGTGGCAGTGCCACCCCTCTGCCCTGGAACGCCATGAGCCTGCAGATCCTCAACACG
CACATCACTGAACTCAATGAGTCCCCGTTCTCAATATCTCAGCCCTCATCGCCCTGAGGATTGAGAAG
AATGAGCTGTCGCGCATCACGCCTGGGGCTTCCGAAACCTGGGCTCGTGCCTATCTCAGCCTCGCC
AACAAAGCTGCAGGTTCTGCCATCGGCCTTCCAGGGCCTGGACAGCCTCGAGTCTCTCCTCTG
TCCAGTAACCAGCTGTTGCAGATCCAGCCGGCCCACTTCTCCAGTGCAGCAACCTCAAGGAGCTGCAG
TTGCACGGCAACCACCTGGAATACATCCCTGACGGAGCCTTCGACCACCTGGTAGGACTCACGAAGCTC
AATCTGGCAAGAATAGCCTCACCCACATCTCACCCAGGGTCTTCCAGCACCTGGCAACCTCCAGGTC
CTCCGGCTGTATGAGAACAGGCTCACGGATATCCCATGGGCACCTTTGATGGGCTTGTTAACCTGCAG
GAACTGGCTCTGCAGCAGAACAGATTGGACTGCTCTCCCTGGTCTTCCACAACAACCACAACCTC
CAGAGACTCTACCTGTCCAACAACCACATCTCCAGCTGCCACCCAGCGTCTTATGCAGCTGCCCCAG
CTCAACCGTCTTACTCTTTGGGAATTCCTGAAGGAGCTCTCTCCGGGGATCTTCCGGGCCATGCCC
AACCTGCGGGAGCTTTGGCTCTATGACAACCACATCTCTTCTACCCGACAATGTCTTCAGCAACCTC
CGCCAGTTGCAGGTCCTGATTCTTAGCCGCAATCAGATCAGTTCATCTCCCCGGGTGCCTTCAACGGG
CTAACGGAGCTTCGGGAGCTGTCCCTCCACACCAACGCACTGCAGGACCTGGACGGGAACGTCTCCGC
ATGTTGGCCAACCTGCAGAACATCTCCCTGCAGAACACCCGCTCAGACAGCTCCAGGGAATATCTTC
GCCAACGTCAATGGCCTCATGGCCATCCAGCTGCAGAACACCCAGCTGGAGAATTGCCCTCGGCATC
TTCGATCACCTGGGAAACTGTGTGAGCTGCGGCTGTATGACAATCCCTGGAGGTGTGACTCAGACATC
CTCCGCTCCGCAACTGGCTCCTGCTCAACCAGCCTAGGTTAGGGACGGACACTGTACCTGTGTGTTTC
AGCCCAGCCAATGTCCGAGGCCAGTCCCTCATTATCATCAATGTCAACGTTGCTGTTCCAAGCGTCCAT
GTCCCGAGGTGCCTAGTTACCCAGAAACACCATGGTACCCAGACACACCCAGTTACCCTGACACCACA
TCCGTCTTCTACCACTGAGCTAACCCAGCCCTGTGGAAGACTACACTGATCTGACTACCATTAGGTC
ACTGATGACCGCAGCGTTTGGGGCATGACCCAGGCCAGAGCGGGCTGGCATTGCCGCCATTGTAATT
GGCATTGTGCCCTGGCCTGCTCCCTGGCTGCCTGCGTCGGCTGTTGCTGCTGCAAGAAGAGGAGCCAA
GCTGTCCTGATGCAGATGAAGGCACCCAATGAGTGTTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTAAACGGCCGGC
  
```

Restriction Sites: Sgfl-Mlul
ACCN: NM_001135057
Insert Size: 1764 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001135057.2](#)

RefSeq Size: 5940 bp

RefSeq ORF: 1764 bp

Locus ID: 131578

UniProt ID: [Q8TF66](#)

Cytogenetics: 3q29

Protein Families: Transmembrane

MW: 65 kDa