

Product datasheet for **SC325227**

ZNF652 (NM_001145365) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF652 (NM_001145365) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF652
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC325227 representing NM_001145365.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGAGCCACACAGCCAGTTCCTTGTGAGGAGCTGGTTGAAAACGTGCTGTGCATGTAGCAGGAATGGCA
CAAGAAGATAGCCGTCGTGGTCAAGTGCCATCTTCTTTTATCATGGTGCCAACCAAGAACTTGACCTG
TCCACCAAAGTGTACAAAAGGGAATCAGGAAGTCCTTATTCTGTGTTAGTGGACCAAGATGAGCAAA
CCGCATCTCCATGAAACAGAAGAACAGCCATATTTTCAGGGAGACAAGAGCAGTGTCTGACGTGCATGCT
GTTAAGGAAGACCGGGAGAATTCTGATGACACAGAGGAGGAAGAGGAAGAAGTCTTTACAAAAGGGAG
CAGATCATAGTGGAGGTAACCTTAATAATCAAACATTAATGTATCTAAAGGGGAAAAGGGTGTCTCT
TCTCAGTCCAAGAGACTCCTGTTCTTAAGACAAGCAGTGAAGGAGGAAGAGGAAGAGTGAAGAGAG
GCCACAGATGACAGCAATGACTATGGAGAGAATGAAAAGCAGAAGAAAAGGAGAAGATAGTAGAGAAA
GTCAGCGTTACACAAAGGAGAACCAGGAGAGCTGCCTCTGTTGCCGCAGTACCCTTCCCTACTCCC
AGAACTACAAGAGTTCGTAGGAAGAGTGTAGAGCCACCTAAGCGTAAGAAGCGGGCCACAAGGAGCCC
AAAGCACCAGTCCAGAAAGCTAAGTGTGAAGAGAAAGAGACTCTGACCTGTGAGAAGTGCCCCAGGGTA
TTTAACACTCGCTGGTACCTGGAGAAGCAGATGAACGTTACTCATAGGCCATGCAGATTTGTGATAAA
TGTGGCAAGAAGTTTGTCTGGAAAAGTGAAGTGTCCCTTACCAGCAAACAGACTGTGAAAAAACATT
CAGTGTGTTTCTGTAAACAAATCGTTCAAGAACTCTGGTCCCTTATGAACATATCAAGATCGTCCAT
GGATATGCAGAAAAGAAATTTTCTGTGAAATTTGTGAGAAGAAATTTCTATACCATGGCTCATGTGCGG
AAACACATGGTTGCACACACAAAAGACATGCCATTTACATGCGAAACCTGTGGAAAATCATTCAAACGC
AGTATGCTACTCAAGGTGCATCCTTGCAGCATTCTGGAGAGAAGCCCTTTAGATGCGGAGAAGTGTGAC
GAAAGGTTTTCAGTACAAGTACCAGTACGCTCCACATGAGCATTATGAGGACAAACAGTTTCATG
TGCCAGTGGTGTGCAAGGATTTCAACATGAAGCAGTACTTCGACGAACACATGAAAACACACACTGGA
GAGAAACCCTTTATCTGTGAAATCTGTGCAAAAAGCTTACCAGCCGCCCAACATGAAGAGACACCGC
AGAAGTACACAGGCGAGAAGCCCTATCCATGTGATGTGTGTGGCCAGCGGTTCCGCTTCTCGAACATG
CTTAAGGCCCAAGGAGAAGTGTCTTCCGGTGACCAGCCCGTGAATGTGCCACCTGCTGTCCAGATC
CCACTTACAACCTCCCCAGCCACCCAGTTCTTCTGTGGTGAACACAGCCACAACCCCAACCCCTCCA
ATCAATATGAATCCTGTAAGCACTTCCCCCTCGGCCATCCCCACCCCTTCTCACACCTGCACATC
CACCCACACCCTACCACCCACACCACCTCCCATCCCTCCAGTCCCTCACCTCCCGCCACCTCCAGCT
CTCTTTAAGAGTGAAGCTTTAAATCATAGAGGCCAGAGTGAAGACAACCTTCTGCGGCACCTGGCAGAG
AAGAACAGTTCAGCACAGCATCATTA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001145365

Insert Size: 1821 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001145365.1](#)

RefSeq Size: 11492 bp

RefSeq ORF: 1821 bp

Locus ID: 22834

UniProt ID: [Q9Y2D9](#)

Cytogenetics: 17q21.32-q21.33

MW: 69.7 kDa

Gene Summary:

Functions as a transcriptional repressor.[UniProtKB/Swiss-Prot Function]

Transcript Variant: This variant (1) represents the longest transcript. Both variants 1 and 2 encode the same protein. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.