

Product datasheet for **SC325229**

Calmegin (CLGN) (NM_001130675) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Calmegin (CLGN) (NM_001130675) Human Untagged Clone
Tag:	Tag Free
Symbol:	CLGN
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC325229 representing NM_001130675.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGCATTTCCAAGCCTTTTGGCTATGTTTGGGTCTTCTGTTTCATCTCAATTAATGCAGAATTTATGGAT
GATGATGTTGAGACGGAAGACTTTGAAGAAAATTCAGAAGAAATGATGTTAATGAAAGTGAACCTTCC
TCAGAGATTAATATAAGACACCTCAACCTATAGGAGAAGTATATTTTGCAGAACTTTTATAGTGGA
AGGTTGGCTGGATGGGTCTTATCAAAGCAAAGAAAGATGACATGGATGAGGAAATTTCAATATACGAT
GGAAGATGGGAAATGAAGAGTTGAAAGAAAACCGGTACCTGGTGACAGAGGACTGGTATTAATCT
AGAGCAAAGCATCATGCAATATCTGCTGTATTAGCAAACCATTCATTTTTGCTGATAAACCTTGATA
GTTCAATATGAAGTAAATTTTCAAGATGGTATTGATTGTGGAGGTGCATACATTAACCTCTAGCAGAC
ACTGATGATTTGATTCTGAAAACCTTTTATGATAAAACATCCTATATCATTATGTTTGGACCAGATAAA
TGTGGAGAAGATTATAAACTTCATTTTATCTTCAGACATAAACATCCAAAACCTGGAGTTTTCGAAGAG
AAACATGCCAAAACCTCCAGATGTAGACCTTAAAAAGTTCTTTACAGACAGGAAGACTCATCTTTATACC
CTTGATGAATCCAGATGACACATTTGAGGTGTTAGTTGATCAAACAGTTGTAACAAAGGAAGCCCTC
CTAGAGGATGTGGTTCCTCCTATCAAACCTCCCAAAGAAATTGAAGATCCCAATGATAAAAAACCTGAG
GAATGGGATGAAAGAGCAAAAATTCCTGATCCTTCTGCCGTCAAACCAGAAGACTGGGATGAAAGTGAA
CCTGCCAAAATAGAAGATTCAAGTGTGTTAAACCTGCTGGCTGGCTTGATGATGAACCAAAATTTATC
CCTGATCCTAATGCTGAAAAACCTGATGACTGGAATGAAGACACGGATGGAGAATGGGAGGCACCTCAG
ATTCTTAATCCAGCATGTCGGATTGGGTGTGGTGTGAGTGGAAACCTCCCATGATAGATAACCAAAATAC
AAAGGAGTATGGAGACCTCCACTGGTCGATAATCCTAACTATCAGGGAATCTGGAGTCTCGAAAAATT
CCTAATCCAGATTATTTGGAAGATGATCATCATTCTTCTGACTTCTTTCAGTGTCTTGGTTTAGAG
CTTTGGTCTATGACCTCTGATATCTACTTTGATAATTTTATTATCTGTTGCGAAAAGGAAGTAGCAGAT
CACTGGGCTGCAGATGGTTGGAGATGGAAAATAATGATAGCAAATGCTAATAAGCCTGGTGTATTAATA
CAGTTAATGGCAGCTGCTGAAGGGCACCCATGGCTTTGGTTGATTTATCTTGTGACAGCAGGAGTGCCA
ATAGCATTAACTTCAATTTTGGTGGCAAGAAAAGTAAAGAAAAAACATAAAGATACAGAGTATAAAA
AAAACCGACATATGTATACCACAAAACAAAAGGAGTACTAGAGCAAGAAGAAAAGGAAGAGAAAGCAGCC
CTGAAAAACCAATGGACCTGGAAGAGGAAAAAAGCAAATGATGGTGAATGCTTAAAAAGAAGAG
GAAAGTGAACCTGAGGAAAAGAGTGAAGAAGAAATGAAATCATAGAAGGGCAAGAAGAAAGTAATCAA
TCAAATAAGTCTGGGTCAGAGGATGAGATGAAAGAAGCAGATGAGAGCACAGGATCTGGAGATGGGCCG
ATAAAGTCAGTACGCAAAGAAGAGTACGAAAGGACTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001130675

Insert Size: 1833 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001130675.1](#)

RefSeq Size: 2828 bp

RefSeq ORF: 1833 bp

Locus ID: 1047

UniProt ID: [O14967](#)

Cytogenetics: 4q31.1

Protein Families: Transmembrane

MW: 70 kDa

Gene Summary: Calmegin is a testis-specific endoplasmic reticulum chaperone protein. CLGN may play a role in spermatogenesis and infertility. [provided by RefSeq, Jul 2008]
Transcript Variant: This variant (2) contains an extra exon in the 5' UTR, compared to variant 1.