

Product datasheet for **SC325255**

ZBTB22 (NM_001145338) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZBTB22 (NM_001145338) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZBTB22
Synonyms:	BING1; fru; fruitless; ZBTB22A; ZNF297; ZNF297A
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC325255 representing NM_001145338.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGAGCCATCTCCTCTGTCTCCAGTGGGGCAGCACTTCCCCTGCCGTGTGCTGGCTCCGCCCCCA
CTACCCCTGCCAGCAGCTGCAGTGGTACATGTGTCTTCCCTGAGGTGACCAGTGCCTCTTGGAGTCC
CTCAATCAGCAGCGTCTGCAGGGCCAGCTCTGCGATGTATCTATCAGAGTGCAGGGCCGGGAGTCCGG
GCTCATCGGGCTGTCTGGCTGCCTCCTCCCTTACTTCCATGATCAGGTCTACTCAAAGGCATGACC
TCCATCTCGTGGCCAGTGTGATGGACCCAGGCGCCTTTGAGACTGTCTAGCCTCCGCTTACACTGGC
CGCCTCAGCATGGCTGCTGCTGACATTGTCAACTTCCTTACAGTGGGTCTGTGCTCCAAATGTGGCAC
ATTGTGGACAAGTGCCTGAACACTCCGAGAAGGCCGGGCTCAGTACCACCACCATCACTACTGCT
GCAGCCACCTCTGCTACTGTCCCTGGTGTGGGTGCCATCCGGGAGTGGGGGACTGTGGCCCTGCT
ACCATGGGCTCTGCGGCTCCCATGCCCTCAGCCGGGCCAGTGAATCAATCTCCAGCAGCAGCAAC
TACTTCAGCCCCAGGGAGTCCACTGATTTCTCATCTTCTCCCAAGAGGCATTTGCAGTCTTGCAGTG
GGCAGTGGGGAGCGTGCAGGAGGTGGCCCTGTATCCAGCCCTGTCTGTGGCAGTGGAGGGGCCACA
TCTGAAAAGTGTGCTGGAGGCAGATGAGCTGTGCGATGATGGTGGGGATGGGAGGGGGCAGTGGTT
CCTGGGGCTGGGCTCCGGAGACCCACCTACACACCCCTAGCATCATGCCACAGAAACACTGGGTATAC
GTGAAGCGAGGTGGTAATTGCCAGCGCAACACCCCTGGTTCCTCAAGACCCAGATCTGGAGGAGGAA
GAGGAGGAGGAAGATCTGGTGTGACCTGTGAGGATGATGAAGATGAAGAACTAGGGGTAGCTCCAGG
GTTCCAGTGGGGGAGGGCCTGAGGCTACCCTCAGCATAAGTGTGTCGTACCTGAGTGAAGCCCCA
GACAAGGGGAGGAGCAGGTCAACTTCTGTGAGTCCCAATGACTTTGGCCATATGAGGGTGGGGT
CCTGTGGCAGTCTTGTGACTCAGGGGGCCAACTCCCTTCTCTATGCCCCCTCCACCCTCCTCGA
CCGCTCCTTCCCTTGGACATGCAGGGCAACCAGATCCTGGTCTTCCCGTGTCTCTTATCCTCATCC
TCACAGGCTCCTGGCAACCACCAGGGAACCAAGCAGAACACGGGGCAGTGACCGTGGGGGACGTCG
GTGGGGAGCCTGGGTGTGCCGGTAGCGTTGGTGGGGTCCCTGGAGGGACTGGCAGTGGGACGGGAAT
AAGATCTTTCTGTGCCATTGTGGGAAGGCCTTCTCCACAAGAGCATGCGGGACCGGCACGTGAACATG
CACCTCAATCTGCGGCCGTTTACTGCCCCGTGTGCAACAAAAGTTCAAGATGAAGCACCATCTGACT
GAGCAGATGAAGACGCACACAGGTCTCAAGCCCTACGAGTGGGAGTCTGCGCCAAGAAGTTCATGTGG
CGAGACAGCTTATGCGCCACCGAGGACTGTGAGCGCCGACCGCCTGGGCGGGTTCGGGGCCGTA
CCTGGGCTGGGACTCCCAGGGGCCATCCTTGGCGTCCAAGAGAGAGTCTCCCGAGTGGGCGGGGGC
AGCGGCGACGAAGCGAGTGCAGGCCACGCCCGTCCAGCAGACGTGTCTGGTCCCCACCCAGAGTCCAC
AAGGTGGAGATGGGCTTCGGTGGAGGTGGAGGAGCAAATGTA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001145338

Insert Size: 1905 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001145338.1](#)

RefSeq Size: 2663 bp

RefSeq ORF: 1905 bp

Locus ID: 9278

UniProt ID: [O15209](#)

Cytogenetics:	6p21.32
Protein Families:	Transcription Factors
MW:	65.6 kDa
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant differs in the 5' UTR compared to variant 1. Both variants 1 and 2 encode the same protein.