

Product datasheet for **SC325268**

ZNF844 (NM_001136501) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF844 (NM_001136501) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF844
Vector:	<u>pCMV6 series</u>



[View online »](#)

Fully Sequenced ORF: >NCBI ORF sequence for NM_001136501, the custom clone sequence may differ by one or more nucleotides

```

ATGGACTTGGTGGCTTTTCGAGGATGTGGCTGTGAACTTCACCCAGGAGGAGTGGTCTTTG
CTGGATCCTTCCCAGAAGAATCTCTACAGAGAAGTGATGCAGGAAACCTTGAGGAATCTG
GCCTCCATAGGAGAAAAATGGAAAGACCAGAACATTGAAGATCAGTACAAAAATCCCAGG
AATAATCTAAGAAGTCTTCTGGGAGAGAGATTGATAAAAATACAGAAGAAAAATCATTGT
GGAGAACTTCTAGCCAGATTCCAGATGACACACTGAACAAAAAACTTCTCCTGGAGTA
AAATCATGTGAAAGCAGTGTGTGGAGAAGTCTTCGTGGGTCTTCTCCCTTAATAGG
CACATTAGAGCTGACACTGCACACAAGCCGTCTGAGTATCAGGAATATGGACAGGAGCCA
TATAAGTGTCAACAACGTAAGAAAGCCTTCAGATGTCACCCCTCCTTTCAATGCAAGAA
AAGGCTCACACTGGAGAAAACTCTATGATTGTAAGAAATGTGGAAAAACCTTCATATCC
CATTCAAGCATTCAAAGACACATGATAATGCACAATGGAGATGGAACCTATAAATGTAAG
TTTTGTGGAAAGCCTGCCCTTGTCTCAGCATATATCTTATACATGAACGAGTTCACACT
GGAGAGAAACCATATAAATGTAACAATGTGGTAAAGCCTTAGTTATTCAACTCCCTT
CAAATACATGAAAGAACTCACACTGGAGAGAAGCCTTATGAATGTAAGGAATGTGGGAAA
GCATTCGGTAGTCCCAATTCCCTTTATGAACATAGAAGAACTCACACTGGAGAGAAGCCA
TATGAATGCAACAATGTGGAAAAGCCTTCAGATGGTTCATTCTTTCAATACATGAA
AGAACTCACAGTGAGGAGAAGGCTTATGAATGTACCAATGTGGGAAAGCATTCAAGTGT
CCCAGTTATCTTTGTAGACATGAAGTGACCCACTCTGGGAAAAAGCCCTGTGAATGTAAA
CAATGTGGGAAAGCATTATCTTAACTTTCAAAGACACATGAAAATGCACACTAGA
ATGAGACCTTATAAATGTAAGACTGTGGAAAAGCCTTTGATTCTCCAGTTCGTTTTGAA
GACATGAAAGAATCACACTGGAGAGAAACCTTATGAATGCAAGCACTGTGGTAAAGCCT
TCAATCGTTCAGTTCCTTTCACTATCATGAAAGGACTCACACTGGAGAGAAACCTATG
AATGTAAGCAGTAGTAAAGCCTTCATTTCTTCCACTTCTTTGATATCATGAAAGGA
CTCACACTGGAGAGAAACCGTATGAGTGTAAAGCAATGTGGGAAAGCCTTCAGATCTGCC
CACACCTTCAAATGCATGGAAGGACTCACACTCAAGAGAAACCTATGAATGTAAGCAGT
GTGGTAAAGCCTTCATTTTTTCCACTTCTTTGATATCATGAAAGGACTCACACTGGAG
AGAAACCTATGAGTGTAAAGCAATGTGGGAAAGCCTTCACATCTGCCTCACACCTTCAA
TGCATGAAAGGACTCACACTGGAAAGCACTGTATGAATCTAAACAATGTAAAAACCT
TTGGATCTGTGAGAAACCTTCAAATTCATGAAAAGACACACCCTGGAGAGAAACCTATA
AGGAATATGGAAAAGCATTCAACAATTTCTTCTTTCAAATACATGCAACAATGCACA
GAGGACAGAATGCCTATGAATGTAAGAGTGTGACAAAGCATTATATCTGCCAAGATCC
TTCCAGTACATGCAAGAACACACCCTGGAGAGAAACCTATGAATGTAAGGAATGCGGAA
AAGCGTCAATTATTTTTCTTCTTTGCGTATACACAAAAGGATGCACACTGGAGAGAAAC
CATATTAATGTAAGGATTGTGGGAAAGCATTAGTTTGCCTGGTTCCTTTGATAGACATA
AAAGGGCTCACACTGGAG

```

Restriction Sites: Please inquire

ACCN: NM_001136501

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001136501.1](#), [NP_001129973.1](#)

RefSeq Size: 2834 bp

RefSeq ORF: 2001 bp

Locus ID: 284391

UniProt ID: [Q08AG5](#)

Cytogenetics: 19p13.2

Gene Summary: May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]
Transcript Variant: This variant (1) represents the longer transcript and encodes the protein.