

Product datasheet for **SC325287**

ZNF234 (NM_001144824) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF234 (NM_001144824) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF234
Synonyms:	HZF4; ZNF269
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC325287 representing NM_001144824.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGACCACATTC AAGGAGGACTGACCTTCAAGGATGTGGCTGTGGTCTTCACTGAGGAGGAGCTGGGG
GTGGGGCATCACCCCTTCAAACATGATGTATTCCCTTTAGAAAAGGAAAAAAGCTTGATATAATGAAG
ACAGCAACTCAAAGAAAAGGAAATCAGCAGAC AAGATCCAAAGTGAGGTGGAGACTGTTCCAGAAGCA
GGACGACATGAAGAGCTTACTGGGGCAAATCTGAAACAGATTGCAAGTGATTTAATCAAGTATGAA
GACTCTATGATAAGTATTCTCGGTTCCCAGACAAGGTGATTTGTCTGCCAGGTAGGGCAGGACTA
TATACAACTCACACAGGACAGAAATTTACCAATGTGATGAGTACAAAAATCCTTCACTGATGTCTTC
AACTTTGATCTTCATCAACAGTTACTCAGGAGAGAAGTCTCATACATGTGATGAGTGTGAAAAAGC
TTCTGTTACATCTCAGCCCTTCATATTCATCAAAGAGTCCACATGGGAGAGAAATGCTATAAGTGTGAC
GTGTGTGGTAAGGAATTTAGTCAGAGCTCACATCTTCAAACATCAGAGAGTCCACACTGTAGAGAAA
CCATTCAAATGTGTGGAATGTGGAAAGGCTTCAGTCGTAGATCAACACTTACTGTACATTGCAAAATTA
CACTCAGGAGAGAAACCTTACAATTGTGAGGAATGTGGAAGGGCCTTCATACATGCTCCCATCTTCAG
GAACATCAGAGAAATTCATACTGGGGAGAAACCATTCAAATGTGATACATGTGGTAAGAACTTCCGTCGT
AGATCAGCACTTAATAATCATTGCATGGTCCACACAGGAGAGAAACCATACAAATGTGAGGACTGTGGT
AAGTGTTCACCTGTAGCTCAAACCTTCGTATCCATCAAAGGGTCCACACAGGAGAGAAACCTTACAAG
TGTGAAGAATGTGGTAAGTGCTTTATTCAGCCTTACAATTT CAGGCCATCGGAGAATCCACACTGGA
GAGAAACCATACGTATGTAAGTGTGTGGTAAGGGTTTCATTTACAGTTCAAGTTTT CAGGCCATCAG
GGAGTCCACACTGGAGAGAAAGCCATACAAATGCAATGAGTGTGGGAAGAGCTTCAGAATGAAAATTCAT
TATCAAGTGCATCTGGTAGTCCACACAGGGGAAAAACCTATAAAATGTGAAGTATGTGGTAAAGCCTTC
CGTCAGAGTTCATATCTTAAAATCCATCTGAAAGCACATAGTGTACAGAAACCTTTTAAAGTGTGAAGAG
TGTGGGCAGGGCTTCAATCAGAGCTCACGACTTCAGATTCACCAGCTGATCCATACCGGTGAGAAACCA
TACAAATGTGAAGAGTGCAGAAAGGGATTTAGTCGTAGAGCAGATCTTAAAATTCATTGTAGGATCCAC
ACAGGGGAGAAACCATATAATTGTGAGGAGTGTGGGAAGGTCTTCAGTCAGGCCTCGCATCTTCTAACC
CATCAGAGAGTTCACAGTGGGAAAAACCTTTAAATGTGAAGAGTGTGGGAAGAGCTTCAGTCGGAGT
GCACACCTTCAAGCCCATCAAAGTCCACACTGGAGAAAAGCCATACAAATGTGGGAGTGTGAAAG
GGCTTCAAGTGGAGCTTGAACTTGACATGCATCAGAGGGTGCACACAGGAGAAAAACCTATACATGT
GGGGAGTGTGGGAAGCACTTCAGTCAGGCCTCAAGTCTCAAACCTTCATCAGAGTGTCCACACAGGAGAG
AAACCATACAAATGTGATGTATGTGGTAAAGTCTTCAGTCGGTCTTCACAATTACAGTATCATAGGCGA
GTTCACTGGGAAAAACCTTACAAATGTGAGATATGTGGTAAGAGGTT CAGCTGGCGATCAAATCTT
GTAAGTCAACAAAAATTCATGCTGCTGGTACATTTTATGAAAATGATGAGAATAGTAAGAACATCAGA
GAGTTGTCAGAGGGAGGAAGTTCTACAAGGTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001144824

Insert Size: 2103 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001144824.1
RefSeq Size:	4482 bp
RefSeq ORF:	2103 bp
Locus ID:	10780
UniProt ID:	Q14588
Cytogenetics:	19q13.31
Protein Families:	Transcription Factors
MW:	80.6 kDa
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (2) differs in the 5' UTR compared to variant 1. Variants 1 and 2 both encode the same protein. Sequence Note: This RefSeq record was created from transcript and genomic sequence data because no single transcript was available for the full length of the gene. The extent of this transcript is supported by transcript alignments.