

## Product datasheet for **SC325288**

### **CAPN8 (NM\_001143962) Human Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	CAPN8 (NM_001143962) Human Untagged Clone
Tag:	Tag Free
Symbol:	CAPN8
Synonyms:	nCL-2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC325288 representing NM\_001143962.  
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCAGCCCAGGCAGCTGGTGTATCTAGGCAGCGGGCAGCCACTCAAGGCTTTGGCTCCAACAAAAC
GCTTTGAAGTACTTGGCCAGGATTTCAAGACCCTGAGGCAACAGTGTCTGGACTCAGGGGTCTATTT
AAGGACCCTGAGTTCCAGCATGTCATCAGCTTTGGGCTACAAGGATCTTGACCAGGCTCTCCGCAA
ACTCAAGGCATCATCTGGAAGCGGCCACGGAGTTGTGTCCAGCCCTCAGTTTATCGTTGGTGGAGCC
ACGCGCACAGACATTTGTCAGGGTGGTCTAGGTGACTGCTGGCTTCTGGTGCCTATGCCTCCCTGACC
CTGAATGAAGAGCTGCTTACCGGGTGGTCCCCAGGGACCAGGACTTCCAGGAGAATAATGCGGGAATC
TTTCACTTTCAGTCTGCGAGTACGGAGAGTGGTGGAGGTGGTCTTACGACAGGCTGCCACCAAG
AATGGACAGCTGCTCTTCTACTCGGAACAAGCAATGAATTCTGGAGTGCCTGCTGGAGAAAGCC
TATGCCAAGCTTAAATGGTTGTTATGAGGCTCTCGCTGGAGTTCCACAGTGGAGGGGTTGAGGATTC
ACAGGTGGCATCTGAGTTTTATGACCTGAAGAAACCACAGCCAATCTATATCAGATCATCCGGAAG
GCCCTCTGTGCGGGTCTCTGCTGGGCTCTCATTGATGTCTCAGTGCAGCCGAGCCGAAAGCCATC
ACCAGCCAGAAGCTGGTTAAGAGTCATGCTACTCTGTCACTGGAGTCAAGAGGTGAATTTCCAGGGC
CATCCAGAGAAGCTGATCAGACTCAGGAATCCATGGGGTGAAGTGGAGTGGTCCGGAGCCTGGAGCGAT
GATGCACCAGAGTGGAAATCAGATAGACCCCGCGGAAGGAAGAACTGGACAAGAAAGTTGAGGATGGA
GAATTCGGATGTCACCTTCAGATTTCTGTGAGGAGTCTCTCGTTGGAGATCTGCAACCTGTCCCCG
GACTCTCTGAGTAGCGAGGAGGTGCACAAATGGAACCTGGTCTGTTCAACGGCCACTGGACCCGGGGC
TCCACAGCTGGGGCTGCCAGAACTACCCAGCCACGACTGGACCAATCCCAGTCAAAAATCCGTTTG
GATGAAGTGGATGAGGACCAGGAGAGAGCATCGGTGAACCTGCTGTACAGTGTCTGGGCCTGATG
CAGAAAAATCGCAGGTGCGGGAAGCGGATAGGACAAGGCATGCTTAGCATCGGCTATGCCGCTACCCAG
GTTCCCAAGGAGCTGGAGAGTCACACGGACGCACACTTGGGCCGGGATTTCTTCTGGCCTACCAGCCC
TCAGCCCGCACCAGCACCTACGTCAACCTGCGGGAGGTCTCTGGCCGGGCCCGGCTGCCCTGGGGAG
TACCTGGTGGTGCATCCACATTTGAACCTTCAAAGACGGCGAGTTCTGCTTGAGAGTGTCTCAGAG
AAGAAGGCCAGGCCCTAGAAATGGGGATGTGGTAGCTGGAACCCATATGAGCCACATCCCAGTGGAG
GTGGATCAGGAAGATGACCAGTTCAGGAGGCTGTTGAGAAGTTGGCAGGGAAGGATTCTGAGATTACT
GCCAATGCACTCAAGATACTTTGAATGAGGCGTTTTCCAAGAGAACAGACATAAAATTCGATGGATTC
AACATCAACACTTGCAGGAAATGATCAGTCTGTTGGATAGCAATGGAACGGGCACCTTTGGGGCGGTG
GAATTAAGACGCTCTGGCTGAAGATTCAGAAGTATCTGGAGATCTATTGGGAAACTGATTATAACCAC
TCGGGCACCATCGATGCCACGAGATGAGGACAGCCCTCAGGAAGGCAGGTTTACCCTCAACAGCCAG
GTGCAGCAGACCATGTCCTGCGGTATGCGTGCAGCAAGCTCGGCATCAACTTTGACAGCTTCGTGGCT
TGTATGATCCGCTGGAGACCCTCTTAAACTATTCAGCCTTCTGGACGAAGACAAGGATGGCATGGTT
CAGCTCTCTGCGCGAGTGGCTGTGCTGCGTGTGGTCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

**Plasmid Map:**


**ACCN:** NM\_001143962

**Insert Size:** 2112 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001143962.1](#)

**RefSeq Size:** 2427 bp

**RefSeq ORF:** 2112 bp

**Locus ID:** 388743

**UniProt ID:** [A6NHC0](#)

**Cytogenetics:** 1q41

<b>Protein Families:</b>	Druggable Genome
<b>MW:</b>	79.1 kDa
<b>Gene Summary:</b>	Calcium-regulated non-lysosomal thiol-protease. Involved in membrane trafficking in the gastric surface mucus cells (pit cells) and may involve the membrane trafficking of mucus cells via interactions with coat protein. Proteolytically cleaves the beta-subunit of coatomer complex (By similarity).[UniProtKB/Swiss-Prot Function]