

## Product datasheet for **SC325322**

### **NELL1 (NM\_201551) Human Untagged Clone**

#### **Product data:**

|                           |  |
|---------------------------|--|
| Product Type:             | Expression Plasmids                    |
| Product Name:             | NELL1 (NM_201551) Human Untagged Clone |
| Tag:                      | Tag Free                               |
| Symbol:                   | NELL1                                  |
| Synonyms:                 | IDH3GL; NRP1                           |
| Mammalian Cell Selection: | Neomycin                               |
| Vector:                   | pCMV6-Entry (PS100001)                 |
| E. coli Selection:        | Kanamycin (25 ug/mL)                   |




[View online »](#)

**Fully Sequenced ORF:** >SC325322 representing NM\_201551.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGCCGATGGATTTGATTTTAGTTGTGTGGTTCTGTGTGCTGACTGCCAGGACAGTGGTGGGCTTTGGG
ATGGACCCTGACCTTCAGATGGATATCGTCACCGAGCTTGACCTTGTAACACCACCCTTGAGTTGCT
CAGGTGTCTGGAATGCACAATGCCAGCAAAGCATTATTTTCAAGACATAGAAAAGAGAGATCCATGCA
GCTCCTCATGTGAGTGAGAAATTAATTCAGCTGTTCCGGAACAAGAGTGAATTCACCATTTTGGCCACT
GTACAGCAGAAGCCATCCACTTCAGGAGTGATACTGTCCATTCGAGAAGTGGAGCACAGCTATTTTGAA
CTGGAGAGCAGTGGCCTGAGGGATGAGATTCGGTATCACTACATACACAATGGGAAGCCAAGGACAGAG
GCACCTTCCTTACCGCATGGCAGATGGACAATGGCACAAGTTGCACTGTGAGTTAGCGCCTCTCATCTC
CTGCTCCATGTCGACTGTAAACAGGATTTATGAGCGTGTGATAGACCCTCCAGATACCAACCTTCCCCA
GGAATCAATTTATGGCTTGGCCAGCGCAACCAAAAGCATGGCTTATCAAAGGGATCATCCAAGATGGG
AAGATCATCTTTATGCCGATGGATATATAACACAGTGTCCAAATCTAAATCACACTTGCCCAACCTGC
AGTGATTTCTTAAGCCTGGTGAAGGAATAATGGATTTACAAGAGCTTTTGGCCAAGATGACTGCAAAA
CTAAATATGCAGAGACAAGACTTAGTCAATTGAAAACTGTCATTGTGAGAAGACTTGTCAAGTGAGT
GGACTGCTCTATCGAGATCAAGACTCTTGGGTAGATGGTGACCATTGCAGGAACTGCACTTGCAAAAAGT
GGTGCCGTGGAATGCCGAAGGATGTCTGTCCCCCTCTCAATTGCTCCCCAGACTCCCTCCCAGTGCAC
ATTGCTGGCCAGTGTCTAAGGTCTGCCGACCAAAATGTATCTATGGAGGAAAAGTTCTTGCAAGGGC
CAGCGGATTTTAACCAAGAGCTGTGGGAAATGCCGAGGTGGAGTTTTAGTAAAAATTACAGAAATGTGT
CCTCCTTTGAACTGCTCAGAAAAGGATCACATTTCTCTGAGAATCAGTGTCTGCCGTGTCTGTAGAGGT
ATAAATTTTGTGCAAGGACCTAAATGTGGTGAAAACCTCAGAGTGCAAAAACCTGGAATACAAAAGCT
ACTTGTGAGTGCAAGAGTGGTTACATCTGTCTCCAGGGAGACTCTGCCTACTGTGAAGATATTGATGAG
TGTGCAGCTAAGATGCATTACTGTGATGCCAATACTGTGTGTGCAACCTTCTGGGTTATATCGTGT
GACTGTGTCCCAGGATACATTCTGTGGATGACTTCTCTGTACAGAACACGATGAATGTGGCAGCGGC
CAGCACAACCTGTGATGAGAATGCCATCTGCACCAACACTGTCCAGGGACACAGCTGCACCTGCAAAACG
GGCTACGTGGGGAACGGGACCATCTGCAGAGCTTTCTGTGAAGAGGGCTGCAGTACGGTGGAACGTGT
GTGGCTCCCAACAAATGTGTCTGTCCATCTGGATTCACAGGAAGCCACTGCGAGAAAGACATTGATGAA
TGTGCCTAAGAACTCACACCTGTTGGAACGATTCGCCTGCATCAACCTGGCAGGGGGCTTTGACTGT
CTCTGCCCTCTGGGCCCTCTGCTCTGGTGACTGTCTCATGAAGGGGGCTGAAGCACAATGGCCAG
GTGTGGACCTTGAAAGAAGACAGGTGTTCTGTCTGCTCCTGCAAGGATGGCAAGATATTCTGCCGACGG
ACAGCTTGTGATTGCCAGAATCCAAGTGTGACCTATTCTGTTGCCCAGAATGTGACACCAGAGTACA
AGTCAATGTTTAGACCAAAATGGTCACAAGCTGTATCGAAGTGGAGACAATTGGACCCATAGCTGTGAG
CAGTGTGGGTGTCTGGAAGGAGAGGTAGATTGCTGGCCACTCACTTGCCCAACTTGAGCTGTGAGTAT
ACAGCTATCTTAGAAGGGGAATGTTGTCCCCGCTGTGTCAGTGACCCCTGCCTAGCTGATAACATCACC
TATGACATCAGAAAACCTGCCTGGACAGCTATGGTGTTCACGGCTTAGTGGCTCAGTGTGGACGATG
GCTGGATCTCCCTGCACAACCTGTAATGCAAGAATGGAAGAGTCTGTTGTTCTGTGGATTTGAGTGT
CTTCAAATAATTGA
ACGCGTACGCGGCCGCTCGAGCAGAAAACCTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:** SgfI-MluI

**Plasmid Map:**


**ACCN:** NM\_201551

**Insert Size:** 2292 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_201551.2](#)

**RefSeq Size:** 3129 bp

**RefSeq ORF:** 2292 bp

**Locus ID:** 4745

**UniProt ID:** [Q92832](#)

**Cytogenetics:** 11p15.1

**Protein Families:** Secreted Protein

**MW:** 84.4 kDa

**Gene Summary:** This gene encodes a cytoplasmic protein that contains epidermal growth factor (EGF)-like repeats. The encoded heterotrimeric protein may be involved in cell growth regulation and differentiation. A similar protein in rodents is involved in craniosynostosis. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Dec 2013]  
Transcript Variant: This variant (2) lacks an alternate in-frame exon compared to variant 1. The resulting isoform (2) has the same N- and C-termini but is shorter compared to isoform 1.