

Product datasheet for **SC325326**

ZNF438 (NM_001143769) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF438 (NM_001143769) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF438
Synonyms:	bA330O11.1
Vector:	<u>pCMV6 series</u>



[View online »](#)

Fully Sequenced ORF: >NCBI ORF sequence for NM_001143769, the custom clone sequence may differ by one or more nucleotides

```

ATGCTGCCATGTCATTACCCATCACGCTCTGATCAGGTGAATCTGGGACCCTCCATCAAC
TCCAAGCTGCTGGGGATGTCCACCCAGAACTATGCCCTGATGCAGTTGCTGGCCAGGAG
GGGACATTTTCTTGTGTCTGCCACATGTTGCCTCAGCTCAGCCAATTCAGAAACCC
AGAATGTCCTACCTGAAAACCTGAAAACCTCCTATTCTAGATATCAACCCCTAGAAAT
AGCAAAGCATCAAGAAAGAAACCCATCCTGATCTTTCCTAAGAGTGGCTGTAGCAAAGCT
CCTGCCCAAACCCAAATGTGTCCTCAGATGTCCCCTTCCCACCTCACCACCTGAACTC
CTGTACAAAACCCAGTCCATTTGAGGAAGTACCATCACTAGAGCAAGCCCCAGCCAGCATT
AGCACAGCTGCGCTGACCAATGGAAGTGACCATGGGGACTTGAGACCACCAGTGACCAAC
ACCCATGGCAGTCTGAACCCTCTGCTACCCCAGCATCATCCACACCAGAGGAGCCTGCC
AAGCAGGACCTCACAGCTCTTTCAGGGAAAGCACACTTTGTAAGCAAGATAACATCTAGT
AAACCTTCTGCTGTTGCCAGTAAAAATTTAAGAACAAGTTGATCTTGCAAAAACCATG
ACCAATTTATACCAACCATTCTGGCAATGCAGTTCAGTTGATCTTTCAGTCCCCAAA
GGGAAACTGCCAATCCCACCTACTCAAGAATGAAGACAATGGAGTTTACAAAATCAAA
TCAGATGCTAACATTGCAGGTTTTTCTTTACCAGGACCTAAGGCCGACTGTGATAAGATA
CCCTCCACCACAGAAGGCTTTAATGCAGCCACCAAGGTGGCAAGCAGGCTACCTGTTCCA
CAAGTGTACAGCAGAGTGCCTGTGAAAGTGCCTTTTGTCCACCCACAAAACCTTGATCTT
AACCACAAAACAAAACCTGAACAGTGGAGCAGCAAAGAGAAAAGGAAGAAAACGGAAGGTA
CCAGATGAAATTTTGGCATTTCAGGGAAAAAGGAGGAAATATATCATTAAAGTGTAGA
GATGGTAAAGAAAGAGTAAAAATGATCCCAAGAATTCAGAGACAAAAGCTGGGGACC
CTGAAAAAATACCGTAGCATTATGCCCAAACCTATCATGGTCATACCCACTTTGGCCTCC
CTGGCTTCTCCAACACTACACTACAGTCCCAGATGCTTGGGGGCTAGGACAGGATGTTTTG
TTAAATAATTCACTCACTCCTAAATATCTTGGCTGTAAGCAAGACAACAGCTCTTCCCCT
AAGCCCAGCTCCGTGTTCAAAAATGGATTCTCTGGCATTAAAGAAGCCTTGGCACAGATGT
CACGTCTGCAACCACCTTCCAGTTCAAACAGCACCTTCGAGACCACATGAATACACAC
ACCAACAGACGCCCTTACAGTTGTCGGATTTGTCGCAAGTCTATGTACGTCCTGGCAGC
CTGAGCACACACATGAACTTCATCATGGTGAGAACCGTCTGAAGAACTCATGTGTTGT
GAGTTTTGTGCAAAAAGTGTGGCCACATCCGAGTCTATTTGGCCATCTGAAAGAAGTG
CATAGGGTTGTGATCAGCACTGAGCCTGCGCCCAGTGAAGTGCAGCCAGGAGACATACCA
AAGAACAGAGACATGAGTGTGCGAGGCATGGAGGGATCATTGGAGAGGGAAAACAAGTCC
AACCTGGAAGAAGACTTCTTCTAAACCAGGCAGACGAAGTCAAATTACAATCAAATGT
GGTCGTTGTCAGATTACTGCTCAGTCTTTTGCAGAAATAAAATTTCAATTAATGATGTT
CATGGAGAGGAAATTGAGGGCAGGCTACAAGAAGGGACCTTCCCAGGAAGCAAGGGGACT
CAGGAAGAGTTGGTGCAGCACGCTAGCCCCGACTGGAAAAGGCATCCTGAGAGAGGGAAG
CCGGAGAAGGTTCAATCTCCTCCGAGGAATCACATGCATGTCCAAGACTGAAAAGGCAG
CTCCACCTTATCAGAATGGCGTGGAAATGCTCATGGAATGAAGGACCCAGTCAGGA
ACCAACAAGCCAAGGGAAACCTGCCAGGGCCCTGAGTGTCTTGGCCTCCACAGTTTCTC
TTGTGGTCCCATTACAGGCTTTAACTGCCTGCTTTGTGCAGAGATGCTGGGACGGAAGAG
GAGCTCTCCACCACTGGAAGCACCAGCATAACTGTGAGGACCCTTCCAAACTGTGGGCT
ATTTTAAATACGGTCTCAACCAGGGAGTGATCGAACTTCCAGTGAAGCTGAGAAA

```

Restriction Sites: Please inquire

ACCN: NM_001143769

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001143769.1 , NP_001137241.1
RefSeq Size:	3431 bp
RefSeq ORF:	2340 bp
Locus ID:	220929
UniProt ID:	Q7Z4V0
Cytogenetics:	10p11.23
Gene Summary:	Isoform 1 acts as a transcriptional repressor.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (5) has multiple differences in the 5' UTR, which results in the primary ORF beginning at a downstream start codon, compared to variant 1. The resulting protein (isoform c) has a shorter N-terminus, compared to isoform a.