

Product datasheet for **SC325328**

KLRAQ (PPP1R21) (NM_001135629) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	KLRAQ (PPP1R21) (NM_001135629) Human Untagged Clone
Tag:	Tag Free
Symbol:	KLRAQ
Synonyms:	CCDC128; KLRAQ1
Vector:	<u>pCMV6 series</u>



[View online »](#)

Fully Sequenced ORF: >NCBI ORF sequence for NM_001135629, the custom clone sequence may differ by one or more nucleotides

```

ATGGCCTCGGCTGAGTTGCAGGGGAAGTACCAGAAGCTGGCTCAGGAGTACTCGAAGCTT
CGGGCTCAGAATCAGGTTCTGAAAAAGGTGTTGTGGATGAACAAGCAAATTCGACGCT
TTAAAGGAGCAACTGAAAATGAAGGATCAGTCATTGAGAAAACACAACAGGAAATGGAC
AGTTTGACATTTGAAAATCTGCAGCTTGCCAAGAGGGTAGAACTACTTCAAGATGAACTA
GCTCTAAGTGAACCACGAGGCAAGAAAAACAAGAAAAGTGGAGAATCTTCTTCTCAGTTG
AGTCAAGAGCAGAAGAGTGCTTTTGATGAAGATCTGCAAAAAGAAGATAGAAGAGAATGAA
CGTTGTCATATACAATTTTTTGAAGCTGATGAGCAGCACAAAGCATGTGGAAGCAGAGCTG
AGGAGTCGACTGGCCACTCTGGAGACAGAAGCAGCCAGCACCAAGCTGTGGTTGACGGT
CTCACCCGGAAGTACATGGAAACCATTGAGAAGCTGCAGAACGACAAGGCTAAACTAGAA
GTGAAATCTCAGACTCTAGAAAAGGAAGCAAGGAATGTCGACTTCGAACGGAAGAATGT
CAATTACAGTTAAAGACTCTTCATGAAGATTTGTCAGGTAGATTAGAGGAATCCTTATCA
ATCATCAATGAAAAAGTACCTTTAATGATACAAAATATAGTCAGTACAACGCTCTGAAC
GTTCCACTCCACAATAGGAGACACCAGCTGAAGATGCGAGATATTGCTGGGCAGGCCCTG
GCTTTTGTTCAGGATCTTGTGACGGCTCTTCTAACTTTTCATACCTACACAGAACAGAGG
ATTCAAATTTTTCTGTTGATTCTGCCATTGACTATATCTCCATTGAATCAGAAGTTC
TCACAATACCTTCATGAAAATGCGTCTATGTCCGCCCTCTTGAGGAAGGAATGCTTCAT
TTATTTGAAAGTATCACTGAGGATACTGTGACTGTCTTGAGAGCAACTGTGAAATGAAA
ACTTTTTCAGAACACTTAACTCCTACATATGTTTTCTTAGGAAGATTCTTCCCTATCAG
TTAAAAAGTTTAGAAGAAGAAATGGAATCCTCTCTTTGCACATCTGCGTTAAGAGCCAGG
AATCTAGAGCTGTCCAGGACATGAAAAAATGACAGCTGTGTTTGAGAAGCTGCAGACT
TACATAGCTCTTCTTGCCTTGCCAAGTACAGAGCCAGATGGACTCCTTCGGACAAAAC
AGTTCGTGTGTTAACAAATGTTGGTGTCTGCTGTCATGGATTTTCATGACGTTATGAAAGAT
ATTTCCAAACATTATAGTCAAAAAGCTGCAATAGAGCATGAACTTCCAACAGCAACACAG
AAGCTGATAACAATAATGACTGTATCCTGTGCATCAGTAGTGGCATTAAACAAATGGAGCA
GGAAAGATTGCATCCTTCTCAGCAACAATTTGGACTACTTCATTGCTTCACTGAGCTAT
GGACCTAAGGCAGCGAGTGGATTATTAGTCTCTTTCAGCTGAATGCATGCTACAGTAT
AAGAAAAAGCTGCTGCCTATATGAAGTCTTTGAGAAAGCCCTCTTGGAGTCTGTGCCT
TATGAAGAAGCACTGGCAAACCCCGCATCCTTCTCAGCTCTACTGAAAGTCGAGAAGGC
CTTGACACAGCAAGTTCAACAGAGTTTGGAAAAGATTTCTAACTGGAGCAGGAAAAAGAA
CATTGGATGTTGGAAGCACAATTAGCCAAAATCAAGCTAGAGAAAAGAAAACAGCGAATT
GCAGATAAGCTGAAGAATACAGGTAGTGCCAGCTGGTTGGGCTGGCCCAGGAAAATGCT
GCTGTGTCAAATACTGTGGCCAGGATGAAGCCACAGCTAAGGCTGTGTTGGAGCCATT
CAGAGCACCAAGTCAATTTGGGACTTTAACCAGGACATCTGACAGTGAGGTTCCAGATGTG
GAATCTCGTGAAGACTTAATTAATACTACTACATGGCAAGGATAGTGGAACTTACGTCT
CAGTTGCAGCTGGCTGACAGTAAGTCAGTGCATTTTTATGCCGAGTGCCGAGCACTGTCT
AAAAGACTGGCCTTGGCTGAAAAGTCTAAGGAAGCATTGACAGAAGAAAATGAAACTTGCC
AGTCAGAACATCAGCAGACTTCAGGATGAGCTGACAACCTACCAAGAGGAGTTACGAGGAT
CAGTTAAGTATGATGAGTGACCACCTGTGCAGCATGAATGAGACATTATCTAAACAGAGA
GAAGAGATTGACACACTAAAGATGTCCAGTAAGGGGAATTCTAAAAAGAACAAGAGTCGA

```

Restriction Sites: Please inquire

ACCN: NM_001135629

Insert Size: 3181 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001135629.1 , NP_001129101.1
RefSeq Size:	3181 bp
RefSeq ORF:	2343 bp
Locus ID:	129285
UniProt ID:	Q6ZMI0
Cytogenetics:	2p16.3
Gene Summary:	Putative regulator of protein phosphatase 1 (PP1) activity (PubMed:19389623). May play a role in the endosomal sorting process or in endosome maturation pathway (PubMed:30520571).[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) encodes the longest protein (isoform 1).