

Product datasheet for **SC325329**

PTPH1 (PTPN3) (NM_001145369) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PTPH1 (PTPN3) (NM_001145369) Human Untagged Clone
Tag:	Tag Free
Symbol:	PTPN3
Synonyms:	PTP-H1; PTPH1
Vector:	<u>pCMV6 series</u>



[View online »](#)

Fully Sequenced ORF: >NCBI ORF sequence for NM_001145369, the custom clone sequence may differ by one or more nucleotides

```

ATGGATATTTGCGAAGGAAGGTTAACCTGCCCTCTTAACTCAGCAGTGGTCTAGCGTCC
TATGCCGTACAATCTCATTTTGGAGACTATAATTCTTCCATACATCATCCAGGCTATCTT
TCCGATAGTCACTTTATACCCGATCAAAATGAGGACTTTTTAACAAAAGTCAATCTCTG
CATGAGCAGCACAGTGGGCTAAAAAATCAGAAGCAGAATCCTGCTATATCAACATAGCG
CGGACCCCTCGACTTCTATGGAGTAGAACTGCACAGTGGTAGGGATCTGCACAATTTAGAC
CTAATGATTGGAATTGCTTCCGCGGGTGTGCTGTGTACCGAAAAATACATTTGCACAAAGT
TTCTATCCTTGGGTGAACATTCTCAAAATTTCTTTCAAAAAGAAAAAGTTCTTCATACAT
CAGCGACAGAAACAGGCTGAATCCAGGGAACATATTGTGGCCTTCAACATGCTGAATTAC
CGATCTTGCAAAAACCTTGGAATCCTGTGTTGAGCACCATACGTTCTTTTCAGGCAAAG
AAGCTACTACCTCAGGAAAAAGATGTTCTGTCTCAGTACTGGACTATGGGCTCTCGGAAC
ACCAAAAAGTCGGTAAATAACCAATATTGCAAAAAGGTGATTGGCGGGATGGTGTGGAAC
CCAGCCATGCGGAGATCCTTATCAGTGGAGCACTTAGAAACCAAGAGTCTGCCTTCTCGT
TCCCCTCCATTACTCCCAACTGGCGAAGTCTCGGCTCCGGCAGAAATCCGAAAGCCA
CGCCACTCTTCTGCAGATAACCTTGCAAATGAAATGACCTACATCACGGAACCGGAAGAT
GTATTTTACACGTACAAGGGCTCTCTGGCCCTCAAGACAGCGATTCTGAAGTTTCTCAG
AACCGAAGCCCGCACCAAGAGAGTTTATCCGAGAACAATCCGGCACAAAGCTACCTGACC
CAGAAGTCATCCAGTTCTGTGTCTCCATCTTCAAATGCTCCAGGCTCCTGCTCACCTGAC
GGCGTTGATCAGCAGCTCTTAGATGACTCCACAGGGTGACCAAAGGGGGCTCCACCGAG
GACGCCAGCCAGTACTACTGTGACAAGAATGATAATGGTGACAGTACTTAGTCTTGATC
CGTATCACACCAGATGAAGATGGAAAAATTTGGATTTAATCTTAAGGGAGGAGTGGATCAA
AAGATGCCTCTTGTGGTATCAAGGATAAACCAGAGTCACTCGGGACACCTGCATTCCT
AAGCTGAACGAAGGGGATCAAATCGTGTTAATCAATGGCCGGGACATCTCAGAACACAGC
CATGACCAAGTGGTGTGTTTATCAAAGCCAGCCGGGAGTCCCACTCACGGGAGCTGGCC
CTGGTGATCAGGAGGAGAGCTGTCCGCTCATTGCTGACTTCAAGTCTGAAGATGAACTG
AACCGACTTTTCCCCGAAGCCATTTTCCCATGTGTCCGGAGGGTGGGGACACTTTGGAG
GGATCCATGGCACAGCTAAAGAAGGGCCTCGAAAGCGGGACGGTGTGATCCAGTTTGAG
CAACTCTACAGAAAAAGCCAGGTTTGGCCATCACGTTTGCAAAGCTGCCTCAAAATTTG
GACAAAAACCGATATAAAGATGTGCTGCCTTATGACACCACCCGGGATTATTGCAGGGA
AATGAAGATTATTAATGCAAGTTACGTGAACATGGAATTCCTGCTGCTAACCTGTG
AACCAAGTACATCGCCACTCAGGGCCCTGCCGATACCTGTGCACAGTTTTGGCAGGTT
GTCTGGGATCAGAAGTTGCACTCATTGTGATGTTGACGACTCTCACAGAACGAGGGCGG
ACCAAAATGTCACCACTACTGGCCAGATCCCCCGACGTCATGAACCACGGCGGCTTTTAC
ATCCAGTGTGACTCAGAGGACTGCACCATCGCCTATGTGTCCCGAGAAATGCTGGTCACA
AACACCCAGACCGGGGAAGAACACACAGTGACACATCTCCAGTACGTCCGATGGCCTGAC
CACGGTGTGCCGATGACTCCTCCGACTTTCTGGAATTTGTAACATGTGAGGTCTCTG
AGAGTGGACAGCGAGCCCGTCTAGTTCACTGCAGTGTGGAATAGTTCGAACCGGTGTG
TTGGTCACTATGGAACAGCCATGTGCCTAACTGAGAGGAACCTGCCATTTACCCACTG
GATATTGTCCGAAAAATGCGAGACCAGCGCCATGATGGTGCAGACATCAAGCCAGTAC
AAGTTTGTGTGGAAGCGATTCTTCGTGTGATGAAGAAGGTTTGTCCAAATGCTGGAT
CCTAGT

```

Restriction Sites: Please inquire

ACCN: NM_001145369

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001145369.1 , NP_001138841.1
RefSeq Size:	8783 bp
RefSeq ORF:	2349 bp
Locus ID:	5774
UniProt ID:	P26045
Cytogenetics:	9q31.3
Protein Families:	Druggable Genome, Phosphatase
Gene Summary:	<p>The protein encoded by this gene is a member of the protein tyrosine phosphatase (PTP) family. PTPs are known to be signaling molecules that regulate a variety of cellular processes including cell growth, differentiation, mitotic cycle, and oncogenic transformation. This protein contains a C-terminal PTP domain and an N-terminal domain homologous to the band 4.1 superfamily of cytoskeletal-associated proteins. P97, a cell cycle regulator involved in a variety of membrane related functions, has been shown to be a substrate of this PTP. This PTP was also found to interact with, and be regulated by adaptor protein 14-3-3 beta. Several alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Feb 2009]</p> <p>Transcript Variant: This variant (3) differs in the 5' UTR, and initiates translation from an in-frame downstream AUG, resulting in an isoform (3) with a shorter N-terminus compared to isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>