

Product datasheet for **SC325330**

SUN1 (NM_001130965) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SUN1 (NM_001130965) Human Untagged Clone
Tag:	Tag Free
Symbol:	SUN1
Synonyms:	UNC84A
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC325330 representing NM_001130965.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGATTTTTCTCGGCTTCACATGTACAGTCTCCCCAGTGTGTGCCGGAGAACACGGGCTACACGTAT
GCGCTCAGTTCAGCTATTCTTCCAGATGCTCTGGATTTTGAGACGGAGCACAAATTGGACCTGTATTT
GATTCTCCACGGATGTCCCGCCGTAGTTTGCGCCTGGCCACGACAGCATGCACCCTGGGGATGGTGAG
GCTGTGGGTGCCGACAGCGGCACCAGCAGCGCTGTCTCCCTGAAGAACCAGCGGCCAGAACAACAAAA
CAGCGCAGAAGCACAAACAAATCAGCTTTTAGTATCAACCACGTGTCAAGGCAGGTACAGTCTCTGGC
GTCAGCCACGGCGGCACTGTGAGCTGCAGGATGCTGTGACTCGACGGCCTCCTGTATTGGACGAGTCT
TGGATTCGTGAACAGACCACAGTGGACCACTTCTGGGTCTTGATGATGATGGTATCTTAAAGTGGA
AATAAAGTGCCTTTCAGGAAACGGGGATGTGGGAGCCGCCGCCACCGCGCACACGGCTTCTCC
TGCAGCAACTGCAGCATGCTGTCCGAGCGCAAGGACGTGCTCACGGCGCACCCCGCGCCCCGGGCC
GTGTCGAGAGTTTATTCTAGGGACAGGAATCAAAAATGTTACTTCTTGTGCAGATTCTGCGCAGGATC
GGAGCTGTGGCCAGGCTGTGTCCAGGACGGCGTGGTCCGCCCTTTGGCTGGCCGTGGTTGCTCCAGGG
AAGGCAGCCTCTGGAGTGTCTGGTGGCTGGGGATTGGATGGTACCAGTTTGTACTTTGATTTCTTGG
CTGAATGTGTTTCTTACCAGGTGCCTTCGAAACATCTGCAAGTTTTAGTCTTGCTCATCCCACTC
TTCCTTTTACTAGCAGTCTCTCCTTACGGGGCCAGGGCAATTTCTTTTCGTTCTTGCCCGTGTGAAC
TGGGAAGCATGCATAGAACACAGCGGGTGGATGACCCCAAGGACGTGTTTAAACCCACGACTTCTCGC
CTGAAGCAGCCTTGCAGGGTGCAGTGAAGCTTTTCCGTGGCATTGGATGAGTGGCGTGGAGCAGCAG
TAGGCCCTCTGTCTGGACAGTCCACCACCATGGTGAAGTCTCCGAGAGCTGACCACTTGTCTACAG
AAGTGCAGGCTCGGGTGGACCAGATGGAAGCGCGCTGCCGGCCGTGAGTTCGGTCCAGAGACGCT
GTGGGACAGCCCCGAGGAGACTGACTTTATGGCCTTTCACCAAGAACATGAAGTGCATGTCCACAC
TTGAAGATATTCTGGGAAAACAGAGAAAAATCTGAGGCCATCCAGAAGGAACTAGAACAGACCAAG
CAAAAAACAATCAGTGCAGTTGGTGAAGCTCCTGCCACAGTGCAGCACCTCCAGCTGGAGCTGGAT
CAGCTAAAGTGCAGCTGTCCAGCTGGCGACACGTGAAGACCGGTGTGAGACAGTGGATGCCGTACAA
GAAAGAGTGGACGTGCAAGTGCAGAGAAATGGTGAAGTCTGTTTTCCGAAGATCAGCAAGGCGTTCT
CTGGAACAGCTGCTGCAGAGTTCTCATCACAGTTTGTGAGCAAAGGCGACTTGCAGACGATGCTGCGA
GACCTGCAGCTGCAGATCCTGCCGAACGTACCCACCACGTTTCCGTGACCAAGCAGCTCCCAACCTCA
GAAGCCGTGGTGTCTGCTGTGAGCGAGGCGGGGGCGTCTGGAATAACAGAGGCGCAAGCAGTGCATC
GTGAACAGCGCCTTGAAGCTGTATCCCAAGATAAGACCGGGATGGTGGACTTTGCTCTGGAATCTGGT
GGTGGCAGCATCTTGAGTACTCGCTGTTCTGAAACTTACGAAACCAAAACGGCGCTGATGAGTCTGTTT
GGGATCCCGCTGTGGTACTTCTCGCAGTCCCGCGCGTGGTATCCAGCCTGACATTTACCCCGGTAAC
TGCTGGGCATTTAAAGGCTCCCAGGGTACCTGGTGGTGAAGGCTCTCCATGATGATCCACCCAGCCGCC
TTCCTCTGGAGCACATCCCTAAGACGCTGTCCCAACAGGCAACATCAGCAGCGCCCCAAGGACTTC
GCCGTCTATGGATTAGAAAATGAGTATCAGGAAGAAGGGCAGTCTCTGGACAGTTCACGTATGATCAG
GATGGGAGTGCCTCCAGATGTTCCAGGCCCTGAAAAGACCCGACGACACAGCTTCCAAATAGTGGAA
CTTCGGATTTTTTCTAACTGGGCCATCCTGAGTATACCTGTCTGTATCGGTTTCAGAGTTCATGGCGAA
CCTGTCAAGTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001130965

Insert Size: 2358 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001130965.2](#)

RefSeq Size: 3977 bp

RefSeq ORF: 2358 bp

Locus ID: 23353

UniProt ID: [O94901](#)

Cytogenetics: 7p22.3

Protein Families: Transmembrane

MW: 87.1 kDa

Gene Summary: This gene is a member of the unc-84 homolog family and encodes a nuclear envelope protein with an Unc84 (SUN) domain. The protein is involved in nuclear anchorage and migration. Alternatively spliced transcript variants have been described. [provided by RefSeq, Jan 2019]
Transcript Variant: This variant (1) encodes isoform a. Variants 1, 72, and 73 all encode the same isoform (a).