

Product datasheet for **SC325343**

ADAM29 (NM_001130705) Human Untagged Clone

Product data:

| | |
|---------------------------|--|
| Product Type: | Expression Plasmids |
| Product Name: | ADAM29 (NM_001130705) Human Untagged Clone |
| Tag: | Tag Free |
| Symbol: | ADAM29 |
| Synonyms: | CT73; svph1 |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pCMV6-Entry (PS100001) |
| E. coli Selection: | Kanamycin (25 ug/mL) |



[View online »](#)


Fully Sequenced ORF: >SC325343 representing NM_001130705.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGAAGATGTTACTCTGCTGCATTGCCTTGGGGTGTTCCTGTCCTGTTCTGGACACATCCAGGATGAG
CACCCCAATATCACAGCCCTCCGGATGTGGTGATTCTGTGAGGATAACTGGCACCACCAGAGGCATG
ACACCTCCAGGCTGGCTCCTATATCCTGCCCTTTGGAGGCCAGAAACACATTATCCACATAAAGGTC
AAGAAGCTTTTGTTCCAAACCTCCCTGTGTTACCTACACAGACCAGGGTGTATCCTTGAGGAC
CAGCCATTTGCCAGAATAACTGCTACTATCATGGTTATGTGGAAGGGGACCAGAATCCCTGGTTTCC
CTCAGTACCTGTTTTGGGGTTTTCAAGGAATATTACAGATAAATGACTTTGCTTATGAAATCAAGCCC
CTAGCATTTTCTACCACGTTTGAACATCTGGTATAACAAGATGGACAGTGAGGAGAAAACATTTTCAACC
ATGAGATCCGGATTTATGCAAAATGAAATAACATGCCAATGGAATTTGAAGAAATTGATAATCCACT
CAGAAGCAAAGTCTTATGTGGGCTGGTGGATCCATTTTAGGATTGTTGAAATTGTAGTCGTCATTGAT
AATTATCTGTACATTCGTTATGAAAGGAACGACTCAAAGTTGCTGGAGGATCTATATGTTATTGTTAAT
ATAGTGGATTCCATTTTGGATGTCATTGGTGTTAAGGTGTTATTATTTGGTTTGGAGATCTGGACCAAT
AAAAACCTCATTGTAGTAGATGATGAAGGAAATCTGTGCACCTGTATTGCAAGTGAAGTCGGAGAAC
ATTACGCCCGGATGCAACATGACACCTCACATCTTTTCACAACTCTAGGATTAAGAGGGTTAAGTGGC
ATAGGAGCTTTTAGAGGAATGTGTACACCACACCGTAGTTGTGCAATTGTTACTTTCATGAACAAAAC
TTGGGCACTTTTTCAATTGCAGTGGCTCATCATCTAGGTCATAATTTGGGCATGAACCATGATGAGGAT
ACATGTCGTTGTTCAACCTAGATGCATAATGCATGAAGGCAACCCACCAATAACTAAATTTAGCAAT
TGTAAGTTATGGTATTTTGGGAATATACTGTAGAGAGGACAAAGTGTGTTGCTTGAACAGTACACACA
AAGGACATCTTTAATGTGAAGCGCTGTGGGAATGGTGTGTTGTTGAAGAAGGAGAAGAGTGTGACTGTGA
CCTTTAAAGCATTGTGCAAAAGATCCCTGCTGTCTGCAAAATGCACTCTGACTGATGGTTCTACTTGT
GCTTTTGGGCTTTGTTGCAAAGACTGCAAGTTCCTACCATCAGGGAAAAGTGTGTAGAAAGGAGGTCAAT
GAATGTGATCTTCCAGAGTGGTGAATGGTACTTCCCATAAGTGCCAGATGACTTTTATGTGGAAGAT
GGAATTCCTGTAAGGAGAGGGGCTACTGCTATGAAAAGAGCTGTGATGACCGCAATGAACAGTGTAGG
AGGATTTTGGTGCAGGCGCAAATACTGCAAGTGAGACTTGCTACAAAGAATTGAACACCTTAGGTGAC
CGTGTGGTCACTGTGGTATCAAAAATGCTACATATAAAGTGAATATCTCAGATGTCCAGTGTGGA
AGAATTCAGTGTGAGAATGTGACAGAAATCCCAATATGAGTGATCATACTACTGTGATTGGGCTCGC
TTCAATGACATAATGTGCTGGAGTACTGATTACCATTTGGGGATGAAGGGACCTGATTTGGTGAAGTG
AAAGATGGAACAGAGTGTGGGATAGATCATATATGCATCCACAGGCACTGTGTCATATAACCATCTTG
AATAGTAATTGCTCACCTGCATTTTGAACAAGAGGGGCATCTGCAACAATAAACATCACTGCCATTGC
AATTATCTGTGGGACCCTCCCAACTGCCTGATAAAAGGCTATGGAGGTAGTGTGACAGTGGCCACCC
CCTAAGAGAAAGAAGAAAAAGAAGTTCTGTTATCTGTGTATATTGTTGCTTATTGTTTTGTTATTTTA
TTATGTTGCTTTTATCGACTTTGTAAAAAAGTAAACCAATAAAAAAGCAGCAAGATGTTCAAATCCA
TCTGCAAAAGAAGAGGAAAAAATTCAGCGTCGACCTCATGAGTTACCTCCCAGAGTCAACCTTGGGTG
ATGCCTCCCAGAGTCAACCTCCTGTGACGCCTCCCAGAGTCATCCTCAGGTGATGCCTCCCAGAGT
CAACCTCCTGTGACACCCTCCCAGAGTCAACCTCGGGTGTGCTTCTCAGAGTCAACCTCCTGTGATG
CCTTCCCAGAGTCATCCTCAGTTGACGCCTTCCCAGAGTCAACCTCCTGTGACACCCTCCCAGAGGCAA
CCTCAGTTGATGCCTTCCCAGAGTCAACCTCCTGTGACGCCTCCTAG
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001130705

Insert Size: 2463 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001130705.1](#)

RefSeq Size: 3138 bp

RefSeq ORF: 2463 bp

Locus ID: 11086

UniProt ID: [Q9UKF5](#)

Cytogenetics: 4q34.1

Protein Families: Druggable Genome, Transmembrane

MW: 92.8 kDa

Gene Summary: This gene encodes a member of the ADAM (a disintegrin and metalloprotease domain) family. Members of this family are membrane-anchored proteins structurally related to snake venom disintegrins, and have been implicated in a variety of biological processes involving cell-cell and cell-matrix interactions, including fertilization, muscle development, and neurogenesis. The protein encoded by this gene is highly expressed in testis and may be involved in human spermatogenesis. Alternative splicing results in multiple transcript variants that encode the same protein. [provided by RefSeq, Jul 2008]
Transcript Variant: This variant (4) differs in the 5' UTR, compared to variant 1. Variants 1, 2, 3, 4, 5, 6, and 7 encode the same protein.