

Product datasheet for **SC325344**

iASPP (PPP1R13L) (NM_001142502) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	iASPP (PPP1R13L) (NM_001142502) Human Untagged Clone
Tag:	Tag Free
Symbol:	PPP1R13L
Synonyms:	IASPP; NKIP1; RAI; RAI4
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC325344 representing NM_001142502.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGACAGCGAGGCATTCAGAGCGCGCGGACTTTCTGGACATGAACCTCCAGTCGCTGGCCATGAAA
CACATGGATCTGAAGCAGATGGAGCTGGACACGGCGCGGCCAAGGTGGATGAACCTGACCAAGCAGCTG
GAGTCGCTGTGGTCAGACTCTCCCGCGCCTCTGGCCCGCAGGCCGGACCCCTTCTAGGCCGCCCGG
TACAGCTCCAGCTCGATCCCTGAGCCCTTCGGCAGCCGAGGGTCCCCCGGAAGCGGCCACCGACGGC
GCAGACACCCCGTTCCGACGATCAGAGAGTGCCCCAACCCCTACACCCCTACAGCCCGCTGTCCCCAAG
GGACGGCCGTCGTCGCCGCGACCCCGCTACCTGCAGCCGGACGCCTACGGCAGCCTGGACCGCGCG
ACCTCGCCCCGGCCCCGCGCCTTCGATGGCGCAGGCAGCTCCCTCGGCCGTGCGCCCTCCCCGGGCC
GGGCCAGGCCCGCTCCGCCAGCAGGGTCCCCCACGCCTTCGACTTCTGGGCCGCGCAGGCTCCCC
CGCGGCAGCCCCCTGGCGAGGGGCCCCAGGCCTTCTCCCGAGCGTGGGCCGTACCGCGCCCCCT
GCCACAGCCTACGACGCCCAGCGTCCGCCTTCGGGAGCTCCCTGCTAGGCTCCGGCGGCAGCGCATT
GCCCGCCTCTGCGCGCGCAAGACGACCTGACGCTGCGCCGGCGGCCCTCGAAAGCTGGAACGAGTCT
GACCTGGACGTGGCGTACGAGAAGAAGCCTTCGACAGACGCGAGCTATGAACGCCTGGACGTCTTCGCA
AGGCCTGCCTCGCCGAGCCTGCAGCTGTTGCCTTGAGGGGAGAGCAGCCTGGATGGACTGGGGGGCACC
GGCAAGGACAACCTCACTAGCGCCACCCTGCCGCGCAATTACAAGGTCTCTCTCTGGCCAGCGACCCG
CGTTCAGACGCGGGCAGCTACCGCGCTCGCTGGGCTCCGCGGGGCCGTGGGCACTTTGCCTCGCAGC
TGGCAGCCCGTCAGCCGATCCCATGCCCCCTCCAGCCCCAGCCCCGCGGGGCCCGCGCCAGCGT
CCCATCCCCCTCAGCATGATCTCAAGCTGCAGAACGCCTTCTGGGAGCACGGGGCCAGCCGCGCATG
CTCCCTGGGTCCCCCTTTACCCGAGCACCCCGCCTAAGCTGCAGCCCCAACCAACACAGCCCC
CAGCCACAATCACAACCACAGCCCCAGCTGCCCCACAGCCCCAGACCAACCCCAACCCCTACCCCA
GCCCCACAGCATCCCAACAGACATGGCCCTGTGAACGAAGGACCCCAACCCCAACCCGAGCTG
GAGCCTGAGCCGGAGATAGAGGGGTGCTGACACCAGTGTGGAGGCTGGCGATGTGGATGAAGGCCT
GTAGCAAGGCCTCTCAGCCCCACGAGGCTGCAGCCAGCACTGCCACCGGAGGCACAGTCGGTGCCCGAG
CTGGAGGAGGTGGCACGGGTGTTGGCGAAATTCGCCGCCCTCAAACGCAGGGGCTCCATGGAGCAG
GCCCTGCTGTGGCCCTGCCCTACCCACAAGAAACAGTACCAGCAGATCATCAGCCGCTCTTCCAT
CGTCATGGGGGGCAGGGCCCGGGGGCCGGAGCCAGAGCTGTCCCATCACTAGGGATCTGAGGCC
AGGGCAGGGCCCCCTGCTCCTGCCCCACAGCTCCCATTCACCCCGGCCCGTCCAGAGCAGCCCA
CCAGAGCAGCCGAGAGCATGGAGATGCCTCTGTGCTGCGGAAGGCGGGCTCCCGCGCAAGGCCCGC
CGCGCGCCCTCAACCCTCTGGTCTCCTCTGGACGCGCGCTGACCGGGAGCTGGAGGTGGTGCAG
CAGGCGGTGAAGGAGATGAACGACCCGAGCCAGCCCAACGAGGAGGGCATCACTGCCTTGCAACAGCC
ATCTGGCGGCCAACTACTCTATCGTGGATTTCTCATACCGCGGGTGCCAATGTCAACTCCCCGAC
AGCCACGGCTGGACACCTTGCACTGCGCGCGCTCGTGAACGACACAGTCACTGCATGGCGCTGGTG
CAGCACGGCGTGCAATCTCGCCACCAGCTCAGCGACGGCGCCACCGCCTTCGAGAAGTGCAGCCCT
TACCGCGAGGGTTATGCTGACTGCGCCACCTACCTGGCAGAGCTCGAGCAGAGTATGGGGCTGATGAAC
AGCGGGCAGTGTACGCTCTCTGGACTACAGCGCCGAGTTCGGGGACGAGCTGTCTTCCGCGAGGGC
GAGTCGGTACCGTGTGTCGAGGAGGACGGGCCGAGGAGACCGACTGGTGGTGGCCGCGCTGCACGGC
CAGGAGGGCTACGTGCCCGGAATACTTCGGGCTGTTCCCGAGGTGAAGCCTCAAAGGAGTAAAGTC
TAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTAAACGGCCGGC
  
```

Restriction Sites: Sgfl-Mlul

Plasmid Map:


ACCN: NM_001142502

Insert Size: 2487 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq:	<u>NM_001142502.1</u>
RefSeq Size:	3122 bp
RefSeq ORF:	2487 bp
Locus ID:	10848
UniProt ID:	<u>Q8WUF5</u>
Cytogenetics:	19q13.32
Protein Families:	Transcription Factors
MW:	89.1 kDa
Gene Summary:	<p>IASPP is one of the most evolutionarily conserved inhibitors of p53 (TP53; MIM 191170), whereas ASPP1 (MIM 606455) and ASPP2 (MIM 602143) are activators of p53.[supplied by OMIM, Mar 2008]</p> <p>Transcript Variant: This variant (1) represents the longer transcript. Variants 1 and 2 encode the same protein.</p>