

Product datasheet for **SC325346**

ZNF438 (NM_001143767) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF438 (NM_001143767) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF438
Synonyms:	bA330O11.1
Vector:	<u>pCMV6 series</u>



[View online »](#)

Fully Sequenced ORF: >NCBI ORF sequence for NM_001143767, the custom clone sequence may differ by one or more nucleotides

```

ATGCAGAACTCTGTATCAGTACCACCAAAAGATGAAGGTGAATCAAACATCCCTTCTGGA
ACAATACAGAGTAGGAAAGGTTTGCAGAATAAGAGTCAGTTTAGGACCATTGCACCAAAA
ATTGTGCCCAAAGTCCTAACGTCCAGAATGCTGCCATGTCATTACCATCACGCTCTGAT
CAGGTGAATCTGGGACCTCCATCAACTCCAAGCTGCTGGGGATGTCCACCCAGAAGTAT
GCCCTGATGCAGGTTGCTGGCCAGGAGGGGACATTTTCTTTGTTGCTCTGCCACATGTT
GCCTCAGCTCAGCCAATTACAGAAACCCAGAATGTCCCTACCTGAAAACCTGAAACTTCCT
ATTCTAGATATCAACCCCTAGAAATAGCAAAGCATCAAGAAAGAAACCATCTCTGATC
TTTCTAAGAGTGGCTGTAGCAAAGCTCCTGCCAAACCCAAATGTGTCCTCAGATGTCC
CCTTCCCCACCTCACCACCTGAACTCCTGTACAAACCCAGTCCATTTGAGGAAGTACCA
TCACTAGAGCAAGCCCCAGCCAGCATTAGCACAGCTGCGTGACCAATGGAAGTGACCAT
GGGGACTTGAGACCACAGTGACCAACACCCATGGCAGTCTGAACCTCCTGCTACCCCA
GCATCATCCACACCAGAGGAGCCTGCCAAGCAGGACCTCACAGCTCTTTCAGGGAAAGCA
CACTTTGTAAGCAAGATAACATCTAGTAAACCTTCTGCTGTTGCCAGTGAAAAATTTAA
GAACAAGTTGATCTTGCAAAAACCATGACCAATTTATCACAACCATTTCTGGCAATGCA
GTTCAAGTTGATCTCTTCACTCCCCAAAGGAAACTGCCAATCCCACCTACTCAAGAATG
AAGACAATGGAGGTTTACAAAATCAAATCAGATGCTAACATTGCAGGTTTTTCTTTACCA
GGACCTAAGGCCGACTGTGATAAGATACCTCCACCACAGAAGGCTTTAATGCAGCCACC
AAGGTGGCAAGCAGGCTACCTGTTCCACAAGTGTACAGCAGAGTGCCTGTGAAAGTGCC
TTTTGTCCACCCACCAAACTTGATCTTAACCACAAAACAAACTGAACAGTGGAGCAGCA
AAGAGAAAAGGAAGAAAACGGAAGGTACCAGATGAAATTTTGGCATTTCAGGGAAAAAGG
AGGAAATATATCATTAAATAGTGTAGAGATGGTAAAGAAAGAGTAAAAATGATCCCCAA
GAATTCAGAGACCAAAAGCTGGGGACCTGAAAAAATACCGTAGCATTATGCCCAAACCT
ATCATGGTCATACCCACTTTGGCCTCCCTGGCTTCTCCTCACTACACTACAGTCCCAGATG
CTTGGGGGCTAGGACAGGATGTTTTGTTAAATAATTCACTCACTCTAAATATCTTGGC
TGTAAGCAAGACAACAGCTCTTCCCTAAGCCCAGCTCCGTGTTCAAGAAATGGATTCTCT
GGCATTAAAGAGCCTTGGCACAGATGTCACGTCTGCAACCACCACTTCCAGTTCAAACAG
CACCTTCGAGACCACATGAATACACACACCAACAGACGCCCTTACAGTTGTCGGATTGT
CGCAAGTCTATGTACGCTCTGGCAGCCTGAGCACACACATGAACTTCATCATGGTGAG
AACCGTCTGAAGAACTCATGTGTTGTGAGTTTTGTGCAAAAGTGTGTTGGCCACATCCGA
GTCTATTTTGGCCATCTGAAAGAAGTGATAGGGTTGTGATCAGCACTGAGCCTGCGCCC
AGTGAAGTGCAGCCAGGAGACATACCAAAGAACAGAGACATGAGTGTGCGAGGCATGGAG
GGATCATTGGAGAGGGAAAAACAAGTCCAACCTGGAAGAAGACTTCCTTCTAAACAGGCA
GACGAAGTCAAATTACAAATCAAATGTGGTCGTTGTGAGATTACTGCTCAGTCTTTTGGC
GAAATAAAATTTCACTTACTTGATGTTTATGGAGAGGAAATGAGGGCAGGCTACAAGAA
GGGACCTTCCCAGGAAGCAAGGGGACTCAGGAAGAGTTGGTGCAGCACGCTAGCCCCGAC
TGAAAAAGGCATCCTGAGAGAGGGAAGCCGGAAGGTTTCACTCCTCCTCGAGGAATCA
CATGCATGTCCAAGACTGAAAAGGCAGCTCCACCTTCATCAGAATGGCGTGGAAATGCTC
ATGGAAAAATGAAGACCCAGTCAGGAACCAACAAGCCAAGGGAACCTGCCAGGGCCCT
GAGTGTCTGGCCTCCACACGTTTCTTGTGGTCCCATTCAAGGCTTTAACTGCCTGCTT
TGTGCAGAGATGCTGGGACGGAAGAGGACCTCCTCCACCACTGGAAGCACCAGCATAAC
TGTGAGGACCTTCCAACTGTGGGCTATTTTAAATACGGTCTCCAACAGGGAGTGATC
GAACTTTCCAGTGAAGCTGAGAAA

```

Restriction Sites: Please inquire
ACCN: NM_001143767

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001143767.1, NP_001137239.1</u>
RefSeq Size:	3119 bp
RefSeq ORF:	2487 bp
Locus ID:	220929
UniProt ID:	<u>Q7Z4V0</u>
Cytogenetics:	10p11.23
Gene Summary:	Isoform 1 acts as a transcriptional repressor.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (3) differs in the 5' UTR, compared to variant 1. Variants 1-4 encode the same isoform (a).