

Product datasheet for **SC325498**

CYB561D1 (NM_001134403) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CYB561D1 (NM_001134403) Human Untagged Clone
Tag:	Tag Free
Symbol:	CYB561D1
Vector:	<u>pCMV6 series</u>
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001134403, the custom clone sequence may differ by one or more nucleotides ATGCAGCCCTGGAGGTAGGTCTGGTCCCGCTCCAGCTGGGGAGCCGAGACTGACCCGC TGGCTGCGGAGAGGCAGTGGGATCTTGGCGCACCTGGTAGCTTTGGGCTTACCATCTTT CTGACAGCGCTGTCCCGCCAGGAACCAAAACAGGTCCCCTGATGGAGGATAGAAGTGAA GGAGGCCGGGCGCGGTGGGTCATGCCTGAAATCCCAGCACTTTGGGAGGCCGACGCGGGT GGATCACTTGAGGTCTTTCTCCTGGCACCCCTGTATTCATGGCCTTGGCGTTCTGCCTCT GCATGGCTGAAGCCATCTACTTCTCACCTGAACACTCCCTGTTCTTCTTCTGCTCCC GAAAAGCACGGATCCGGTCCACTGGGCAGGGCAGACCC
Restriction Sites:	Please inquire
ACCN:	NM_001134403
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.



[View online »](#)

RefSeq:	<u>NM_001134403.1, NP_001127875.1</u>
RefSeq Size:	5048 bp
RefSeq ORF:	402 bp
Locus ID:	284613
Cytogenetics:	1p13.3
Protein Families:	Transmembrane
Gene Summary:	<p>Probable transmembrane reductase that may use ascorbate as an electron donor and transfer electrons across membranes to reduce monodehydro-L-ascorbate radical and iron cations Fe(3+) in another cellular compartment.[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (4) contains an alternate coding region segment, which changes the reading frame, compared to variant 1. The resulting protein (isoform 4) has a shorter and distinct C-terminus when it is compared to isoform 1. Sequence Note: The RefSeq transcript and protein were derived from transcript and genomic sequence to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on alignments.</p>