

Product datasheet for **SC325515**

FKBP11 (NM_001143782) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	FKBP11 (NM_001143782) Human Untagged Clone
Tag:	Tag Free
Symbol:	FKBP11
Synonyms:	FKBP19
Vector:	<u>pCMV6 series</u>
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001143782, the custom clone sequence may differ by one or more nucleotides <pre> ATGACCCTGCGCCCTCACTCCTCCCGCTCCATCTGCTGCTGCTGCTGCTCAGTGCG GCGGTGTGCCGGGCTGAGGCTGGGCTCGAAACCGAAAGTCCCGTCCGACCCTCCAAGT GAGACCCTGGTGGAGCCCCAGAACCATGTGCCGAGCCCGCTGCTTTGGAGACACGCTT CACATACACTACACGGGAAGCTTGGTAGATGGACGTATTATTGACACCTCCCTGACCAGA GACCTCTGGTTATAGAACTTGGCCAAAAGCAGGTGATTCCAGGTCTGGAGCAGAGTCTT CTCGACATGTGTGTGGGAGAGAAGCGAAGGGCAATCATTCTTCTCACTTGGCCTATGGA AAACGGGGATTTCACCATCTGTCCCAGGGACTAAAGACAACCTGATGAGGCCACCTGGC ATGACCTCCAGCAGCCAG </pre>
Restriction Sites:	Please inquire
ACCN:	NM_001143782
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).


[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: NM_001143782.1, NP_001137254.1

RefSeq Size: 986 bp

RefSeq ORF: 441 bp

Locus ID: 51303

UniProt ID: Q9NYL4

Cytogenetics: 12q13.12

Protein Families: Druggable Genome, Transmembrane

Gene Summary: FKBP11 belongs to the FKBP family of peptidyl-prolyl cis/trans isomerases, which catalyze the folding of proline-containing polypeptides. The peptidyl-prolyl isomerase activity of FKBP proteins is inhibited by the immunosuppressant compounds FK506 and rapamycin (Rulten et al., 2006 [PubMed 16596453]).[supplied by OMIM, Mar 2008]
 Transcript Variant: This variant (3) differs in the 3' UTR and coding sequence compared to variant 1. The resulting isoform (3) has a shorter and distinct C-terminus compared to isoform 1.