

Product datasheet for SC325543

PPIAL4C (NM_001135789) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PPIAL4C (NM_001135789) Human Untagged Clone
Tag:	Tag Free
Symbol:	PPIAL4C
Synonyms:	COAS-2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC325543 representing NM_001135789. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGTCAACTCCGTCGTCTTTTTGACATCACCGTCGACGGCAAGCCCTTGGGCCGATCTCCATCAAA
CTGTTTGACAGACAAGATTCCAAAGACAGCAGAAAACCTTCGTGCTCTGAGCACTGGAGAGAAAGGATTT
CGTTATAAGGGTTCCTGCTTTCACAGAATTATTCCAGGGTTTATGTGTCAGGGTGGTGACTTCACACGC
CCTAATGGCACCAGTCCATCTATGGGGAGAAATTTGATGATGAGAACCTCATCCGAAAGCAT
ACAGGTTCTGGCATCTTGCCATGGCAAATGCTGGACCAACACAAATGGTCCAGTTTTTCATCTGC
ACTGCAAGACTGAGTGGTTGGATGGCAAGCATGTGGCCTTTGGCAAGGTGAAAGAACGTGTGAATATT
GTGGAAGCCATGGAGCACTTTGGGTACAGGAATAGCAAGACCAGCAAGAAGATCACCATTGCTGACTGT
GGACAATTCTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```

Restriction Sites: Sgfl-Mlul



[View online »](#)

Plasmid Map:



ACCN: NM_001135789

Insert Size: 495 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001135789.3](#)

RefSeq Size: 785 bp

RefSeq ORF: 495 bp

Locus ID: 653598

UniProt ID: [A0A0B4J2A2](#)

Cytogenetics: 1q21.2

MW: 18.2 kDa

Gene Summary: PPIases accelerate the folding of proteins. It catalyzes the cis-trans isomerization of proline imidic peptide bonds in oligopeptides (By similarity).[UniProtKB/Swiss-Prot Function]