

Product datasheet for **SC325628**

FGF13 (NM_001139502) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	FGF13 (NM_001139502) Human Untagged Clone
Tag:	Tag Free
Symbol:	FGF13
Synonyms:	DEE90; FGF-13; FGF2; FHF-2; FHF2; LINC00889
Vector:	<u>pCMV6 series</u>
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001139502, the custom clone sequence may differ by one or more nucleotides ATGTTACGACAAGATTCCATCCAATCTGCGGAATTAAGAAAAAAGAGTCCCCCTTCGT GCTAAGTGTACGAAATCTTCTGCTGCCGCTGAAGCAAGTACACCACAAAGAGAACACA GAGCCGGAAGAGCCTCAGCTTAAGGGTATAGTTACCAAGCTATACAGCCGACAAGGCTAC CACTTGCAGCTGCAGGCGGATGGAACCATTGATGGCACCAAGATGAGGACAGCACTTAC ACTCTGTTTAACTCATCCCTGTGGTCTGCGAGTGGTGGCTATCCAAGGAGTTCAAACC AAGCTGTACTTGGCAATGAACAGTGAGGGATACTGTACACCTCGGAACCTTTTCACACCT GAGTGCAAATTCAAAGAATCAGTGTGTTGAAAATTATTATGTGACATATTCATCAATGATA TACCGTCAGCAGCAGTCAGGCCGAGGGTGGTATCTGGGTCTGAACAAAGAAGGAGAGATC ATGAAAGGCAACCATGTGAAGAAGAACAAGCCTGCAGCTATTTTCTGCCTAAACCACTG AAAGTGGCCATGTACAAGGAGCCATCACTGCACGATCTCACGGAGTTCTCCCGATCTGGA AGCGGGACCCCAACCAAGAGCAGAAGTGTCTCTGGCGTGTGAACGGAGGCAAATCCATG AGCCACAATGAATCAACG
Restriction Sites:	Please inquire
ACCN:	NM_001139502
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001139502.1](#), [NP_001132974.1](#)

RefSeq Size: 2093 bp

RefSeq ORF: 681 bp

Locus ID: 2258

UniProt ID: [Q92913](#)

Cytogenetics: Xq26.3-q27.1

Protein Families: Secreted Protein

Protein Pathways: MAPK signaling pathway, Melanoma, Pathways in cancer, Regulation of actin cytoskeleton

Gene Summary: The protein encoded by this gene is a member of the fibroblast growth factor (FGF) family. FGF family members possess broad mitogenic and cell survival activities, and are involved in a variety of biological processes, including embryonic development, cell growth, morphogenesis, tissue repair, tumor growth, and invasion. This gene is located in a region on chromosome X, which is associated with Borjeson-Forssman-Lehmann syndrome (BFLS), making it a possible candidate gene for familial cases of the BFLS, and for other syndromal and nonspecific forms of X-linked cognitive disability mapping to this region. Alternative splicing of this gene at the 5' end results in several transcript variants encoding different isoforms with different N-termini. [provided by RefSeq, Nov 2008]
Transcript Variant: This variant (5) contains two additional exons at the 5' end compared to transcript variant 1. This results in a slightly shorter isoform (3) with a different N-terminus compared to isoform 1. Variants 3 and 5 encode the same isoform.