

Product datasheet for **SC325651**

TOR2A (NM_001134430) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: TOR2A (NM_001134430) Human Untagged Clone
Tag: Tag Free
Symbol: TOR2A
Synonyms: TORP1
Vector: pCMV6 series

Fully Sequenced ORF: >NCBI ORF sequence for NM_001134430, the custom clone sequence may differ by one or more nucleotides

```
ATGGCGGCTGCGACGCGGGCTGCCGGCCCTGGGGCTCGCTCCTCGGGCTGCTCGGGCTG
GTCTCGGCCGCGGCCGCCGCTGGGACCTGGCTTCCCTGCGCTGCACCTGGGCGCCTTT
TGCGAATGCGACTCCGGCCCGACTTGGCGGGTCTGGAGTGTGACCTGGCTCAGCACCTG
GCCGGCCAGCATCTGGCCAAGGCGCTGGTGGTGAAGGCGCTGAAGGCCTTTGTGCGGGAC
CCAGCCCCACCAAGCCGCTGGTCCTCTCCCTGCACGGCTGGACCGGCACCGGCAAATCC
TATGTCAGCTCCCTGCTGGCGCACTACCTTCCAGGGCGGCCTCCGAGCCCCCGCTG
CACCACCTTTCTCCCGTCTCCACTTCCCCACCCAGCCACATCGAGCGCTACAAGAAG
GATCTGAAGAGCTGGGTCCAAGGGAACCTCACTGCCTGTGGCCGCTCCCTTCTCCTTTC
GATGAGATGGACAAGATGCCCCAGGCCTGATGGAAGTCTGCGGCCCTTCTGGGCTCC
TCTGGGTGGTATACGGGACCAATTACCGCAAAGCCATCTTCATCTTCATCAGATGGCTT
CTCAAACCTCGGGCATCATGGAAGAGCGCCTCTAGACGAGTGGTGGCCCTTCTCCCGCT
CCAGCGGCACCACGTCCGGCACTGCGTGCTCAACGAGCTGGCCAGCTGGGCCTGGAGCC
AAGGGA
```

Restriction Sites: Please inquire

ACCN: NM_001134430

Insert Size: 1460 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001134430.1</u> , <u>NP_001127902.1</u>
RefSeq Size:	1460 bp
RefSeq ORF:	729 bp
Locus ID:	27433
UniProt ID:	<u>Q5JU69</u>
Cytogenetics:	9q34.11
Protein Families:	Secreted Protein, Transmembrane
Gene Summary:	<p>This gene encodes a member of the AAA family of adenosine triphosphatases with similarity to Clp proteases and heat shock proteins. Alternative splicing at this locus results in the translation of multiple isoforms of the encoded protein, some of which contain salusin peptides in the C-terminal region. These peptides may play roles in hypotension, myocardial growth and the induction of mitogenesis, and may also be involved in the pathogenesis of atherosclerosis. The antimicrobial peptide salusin-beta has antibacterial activity. [provided by RefSeq, Nov 2014]</p> <p>Transcript Variant: This variant (3) lacks an exon in the 3' coding region, which results in a frameshift, compared to variant 1. The encoded isoform (c) is shorter and has a distinct C-terminus, compared to isoform a.</p>