

Product datasheet for **SC325655**

AZ12 (NM_001134432) Human Untagged Clone

Product data:

| | |
|----------------------|---|
| Product Type: | Expression Plasmids |
| Product Name: | AZ12 (NM_001134432) Human Untagged Clone |
| Tag: | Tag Free |
| Symbol: | AZ12 |
| Synonyms: | AZ2; NAP1; TILP |
| Vector: | <u>pCMV6 series</u> |
| Fully Sequenced ORF: | >NCBI ORF sequence for NM_001134432, the custom clone sequence may differ by one or more nucleotides ATGGATGCACTGGTAGAAGATGATATCTGTATTCTGAATCATGAAAAAGCCATAAGAGA GATACAGTGACTCCAGTTTCAATATATTCAGGAGATGAATCTGTTGCTTCCCATTGCT CTTGTCACATGATGAAGACATCAAAAAACGACTTAAGGATTCAGAGAAAGAGAACTCT TTGTTAAAGAAGAGAATAAGATTTTTGGAAGAAAAGCTAATAGCTCGATTTGAAGAAGAA ACAAGTCCGTGGGACGAGAACAAGTAAATAAGGCCTATCATGCATATCGAGAGGTTTGC ATTGATAGAGATAATTTGAAGAGCAAACCTGGACAAAATGAATAAAGACAACCTCTGAATCT TTGAAAGTATTGAATGAGCAGCTACAATCTAAAGAAGTAGAACTCTCCAGCTGAGGACA GAGGTGGAACCTCAGCAGGTGATGAGGAATTTAAATCCACCTTCATCAAACCTGGGAGGTG GAAAAGTTGAGCTGTGACCTGAAGATCCATGGTTTGAACAAGAGCTGGAACCTGATGAGG AAAGAATGTAGCGATCTCAAATAGAACTACAGAAAGCCAAACAAACGGATCCATATCAG GAAGACAATCTGAAGAGCAGAGATCTCCAAAACCTAAGCATTTCAAGACAGAGTATCTGT AGCACAGCTGGAGTGCAGTGGCGCAATCATGGGATAGCATGAATTTTGAGGACTGTTCC CTCTTCGCA |
| Restriction Sites: | Please inquire |
| ACCN: | NM_001134432 |
| OTI Disclaimer: | Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP). |
| OTI Annotation: | This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA. |
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001134432.1](#), [NP_001127904.1](#)

RefSeq Size: 1460 bp

RefSeq ORF: 732 bp

Locus ID: 64343

UniProt ID: [Q9H6S1](#)

Cytogenetics: 3p24.1

Protein Pathways: RIG-I-like receptor signaling pathway

Gene Summary: AZI2, or NAP1, contributes to the activation of NFkB (see MIM 164011)-dependent gene expression by activating IKK-related kinases, such as NAK (TBK1; MIM 604834) (Fujita et al., 2003 [PubMed 14560022]).[supplied by OMIM, Mar 2008]

Transcript Variant: This variant (2) differs in the 3' coding region and 3' UTR, compared to variant 1. The resulting isoform (b) has a distinct C-terminus and is shorter than isoform a.