

Product datasheet for SC325660

ADTRP (NM_001143948) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ADTRP (NM_001143948) Human Untagged Clone
Tag:	Tag Free
Symbol:	ADTRP
Synonyms:	ALG1L; C6orf105; dJ413H6.1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC325660 representing NM_001143948. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
 GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
 ATGACGAAGACTTCTACATGCATATACCACTTCCTTGTCTGAGCTGGTATACTTTCCTCAATTATTAC
 ATCTCACAGGAAGGAAAAGACGAGGTGAAACCCAAATCTTGGCAATGGTGCAAGGTGGAATATATG
 ACGCTGCTTAATCTGCTCAAGAACAGGACTGCTGGGTTTGACATCTACCAGCCAGGAAGCTTTAGGCAG
 CTCTTGACAGACCATTTTCTACGGGGTCACCTGCCTGGATGATGTGCTGAAAAGAACCAAGGGGAAAA
 GACATTAAGTTCCTAAGTGCCTTCAGAGACCTGCTTTTCACTCTGGCTTTTCTGTATCCACGTTT
 GTATTTTGGCATTCTGGATCCTTTTCTACAATCGAGATCTCATTTACCCCAAGGCTCTAGATACT
 GTCATCCCGTGTGGCTGAATCATGCAATGCACACTTTTATATTTCCCATCACATTGGCTGAAGTCGTC
 CTCAGGCCTCACTCCTATCCATCAAAGAAGACAGGACTCACCTTGCTGGCTGCTGCCAGCATTGCTTAC
 ATCAGCCGCATCCTATGGCTCTACTTTGAGACGGGTACCTGGGTGTATCCTGTGTTGCCAACTCAGC
 CTCTTGGGTCTAGCAGCTTTCTTCTCTCAGCTACGTTTCATCGCCAGCATCTACCTACTTGGAGAG
 AAGCTCAACCACTGGAATGGGGTGACATGAGGCAGCCACGGAAGAAGAGGAAGTAA
 ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
 TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC

Restriction Sites:	SgfI-MluI
ACCN:	NM_001143948
Insert Size:	747 bp


[View online »](#)

OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001143948.1
RefSeq Size:	1866 bp
RefSeq ORF:	747 bp
Locus ID:	84830
UniProt ID:	Q96IZ2
Cytogenetics:	6p24.1
Protein Families:	Transmembrane
MW:	28.9 kDa
Gene Summary:	<p>Regulates the expression and the cell-associated anticoagulant activity of the inhibitor TFPI in endothelial cells (in vitro).[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (1) represents the longer transcript and encodes the longer isoform (1).</p>