

Product datasheet for **SC325668**

CYB561D1 (NM_001134400) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: CYB561D1 (NM_001134400) Human Untagged Clone
Tag: Tag Free
Symbol: CYB561D1
Vector: pCMV6 series
Fully Sequenced ORF: >NCBI ORF sequence for NM_001134400, the custom clone sequence may differ by one or more nucleotides

```
ATGCAGCCCTGGAGGTAGGTCTGGTCCCGCTCCAGCTGGGGAGCCGAGACTGACCCGC  
TGGCTGCGGAGAGGCAGTGGGATCTTGGCGCACCTGGTAGCTTTGGGCTTACCATCTTT  
CTGACAGCGCTGTCCCGCCAGGAACAAAACAGGTCCCCTGATGGAGGATAGAAGTGAA  
GGAGGCCGGGCGGTCGGTTCATGCCTGAAATCCCAGCACTTTGGGAGGCCGACGCGGGT  
GGATCACTTGAGTTCTGCCTCTGCATGGCTGAAGCCATCCTACTTTCTCACCTGAACAC  
TCCCTGTTCTTCTTCTGCTCCCAGAAAAGCACGGATCCGGCTCCACTGGGCAGGGCAGACC  
CTAGCCATCCTCTGTGCAGCTCTGGGCCTGGGCTTCATCATCTCCAGCAGGACCCGAGT  
GAGCTGCCTCATCTGGTGTCTTGGCACAGCTGGGTGGGAGCCCTGACACTGCTGGCCACT  
GCTGTCCAGGCACTGTGTGGGCTCTGCCTCCTTTGCCCCGGGCAGCCAGGGTCTCAAGG  
GTGGCTCGCCTCAAGCTCTACCATCTGACATGTGGACTGGTGGTCTACCTGATGGCTACA  
GTAACGGTGCTTCTGGGCATGTAAGTATGGTTCAGGCCAGATCAAAGGTGCGGCC  
TGGTACCTGTGCCTGGCACTGCCCGTCTATCCAGCCCTGGTGATCATGCACCAGATTTCC  
AGATCCTACTTGCCGAGGAAGAAAATGGAATG
```

Restriction Sites: Please inquire
ACCN: NM_001134400
OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001134400.1](#), [NP_001127872.1](#)

RefSeq Size: 5010 bp

RefSeq ORF: 756 bp

Locus ID: 284613

UniProt ID: [Q8N8Q1](#)

Cytogenetics: 1p13.3

Protein Families: Transmembrane

Gene Summary: Probable transmembrane reductase that may use ascorbate as an electron donor and transfer electrons across membranes to reduce monodehydro-L-ascorbate radical and iron cations Fe(3+) in another cellular compartment.[UniProtKB/Swiss-Prot Function]
Transcript Variant: This variant (1) represents the longest transcript and it encodes the longest protein (isoform 1). Sequence Note: The RefSeq transcript and protein were derived from transcript and genomic sequence to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on alignments.