

Product datasheet for **SC325677**

ABHD11 (NM_001145364) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ABHD11 (NM_001145364) Human Untagged Clone
Tag:	Tag Free
Symbol:	ABHD11
Synonyms:	PP1226; WBSCR21
Vector:	<u>pCMV6 series</u>
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001145364, the custom clone sequence may differ by one or more nucleotides ATGCGAGCCGGCCAACAGCTTGCAAGCATGCTCCGCTGGACCCGAGCCTGGAGGCTCCCG CGTGAGGGACTCGGCCCCACGGCCCTAGCTTCGCGAGGGTGCCTGTCGCACCCAGCAGC AGCAGCGCGGCCGAGGGGGCGCCGAGCCGAGGCCGCTTCCGCTTCTACAGGCTTCTG GACGGGGAGGCAGCCCTCCCGGCCGTCGTCTTTTGCACGGGCTTTCGGCAGCAAACT AACTTCAACTCCATCGCCAAGATCTTGGCCAGCAGACAGGCCGTAGGGTGTGACGGTG GATGCTCGTAACCACGGTGACAGCCCCACAGCCAGACATGAGCTACGAGATCATGAGC CAGGACCTGCAGGACCTTCTGCCCCAGCTGGGCCTGGTGCCTGCGTCGTTGGCCAC AGCATGGGAGGAAAGACAGCCATGCTGCTGGCACTACAGAGGGACATGGCCGTGCCGAG CACCTGCTCACTAACCTGGTAGAGGTAGACGGGCGCTTCGTGTGGAGGGTGAACCTGGAT GCCCTGACCCAGCACCTAGACAAGATCTTGGCTTTCCACAGAGGCAGGAGTCTACCTC GGGCCAACACTCTTTCTCCTTGGTGAAACTCCCAGTTCGTGCATCCCAGCCACCACCT GAGATTATGCGGCTCTTCCCTCGGGCCAGATGCAGACGGTGCCGAACGCTGGCCACTGG ATCCACGCTGACCGCCACAGGACTTCATAGCTGCCATCCGAGGCTTCTCTGGTC
Restriction Sites:	Please inquire
ACCN:	NM_001145364
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001145364.1](#), [NP_001138836.1](#)

RefSeq Size: 1300 bp

RefSeq ORF: 777 bp

Locus ID: 83451

UniProt ID: [Q8NEV4](#)

Cytogenetics: 7q11.23

Gene Summary: This gene encodes a protein containing an alpha/beta hydrolase fold domain. This gene is deleted in Williams syndrome, a multisystem developmental disorder caused by the deletion of contiguous genes at 7q11.23. [provided by RefSeq, Mar 2016]
Transcript Variant: This variant (8) lacks a coding exon, compared to variant 1, resulting in a shorter protein (isoform 8).