

## **Product datasheet for SC325683**

## HSD17B13 (NM 001136230) Human Untagged Clone

**Product data:** 

**Product Type:** Expression Plasmids

Product Name: HSD17B13 (NM\_001136230) Human Untagged Clone

Tag: Tag Free
Symbol: HSD17B13

Synonyms: HMFN0376; NIIL497; SCDR9; SDR16C3

Mammalian Cell Neomycin

Selection:

**Vector:** pCMV6-Entry (PS100001)

E. coli Selection: Kanamycin (25 ug/mL)

Fully Sequenced ORF: >SC325683 representing NM\_001136230.

Blue=Insert sequence Red=Cloning site Green=Tag(s)

GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC

ATGAACATCATCCTAGAAATCCTTCTGCTTCTGATCACCATCATCTACTCTACTTGGAGTCGTTGGTG
AAGTTTTTCATTCCTCAGAGAGAAAATCTGTGGCTGGGGAGATTGTTCTCATTACTGGAGCTGGGCAT
GGAATAGGCAGGCAGACTACTTATGAATTTGCAAAAACGACAGAGCATATTGGTTCTGTGGGATATTAAT
AAGGTGAAGAAAGAAGTGGGTGATGTAACAATCGTGGTGAATAATGCTGGGACAGTATATCCAGCCGAT
CTTCTCAGCACCAAGGATGAAGAGATTACCAAGACATTTGAGGTCAACATCCTAGGACATTTTTGGATC
ACAAAAGCACTTCTTCCATCGATGATGGAGAGAAATCATGGCCACATCGTCACAGTGGCTTCAGTGTGC
GGCCACGAAGGGATTCCTTACCTCATCCCATATTGTTCCAGCAAATTTGCCGCTGTTGGCTTTCACAGA
GGTCTGACATCAGAACTTCAGGCCTTGGGAAAAACTGGTATCAAAACCTCATGTCTCTGCCCAGTTTTT
GTGAATACTGGGTTCACCAAAAAATCCAAGCACAAGATTATGGCCTGTATTGGAGACAGATGAAGTCGTA
AGAAGTCTGATAGATGGAATACTTACCAATAAGAAAATGATTTTTGTTCCATCGTATATCAATATCTTT
CTGAGACTACAGAAGTTTCTTCCTGAACGCCCTCAGCGATTTTAAATCGTATGCAGAAATTCAATTT

GAAGCAGTGGTTGGCCACAAAATCAAAATGAAA<mark>TGA</mark>

**ACGCGTACGCGCCCCTC**GAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT

TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC

**Restriction Sites:** Sgfl-Mlul

**ACCN:** NM 001136230

**Insert Size:** 795 bp



**OriGene Technologies, Inc.** 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com



**OTI Disclaimer:** 

Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at <a href="mailto:customercom">customercom</a> or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <u>More info</u>

**OTI Annotation:** 

This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components:

The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:** 

- 1. Centrifuge at 5,000xg for 5min.
- 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
- 3. Close the tube and incubate for 10 minutes at room temperature.
- 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
- 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** <u>NM 001136230.2</u>

 RefSeq Size:
 2289 bp

 RefSeq ORF:
 795 bp

 Locus ID:
 345275

 UniProt ID:
 Q7Z5P4

 Cytogenetics:
 4q22.1

**Protein Families:** Druggable Genome, Secreted Protein, Transmembrane

**MW:** 29.6 kDa